



Jan 30, 2024

Biotinylation by antibody recognition

DOI

dx.doi.org/10.17504/protocols.io.14egn3qq6l5d/v1

Bryan_Killinger¹

¹Rush University Medical Center



Bryan Killinger

Rush University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn3qq6l5d/v1>

Protocol Citation: Bryan_Killinger 2024. Biotinylation by antibody recognition. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.14egn3qq6l5d/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 16, 2024



Last Modified: January 30, 2024

Protocol Integer ID: 94351

Keywords: biotinylation by antibody recognition, biotinylation, antibody recognition, antibody

Funders Acknowledgements:

NCI

Grant ID: CCSG P30 CA060553

NIH Office of Director

Grant ID: S10OD025194

National Resource for Translational and Developmental Proteomics

Grant ID: P41 GM108569

EraPerMed DEEPEN-iRBD project

Grant ID: ANR-22-PERM-0006

Michael J. Fox Foundation

Grant ID: ASAP-000458

NINDS

Grant ID: 1R01NS128467

Michael J. Fox Foundation

Grant ID: ASAP-021030

NIH

Grant ID: R21 NS109871

NINDS

Grant ID: K23-NS097625-06

Abstract

This protocol details the biotinylation by antibody recognition.

Attachments



[956-2485.docx](#)

15KB

Materials

Crosslink reversal buffer

	A	B
	SDS	5%
	Tris-HCl pH 8.0	500 mM
	NaCl	150 mM
	EDTA	2 mM

Modified TBST

	A	B
	Tris-HCl	20 mM
	NaCl	200 mM
	EDTA	2 mM
	Triton X-100	0.5%

Stringent wash buffer

	A	B
	Tris-HCl pH 7.6	20 mM
	NaCl	200 mM
	SDS	0.1%
	EDTA	2 mM

TBST

	A	B
	Tris-HCl pH 7.6	20 mM
	NaCl	150 mM




	A	B
	Tween-20	0.1%

High stringency wash buffer

	A	B
	Tris-HCl pH 7.6	20 mM
	NaCl	400 mM
	Tween-20	0.1%

 Trypsin **Promega Catalog #V5111**

 Clarity Western ECL Substrate **Bio-Rad Laboratories Catalog #1705060**

Troubleshooting



Biotinylation by antibody recognition

7h 45m

- 1 Collect the brain sections at 240-micron intervals across the neuroaxis, place them into a net well (Brain research laboratories) and wash 3 times for 01:00:00 each in TBST. 1h
- 2 Place the sections in 0.3% hydrogen peroxide and 0.1% sodium azide diluted in blocking buffer for 01:00:00 at Room temperature to quench endogenous peroxidases. 1h
- 3 Rinse the sections briefly in TBST and incubate in anti-PSER129 antibody EP1536Y diluted 1:50,000 in blocking buffer Overnight at 4 °C with gentle agitation. 1h
- 4 The following day, wash the sections 3 times in TBST, then incubate with biotinylated anti-rabbit antibody diluted 1:200 in blocking buffer for 01:00:00 at Room temperature 1h
- 5 Wash the sections 3 times in TBST, incubate with ABC reagent for 01:00:00 , and wash off with borate buffer. 1h
- 6 Incubate the sections with borate buffer containing biotinyl tyramide as described above.
- 7 Wash the sections Overnight with TBST, gather in a 1.5mL Eppendorf tube, 3000 x g, 00:15:00 to pellet floating sections, and discard the supernatant. 1h 15m
- 8 Briefly sonicate each sample in 1 mL of crosslink reversal buffer (refer materials section) and heat for 00:30:00 at 98 °C followed by 01:00:00 at 90 °C . 1h 30m
- 9 Centrifuge the sample 20000 x g, 00:20:00 of the samples and then dilute the supernatant 1:10 in modified TBST (refer materials section). 20m











- 10 Incubate each sample with 40 mg of streptavidin magnetic beads (Thermofisher Scientific) for 02:00:00 at Room temperature with constant mixing. 2h
- 11 Collect the beads using a magnetic stand (Thermofisher Scientific), wash the beads 3 times in modified TBST, and then Overnight in 10 mL of stringent wash buffer (refer materials section). 2h
- 12 The following day, collect the beads using magnetic stand and resuspend in 100 μ L 1 X Bolt LDS sample buffer with reducing agent (Thermofisher) then heat for 00:10:00 at 98 $^{\circ}$ C . 10m
- 13 Vortex the samples vigorously and remove the beads using magnetic stand.
- 14 Subject 70 μ L of the sample to electrophoresis approximately 2 cm into a Bolt gel (ThermoFisher).
- 15 Fix the gel in 50% ethanol and 10% acetic acid for 01:00:00 . 1h
- 16 Wash the gel several times in dH₂O, and stain the proteins with colloidal Coomassie blue.
- 17 Then excise the entire sample for trypsin digestion and mass spectrometry.
- 18 Wash the gel pieces with 100 millimolar (mM) ammonium bicarbonate (AmB)/acetonitrile (ACN) and reduce with 10 millimolar (mM) dithiothreitol (DTT) at 50 $^{\circ}$ C for 00:45:00 . 45m
- 19 Alkylate the cysteines using 100 millimolar (mM) iodoacetamide in the dark for 00:45:00 at Room temperature (RT). 45m




- 20 Wash the gel bands in [M] 100 millimolar (mM) AmB/ACN prior to adding 1 µg trypsin (Promega #V5111) for Overnight incubation at 37 °C . 45m
- 21 Collect the peptide containing supernatants into a separate tube.
- 22 Wash the gel pieces with gentle shaking in 50% ACN/1% FA at Room temperature for 00:10:00 , and collect the supernatant in the previous tubes. 10m
- 23 Do the final peptide extraction step with 80% ACN/1% FA, and 100% ACN, and collect all supernatant.
- 24 Dry the peptides in a speedvac and reconstitute with 5% ACN/0.1% FA in water before injecting into LC-MS/MS.
- 25 Analyse the peptides by LC-MS/MS using a Dionex UltiMate 3000 Rapid Separation nanoLC coupled to an Orbitrap Elite Mass Spectrometer (Thermo Fisher Scientific Inc.).
- 26 Load the samples onto the trap column, which is 150 µm x 3 cm in-house packed with 3 µm ReproSil-Pur® beads.
- 27 The analytical column is a 75 µm x 10.5 cm PicoChip column packed with 3 µm ReproSil-Pur® beads (New Objective, Inc. Woburn, MA).
- 28 Keep the flow rate at 300 nL/min.
- 29 Ellute all the fractions from the analytical column at a flow rate of 300 nL/min using an initial gradient elution of 5% B from 00:00:00 to 00:05:00 , transition to 40% over 01:40:00 , 60% for 00:04:00 , ramping up to 90% B for 00:03:00 , holding 90% B for 00:03:00 , followed by re-equilibration of 5% B at 00:10:00 with a total run time of 02:00:00 . 4h 5m
- 30 Record the mass spectra (MS) and tandem mass spectra (MS/MS) in positive-ion and high-sensitivity mode with a resolution of ~60,000 full-width half-maximum.



- 31 Select the 15 most abundant precursor ions in each MS1 scan for fragmentation by collision-induced dissociation (CID) at 35% normalized collision energy in the ion trap.
- 32 Dynamically excluded the previously selected ions from re-selection for  00:01:00 . 1m
Store the collected raw files spectra in. raw format.
- 33 Identify the proteins from the MS raw files using the Mascot search engine (Matrix Science, London, UK. version 2.5.1). 
- 34 Search the MS/MS spectra against the SwissProt mouse database. 
- 35 Include carbamidomethyl cysteine as a fixed modification and oxidized methionine, deamidated asparagine and aspartic acid, and acetylated N-terminal as variable modifications in all searches.
- 36 Allow three missed tryptic cleavages. Apply a 1% false discovery rate cutoff at the peptide level.
- 37 Consider only proteins with a minimum of two peptides above the cutoff for further study.
- 38 Visualize the identified peptides/protein by Scaffold software (version 5.0, Proteome Software Inc., Portland, OR). 
- 39 To estimate BAR enrichment, apply  1 μ L of bead eluent to a methanol activated polyvinylidene difluoride (PVDF) membrane and then allow to dry completely.
- 40 Reactivate the membrane then in methanol, rinse with water, and post-fix in 4% PFA for  00:30:00 . 30m
- 41 Rinse the blots with TBST (refer materials section) and block with buffer containing either BSA (TBST and 5% BSA) or non-fat milk (TBST and 5% non-fat milk) for detection of biotin or α syn , respectively.
- 42 Detect the biotinylated proteins by ABC (VectorLabs) diluted 1:10 in BSA blocking buffer for  01:00:00 at  Room temperature . 1h



- 43 Asyn can be detected using SYN1 (BD Biosciences) diluted 1:2,000 and PSER129 detected using EP1536Y diluted 1:50,000 both diluted in non-fat milk blocking buffer.
- 44 Detect the primary antibodies by incubating blots for  01:00:00 in secondary anti-mouse HRP conjugate diluted 1:6,000 or secondary anti-rabbit HRP conjugate (Cell signaling) diluted in milk blocking buffer.
- 45 Following secondary antibody, wash the membranes in high stringency wash buffer (Refer materials section) and image using enhanced chemiluminescence (ECL) substrate (Biorad, product # 1705060) and Chemidoc imager (Biorad).

1h

