

Oct 23, 2025

BL-21 Competency CaCl₂ Treatment

DOI

<https://dx.doi.org/10.17504/protocols.io.kqdg3119el25/v1>

mpalacios¹

¹Liberty University



mpalacios

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg3119el25/v1>

Protocol Citation: mpalacios 2025. BL-21 Competency CaCl₂ Treatment. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.kqdg3119el25/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 22, 2025



Last Modified: October 23, 2025

Protocol Integer ID: 230520

Keywords: Competency, competency cac12 treatment protocol for calcium chloride treatment, competency cac12 treatment protocol, calcium chloride treatment, competency

Abstract

Protocol for Calcium Chloride treatment to induce competency in BL-21 cells.

Materials

PER COLONY -

55 1.5mL sterile microcentrifuge tubes

5mL of sterile 100mM CaCl₂

2mL of sterile 100mM CaCl₂ + Glycerol

1 plate of non-selective media (LB Agar or TSA)

60mL of sterile LB (10mL in a falcon, 50mL in a 250mL shaking flask)

Microcentrifuge with temperature control (pre-chilled to 0°C)

Cold room is helpful

Ice Water bath

Safety warnings

 Due to working in non-selective media, ensure to follow sterile protocol at all steps.



Cell Growth

4h 15m

- 1 Plate your BL-21 on non-selective plates (LB Agar or TSA) 20m
- 2 Incubate overnight at 37 °C
- 3 Using a sterile pipette tip, pick a colony and add to 10mL of non-selective LB 20m
- 4 Incubate overnight at 140 rpm, 37°C, 12:00:00
- 5 Gather 55 sterile 1.5mL tubes and add to freezer to chill
- 5.1 Consider labeling the tubes at this point to save time later
- 6 Add 500 µL of overnight culture to 50 mL of non-selective LB 5m
- 7 Incubate 210 rpm, 37°C, 02:00:00 2h
- 7.1 You can also label tubes now.
- 8 Measure OD600 every 00:30:00 until early to mid-log phase (0.35-0.55) 1h 30m
- 9 Take colony and add to wet ice to halt growth.
- 9.1 Everything from this point forward needs to be as close to 0°C as possible.



CaCl₂ Treatment

3h 50m

- 10 Gather 5 prechilled 1.5mL tubes
- 11 Add 1mL of culture to each tube
- 12 Centrifuge the tubes at  8000 rpm, 0°C, 00:05:00 5m

- 13 Discard supernatant
- 14 Repeat steps 11-13 9 times until a total of 10mL worth of cells is added to each tube 1h
- 15 Resuspend each tube in 1mL of  100 millimolar (mM) CaCl₂ by pipetting 10m
- 16 Incubate on ice for  02:30:00 2h 30m

- 17 Centrifuge the tubes at  8000 rpm, 0°C, 00:05:00 5m

- 18 Discard supernatant and resuspend pellet in 500mL of  100 millimolar (mM) CaCl₂ +
 15 % volume glycerol
- 18.1 Last time pipetting up and down is possible, from here on out mix by flicking
- 19 Aliquot  50 µL into pre-chilled tubes.