



Jan 03, 2025

Bleach life-stage synchronization of *C. elegans*

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vzj91xlx1/v1>

Justin Donnelly¹

¹Arcadia Science

Arcadia Science



Arcadia Science

Arcadia Science

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzj91xlx1/v1>

Protocol Citation: Justin Donnelly 2025. Bleach life-stage synchronization of *C. elegans* . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.rm7vzj91xlx1/v1>

Manuscript citation:

Avasthi P, Borges AL, Cheveralls K, Donnelly J, Mets DG, Reiter T. (2025). An experimental and computational workflow to characterize nematode motility behavior. <https://doi.org/10.57844/arcadia-b89a-7e76>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 23, 2024

Last Modified: March 14, 2025

Protocol Integer ID: 106302

Keywords: *C. elegans*, bleach, synchronization, life-stage, L1, nematode, bleach life, bleach, bleaching gravid adult, hatched worm, greater resistance to bleach, stage synchronization, necessary for cryopreservation, cryopreservation

Abstract

This protocol describes the procedure for synchronizing the life stage of *C. elegans* worms by bleaching gravid adults so that they release their eggs, which have greater resistance to bleach than hatched worms. This procedure can be useful for eliminating the variable of life stage from downstream experiments, and is also necessary for cryopreservation of this organism.

Materials

 Potassium hydroxide **P212121**

 Ethanol 200 Proof **Decon Labs Catalog #2716**

Concentrated bleach

Autoclaved ultrapure water

M9 buffer

15 mL centrifuge tubes

NGM(OP50) plates

Clinical centrifuge

Tube rocker for 15 mL centrifuge tubes

Dissection scope

Protocol materials

 Potassium hydroxide **P212121**

 Ethanol 200 Proof **Decon Labs Catalog #2716**



Preparation of plates

- 1 Chunk worms onto at least six NGM (OP50) plates: Dip a spatula in 100%  Ethanol 200 Proof **Decon Labs Catalog #2716** . Flame-sterilize. Once the flame is extinguished, cut a chunk from the plate with an edge length of ~  0.5 cm . Select a chunk with several larvae. Transfer this chunk to a new plate under flame, with the worm side facing down.
- 2 Incubate worms for 3–5 days, until there are many gravid adults.

Note

Synchronization works best on young adults, which still have many fertilized eggs. The older adults become, the fewer L1s there will be at the end of the protocol.

Synchronization via alkaline bleach treatment

12m

- 3 Under a flame, add  0.5-1 mL autoclaved water to a plate. Swirl and tilt plate to ensure that water passes over the entirety of the plate, then use a pipette to transfer the water to a sterile 15 mL centrifuge tube. This process should remove the vast majority of worms from the plate.
- 4 Repeat step 3 with remaining plates. The total volume should be below  8 mL .
- 5 Once you've transferred the worms, add autoclaved water to the conical tube to a total volume of  8 mL . If synchronizing multiple strains, synchronize one at a time.
- 6 Add  1.2 mL concentrated bleach and  0.6 mL  5 Molarity (M)  Potassium hydroxide **P212121** . Mix solution by gently inverting.
- 7 Place the tube on a rocker. Monitor under a dissecting microscope to observe lysis of adults and release of fertilized eggs. This process should take  00:04:00 –  00:08:00 .

12m



Note

It's important that you don't over-bleach the sample. Extended exposure to alkaline-bleach conditions can damage eggs and prevent hatching. In our hands, our most successful bleach treatments were under 5 minutes. For best results, bias to incomplete lysis over excessive lysis.

- 8 Fill tube to  15 mL with autoclaved water.
- 9 Centrifuge in a swinging bucket centrifuge at  768.6 x g for  00:01:00 –  00:02:00 . 3m
- 10 Decant supernatant. Perform two additional washes with  15 mL autoclaved water.
- 11 Decant supernatant. Perform a third wash with  15 mL M9 buffer.
- 12 Resuspend worms in  1-2 mL M9 buffer. Place on rocker  Overnight at  20 °C to allow worms to hatch into L1s. Worms can remain on rocker for an additional day without negative effects. 8m

Plating and counting of L1s

- 13 *Optional:* To estimate the concentration of L1s, plate several volumes (e.g.,  5 µL ,  10 µL ,  20 µL) on plates and count the number of L1 larvae.
- 14 Under a flame, plate desired quantity of worm larvae on desired number of NGM (OP50) plates by pipetting appropriate volume of buffer onto the surface of the plate. Allow buffer to evaporate before returning to incubator.

Protocol references

M9 buffer preparation:

<https://www.protocols.io/view/preparation-of-m9-worm-buffer-cckysuxw.html>