

Aug 26, 2025

## Blood DNA extraction protocol compatible with LAMP-PCR for *T. cruzi* detection.

 [Microbiology Spectrum](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vz9mnrgx1/v1>

Sneider Gutierrez<sup>1</sup>, Jessy Condori<sup>2</sup>

<sup>1</sup>Jons Hopkins University; <sup>2</sup>Universidad Peruana Cayetano Heredia



**Sneider Gutierrez**

Johns Hopkins Bloomberg School of Public Health

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.rm7vz9mnrgx1/v1>

**External link:** <https://doi.org/10.1128/spectrum.00172-25>

**Protocol Citation:** Sneider Gutierrez, Jessy Condori 2025. Blood DNA extraction protocol compatible with LAMP-PCR for *T. cruzi* detection.. protocols.io <https://dx.doi.org/10.17504/protocols.io.rm7vz9mnrgx1/v1>

**Manuscript citation:**

Guarnizo SAG, Condori BJ, Basma L, Equilia S, Malaga E, Defazio S, Arteaga E, Velarde JK, Obregón M, Takyar A, Duque C, Hakim J, Tinajeros F, Gilman RH, Bowman N, Mugnier MR A specific, stable, and accessible LAMP assay targeting the HSP70 gene of *Trypanosoma cruzi*. Microbiology Spectrum 13(11). doi: [10.1128/spectrum.00172-25](https://doi.org/10.1128/spectrum.00172-25)

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development

**We are still developing and optimizing this protocol**

**Created:** June 28, 2025

**Last Modified:** August 26, 2025

**Protocol Integer ID:** 221226

**Keywords:** Blood DNA extraction, cost-effective method, *T. cruzi* detection, compatible with LAMP-PCR and qPCR, blood dna extraction protocol compatible with lamp, blood dna extraction protocol, isolated dna, quality nucleic acids from fresh peripheral blood, pcr, fresh peripheral blood, reliable results for downstream molecular application, ctab lysis buffer, qpcr

## Abstract

This cost-effective protocol uses CTAB lysis buffer and isopropanol precipitation to purify high-quality nucleic acids from fresh peripheral blood collected in EDTA. The protocol is about four times cheaper when compared with commercial kits. The isolated DNA is compatible with both LAMP-PCR and qPCR assays, providing reliable results for downstream molecular applications.

## Guidelines

Compared to commercial kits, this protocol yields a comparable average DNA amount but with moderately higher variability.

## Materials

### Reagents

|  | A                                                                 | B                 | C                  |
|--|-------------------------------------------------------------------|-------------------|--------------------|
|  | Reagent                                                           | Company           | Catalog number/SKU |
|  | 2-propanol                                                        | sigma aldrich     | 1096341011         |
|  | Hexadecyltrimethylammonium bromide                                | sigma aldrich     | H5882-100G         |
|  | Tris-HCl, 1M Solution, pH 8.0, Molecular Biology Grade, Ultrapure | Fisher scientific | AAJ22638AE         |
|  | Ethylenediaminetetraacetic acid solution                          | MERK              | E8008-100ML        |
|  | Absolute ethanol                                                  | Thermo scientific | T038181000         |
|  | Nuclease free water                                               | Milipore          | 4.86505.1000       |
|  | NaCl (5 M), RNase-free                                            | Milipore          | S9625-1Kg          |

### Disposable materials

|  | A                                                        | B              | C                  |
|--|----------------------------------------------------------|----------------|--------------------|
|  | Disposable material                                      | Company        | Catalog number/SKU |
|  | Conical Base 2.0 mL Microcentrifuge Tubes                | USA Scientific | 14202700           |
|  | 1000 uL tips (TipOne <sup>®</sup> Pipette Tips in Racks) | USA Scientific | 1111-0806          |
|  | 200 uL tips (TipOne <sup>®</sup> Pipette Tips in Racks)  | USA Scientific | 1111-0806          |

### Permanent materials

|  | A                             | B   |
|--|-------------------------------|-----|
|  | Permanent materials/equipment | SKU |



| A                                                                                              | B         |
|------------------------------------------------------------------------------------------------|-----------|
| Vortex mixer                                                                                   | CLS6885DB |
| Centrifuge (capable of 12,000 × g) (Benchmark MC-12 High Speed Microcentrifuge)                | BMS:C1612 |
| Heater (prewarmed to 70 °C)<br>(Corning <sup>®</sup> LSE <sup>®</sup> digital dry bath heater) | CLS6885DB |
| Plastic Test Tube Rack Cube                                                                    | 76263-262 |

## Safety warnings

- ! It is recommended to receive training before processing important samples and to perform extractions in replicates when possible.

## Before start

Collect peripheral blood in an EDTA vacutainer and proceed with the checklist section.



## Check List (Ensure you have the following materials ready before running the experiment)

10m 30s

### 1 Equipment

30s

- Vortex mixer
- Centrifuge (*capable of 12,000 × g*)
- Heater (*prewarmed to 70 °C*)
- Timer

### 2 Cold Storage

2m

- A bucket with ice

### 3 Reagents & Buffers

30s

- CTAB 2% lysis buffer (*CTAB: 2% (w/v), Tris-HCl (pH 8.0): 100 mM, EDTA: 20 mM, NaCl: 1.4 M.*)
- 5 M NaCl (*cold/onice*)
- Isopropanol (*cold/onice*)
- 70% ethanol
- PCR-grade nuclease-free water (*incubate at 70 °C*)

### 4 Consumables & Labware

2m

- Rack for 2 mL tubes
- 2 mL microcentrifuge tubes (*Catalog 1420-2700 USA Scientific; two per sample*)  
Do not use low-binding tubes.
- 1 mL pipette + tips
- 100 µL pipette + tips
- Paper towel
- Dampen a paper towel with 10% bleach to clean gloves between samples

### 5 Labeling & Organization

5m

- Markers for labeling
- Label the tubes as follows:
  - **Left side:** Date and operator's initials
  - **Top:** Sample order (e.g., 1, 2, 3) + "NA" =nucleic acids + full sample name
  - **Right side:** Leave blank for DNA concentration in case you have the equipment to quantify

### 6 Are you running the LAMP-PCR right after? If so,...

30s



- Make sure you have enough master mix vials pre-stored at -20 °C.
- Thaw the SYBR Green at 4 °C.



## Cell lysis

12m 25s

7

Add **300mL of fresh whole blood in an EDTA** to 2mL microcentrifuge tube

10s



8

Add **500 µL of CTAB lysis** buffer to the whole blood.

10s



9

**5 seconds Vortex** thoroughly to mix.

5s



10

Incubate at **65°C for 10** minutes

10m 30s



11

**1 min on ice** to improve protein precipitation

1m 30s



## Protein precipitation

3m 10s

12

Add **500 µL of 5 M** NaCl to the lysed sample by pipetting three times.

10s



13

**5 seconds Vortex**

10s



14

Centrifuge at **12,000x g for 2** minutes.

2m 30s



15

Carefully **pipette 1 mL of the supernatant** from top to bottom, ensuring that the pellet is not touched.

10s



16

**Transfer** the 1 mL supernatant to a **new, pre-labeled 2 mL tube**.

10s



## Nucleic acids precipitation

7m 40s

- 17 Add **1mL of cold isopropanol** to the recovered supernatant

10s



- 18 **Gently invert** to mix 20 times

1m 30s

**Note:** At this point, you should see debris floating in the solution.



***Trick:** If you are working with multiple samples, place the tubes in the rack. Use the palms of your hands to make a sandwich, and gently rotate from right to left to mix the isopropanol*

- 19 ***Trick:** If you did not see precipitation during the previous step, incubate the sample on ice or -20 °C for 10 minutes before proceeding.*



- 20 Centrifuge at **12,000xg for 5 minutes**

5m 30s



- 21 Carefully **remove the supernatant** by pipetting, making sure not to disturb the pellet stuck to the side at the bottom of the tube. Since the DNA is concentrated in the pellet, you can safely discard any blood debris present in the upper phases.

30s

## Nucleic acids clean up and DNA quantification

18m 43s

- 22 Add **1 mL of 70% ethanol** to wash the pellet

10s

- 23 **Vortex** for 1 second.

3s

- 24 Centrifuge at **12,000x g for 2 minutes**

2m 30s

- 25 **Discard** the ethanol by inverting the tube and **tapping** it repeatedly on a paper towel until any remaining ethanol drains down the side.



**Note:** It is recommended to invert the tube upside down only once. While inverting, gently rotate the tube to prevent the liquid from putting pressure on one side of the tube.

- 26 Add **1 mL of 100% ethanol** to wash the pellet 10s
- 27 **Discard** the ethanol by inverting the tube and **tapping** it repeatedly on a paper towel until any remaining ethanol drains down the side. 20s
- 28 Dry the DNA for **3 minutes at 70 °C**. 3m 10s
- 29 **Resuspend** the dried pellet in preheated 25-50 µL by pipetting 20 times. 30s
- 30 **Incubate at 70 °C for 5 minutes** to improve DNA dissolution 5m 10s
- 31 **Vortex** for 5 seconds 10s
- 32 **Centrifuge** at 12000 for 1 minute 1m 10s
- 33 Quantify DNA by **NanoDrop**, taking DNA from the top of the solution. 5m

## Sample storage

- 34 **Store** at 4°C for short-term (< 1week) or -20°C for long-term storage (>1 week). 1m

## Running LAMP- PCR

- 35 Take out the LAMP-PCR mix stored at -20°C and spin it down. **Note:** include a negative and positive control. 1m



36 Add **2  $\mu$ L of SYBR Green** to the lid.

10s



37 Load 5  $\mu$ L of pre-extracted DNA and pipette 5 times.

20s



38 Incubate for 1 hour at 65 °C.

1h 20m



39 Spin it down and record the result: brown = negative, green = positive.

30s

