



Oct 02, 2025

Blunt-End-Single-Tube (BEST): Adapter Hybridization

 Ecology and Evolution

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlkym4wg5r/v1>

Michael V. Westbury¹, Alba Rey-Iglesia¹, Deon de Jager¹, Vanssy Li¹, Eline D. Lorenzen¹

¹University of Copenhagen

Lorenzen Lab



Vanssy Li

Copenhagen university

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlkym4wg5r/v1>

External link: <https://doi.org/10.1111/2041-210X.12871>

Protocol Citation: Michael V. Westbury, Alba Rey-Iglesia, Deon de Jager, Vanssy Li, Eline D. Lorenzen 2025. Blunt-End-Single-Tube (BEST): Adapter Hybridization. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.j8nlkym4wg5r/v1>

**Manuscript citation:**

Carøe, C., Gopalakrishnan, S., Vinner, L., T. Mak, S. S., S. Sinding, M. H., Samaniego, J. A., Wales, N., Sicheritz-Pontén, T., & P. Gilbert, M. T. (2018). Single-tube library preparation for degraded DNA. *Methods in Ecology and Evolution*, 9(2), 410-419. <https://doi.org/10.1111/2041-210X.12871>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 18, 2025

Last Modified: October 02, 2025

Protocol Integer ID: 227597

Keywords: Ancient DNA, Library Build, Adapter Hybridization, BEST, tube library preparation for degraded dna, tube library preparation, adapter hybridization this protocol, process of the hybridization, degraded dna, adapter hybridization, hybridization, p7 adapters for the blunt, dna, tube, p7 adapter

Abstract

This protocol describes the process of the hybridization of P5 and P7 adapters for the Blunt-End-Single-Tube (BEST) library build protocol. Following this protocol produces the P5+P7 adapter complex for use in the BEST library build.

This protocol is based on: Carøe, C., Gopalakrishnan, S., Vinner, L., T. Mak, S. S., S. Sinding, M. H., Samaniego, J. A., Wales, N., Sicheritz-Pontén, T., & P. Gilbert, M. T. (2018). Single-tube library preparation for degraded DNA. *Methods in Ecology and Evolution*, 9(2), 410-419. <https://doi.org/10.1111/2041-210X.12871>

Always **cite the original paper** when using this protocol, in addition to citing this protocols.io protocol.

Protocol materials

⊗ Molecular grade H2O Merck MilliporeSigma (Sigma-Aldrich) Catalog #W4502

⊗ 5M NaCl Invitrogen - Thermo Fisher Catalog #AM9760G

⊗ 1M Tris-HCl, pH 8.0 Invitrogen - Thermo Fisher Catalog #15568025

⊗ EB buffer Qiagen Catalog #19086

⊗ 0.5M EDTA Invitrogen Catalog #AM9260G

⊗ TE Buffer, RNase Free Invitrogen Catalog #AM9858



Before start

- UV-treat aliquots (ddH₂O, etc.) for 10 minutes before use.
- Mix all the reagents well, before use.
- Do not UV Eppendorf Tubes/PCR tubes
- Clean all surfaces with 5% bleach solution followed by 70% ethanol before and after use.
- Clean all equipment with 70% ethanol before and after use.

Buffer Preparation

1 Prepare the following reagents:

Adapters (sequences in table below):

- IS1 (stock: 500 μ M)
- IS2 (stock: 500 μ M)
- IS3 (stock: 500 μ M)

Other reagents:

-  Molecular grade H₂O **Merck MilliporeSigma (Sigma-Aldrich) Catalog #W4502**
-  5M NaCl **Invitrogen - Thermo Fisher Catalog #AM9760G**
-  1M Tris-HCl, pH 8.0 **Invitrogen - Thermo Fisher Catalog #15568025** (identical to  EB buffer **Qiagen Catalog #19086**)
-  0.5M EDTA **Invitrogen Catalog #AM9260G**
-  TE Buffer, RNase Free **Invitrogen Catalog #AM9858**

Table 1 Adapter oligo sequences in 5'-3' direction. Sequences sourced from Carøe et al. (2018) Table S1.

Oligo	Sequence
IS1	5'OMedT*ACACTCTTTCCCTACACGACGCTCTTCCGATCT*A
IS2	5'OMedT*GTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT*A
IS3	AGATCGGAAGAGC[C3spacer]

*=Phosphothioate linkage

All oligos can be synthesized at Biomers, Germany, HPLC purified, shipped dry, and resuspended in 10 mM Tris-HCl.

2 Prepare the oligo hybridization buffer (10X)

Reagent	Volume	Final Concentration
5M NaCl	50 μ L	500 mM
1M Tris-HCl, pH 8.0	5 μ L	10 mM
0.5M EDTA	1 μ L	1 mM



	Reagent	Volume	Final Concentration
	ddH ₂ O	444 µL	/
	Total	500 µL	

Preparation of double stranded adapter

3 P5 Adapter (in PCR tube)

	Reagent	Volume
	IS1 (500µM)	40 µL
	IS3 (500µM)	40 µL
	ddH ₂ O	10 µL
	Oligo hybridization buffer	10 µL
	Total	100 µL

4 P7 Adapter (in a different PCR tube to P5)

	Reagent	Volume
	IS2 (500µM)	40 µL
	IS3 (500µM)	40 µL
	ddH ₂ O	10 µL
	Oligo hybridization buffer	10 µL
	Total	100 µL

5 Vortex both tubes & spin down. Make sure they are well mixed.

6 Incubate in a thermal cycler with a heated lid (105°C) at 95°C for 10s, then ramp down to 12°C at a rate of 0.1°C per second (see table below).



	Temperature	Time
	95°C	10 sec
	Ramp down to 12°C	0.1°C per sec
	10°C	Hold

- 7 After cooldown, mix the two adapters in a 1.5 mL Eppendorf tube. This is your final stock tube, so label properly with name (e.g. BEST P5+P7 adapter mix), date, and concentration.
This gives a P5+P7 solution of **200 µM**.
- 8 Dilute with TE buffer to the desired working concentration, e.g., 20 µM, 10 µM, 5 µM, 3 µM, 2 µM.
The working concentration depends on the input DNA amount, which will be specified in the BEST library build protocol.

Protocol references

Carøe, C., Gopalakrishnan, S., Vinner, L., T. Mak, S. S., S. Sinding, M. H., Samaniego, J. A., Wales, N., Sicheritz-Pontén, T., & P. Gilbert, M. T. (2018). Single-tube library preparation for degraded DNA. *Methods in Ecology and Evolution*, 9(2), 410–419. <https://doi.org/10.1111/2041-210X.12871>