

Jun 13, 2019

Bold SW 1 NV modified medium

DOI

<https://dx.doi.org/10.17504/protocols.io.s23eggn>



Roscoff Culture Collection¹

¹Station Biologique, CNRS-Sorbonne Université, Roscoff France

Roscoff Culture Collection



Daniel Vaultot

Station Biologique, Roscoff, France , Nanyang Technologica...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.s23eggn>

Protocol Citation: Roscoff Culture Collection 2019. Bold SW 1 NV modified medium. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.s23eggn>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 28, 2018

Last Modified: June 13, 2019

Protocol Integer ID: 15163

Keywords: eukaryote, modified medium medium, bold sw, medium medium, medium

Abstract

Medium used for some eukaryotes

Adapted from CCAP web site : <https://www.ccap.ac.uk/media/documents/BB.pdf>

See also : <http://www.ife.ac.uk/ccap/MediaRecipes.html#bb>

Materials

MATERIALS

 Red Sea Salt **Red Sea**

Before start

Please refer to our general recommendations to grow cultures :

<https://www.protocols.io/private/A48906DC1374AD6281495CB86A8F092F>



1 Change from standard Bold medium

	Volume for 500 ml	Component	Concentration Stock Solutions		
	5 ml	CaCl ₂	CaCl ₂ ·2H ₂ O	0.5 gr	20 0 ml
	5 ml	NaNO ₃	NaNO ₃	16 gr	20 0 ml
	5 ml	MgSO ₄	MgSO ₄ ·7H ₂ O	10 gr	20 0 ml
	10 ml	PO ₄	From PCRS11 media		
	5 ml	NH ₄ Cl (K)	From K media		
	3 ml	K trace metal	From K media		
	50 µl	Vitamines			
	467 ml	SW sterile			

.