

Apr 21, 2023

Bone Decalcification Protocol Using 14% EDTA Buffer Solution pH 7.2 - 7.4

DOI

dx.doi.org/10.17504/protocols.io.4r3l2718qg1y/v1

Donna Peehl¹, Renuka Sriram¹, Rosalie Nolley¹

¹University of California, San Francisco



Rosalie Nolley

UCSF

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l2718qg1y/v1>

Protocol Citation: Donna Peehl, Renuka Sriram, Rosalie Nolley 2023. Bone Decalcification Protocol Using 14% EDTA Buffer Solution pH 7.2 - 7.4. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.4r3l2718qg1y/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 09, 2023



Last Modified: April 21, 2023

Protocol Integer ID: 78391

Keywords: bone decalcification protocol, decalcification, successful decalcification, calcified tissue, minerals from bone, removing mineral, routine processing to paraffin, mineral, bone, quality paraffin section

Funders Acknowledgements:

NIH

Grant ID: U24CA253377

Abstract

This protocol describes the steps required for successful decalcification. Decalcification describes the technique for removing minerals from bone or other calcified tissue so that good-quality paraffin sections can be prepared that will preserve all the essential microscopic elements. Decalcification is carried out after the specimen has been thoroughly fixed and prior to routine processing to paraffin.









Materials

	A	B	C
	10% Buffered Formalin	Fisher	STL286001
	Ethylenediaminetetraacetic (Acid free EDTA)	VWR	97061-404
	Ammonium hydroxide	Fisher	1336216
	Phosphate Buffered Saline	Gibco	10010-023




Troubleshooting

Preparation of Ethylenediaminetetraacetic acid (EDTA)

- 1 Make 14% EDTA solution fresh
 - a. Add  140 g free acid EDTA to  700 mL distilled H₂O
 - b. On stir plate in the fume hood, add ammonium hydroxide,  30 mL at a time, until solution clears (about  90 mL total)
 - c. Add H₂O to almost  1 L. Check pH and adjust with ammonium hydroxide dropwise up to pH 7.4, then adjust final volume to  1 L

Procedure

3d






- 2 Dissect bone and remove as much soft tissue as possible.
- 3 After appropriate fixation in 10% buffered formalin -  24:00:00 -  48:00:00 hours, wash tissue in distilled H₂O or used EDTA solution. 3d
- 4 Place tissue in 14% EDTA solution at  4 °C C on a stirring device to circulate the EDTA. Use enough solution to saturate tissue (fluid volume to tissue ratio is critical for the decalcification process). The EDTA solution should be at least 20X the volume of the tissue to ensure proper decalcification.

Note

Fluid volume to tissue ratio is critical for the decalcification process.
- 5 Periodically check bone for adequate decalcification. Refresh EDTA solution daily for first 5 days, then may leave in same EDTA solution without changing.
- 6 Decalcification is complete when bone is soft and pliable. This may take 10 days or more depending on tissue size. Can check by X-ray, or probe with a needle and/or bend tissue



to determine if tissue is soft enough to section. Over-decalcification will cause tissue or cells to lose affinity for certain stains.

- 7 Rinse once with Phosphate Buffered Saline (PBS) for  00:30:00 minutes 
- 8 Rinse twice in ddH₂O for  00:30:00 minutes 
- 9 Place in 70% Etanol (EtOH). Store at  4 °C (no more than 2 weeks)
- 10 Proceed with tissue processing