



Aug 04, 2023

Bone Marrow Derived Macrophage (BMDM) differentiation and maintenance

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlkw7bwl5r/v1>

Narayana Yadavalli^{1,2,3,4,5}, Shawn M. Ferguson^{1,2,3,4,6,5}

¹Department of Cell Biology, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

²Neuroscience, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

³Program in Cellular Neuroscience, Neurodegeneration and Repair;

⁴Wu Tsai Institute Yale University School of Medicine, New Haven, Connecticut 06510, USA;

⁵Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815, USA;

⁶Kavli Institute for Neuroscience, Yale University School of Medicine, New Haven, Connecticut 06510, USA



Berrak Ugur

Yale University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlkw7bwl5r/v1>

Protocol Citation: Narayana Yadavalli, Shawn M. Ferguson 2023. Bone Marrow Derived Macrophage (BMDM) differentiation and maintenance. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.j8nlkw7bwl5r/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 24, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 82373

Keywords: Macrophages, L929 media, Bone marrow, ASAPCRN, bone marrow derived macrophage, macrophage, isolation of bone marrow, bmdm, bone marrow

Funders Acknowledgements:

ASAP

Grant ID: 000580

Abstract

This protocol describes the isolation of bone marrow derived macrophages.

Attachments



[724-1815.docx](#)

17KB

Materials

Reagents required

- L929 conditioned media
-  Penicillin/Streptomycin **Thermo Fisher Scientific Catalog #Invitrogen 15140-122**
-  Gibco™ Fetal Bovine Serum **Fisher Scientific Catalog #16-140-071**
-  Gibco™ DMEM/F-12 HEPES **Thermo Fisher Scientific Catalog #11330032**
-  GlutaMAX™ Supplement **Gibco - Thermo Fisher Scientific Catalog #35050061**



1

Note

Note: For mouse bone marrow-derived macrophage (BMDMs) primary cultures, each experiment involved age and sex matched C57BL/6 mice between 3-9 months of age. *Lrrk2*KO (*Lrrk2*^{tm1.1Mjff}) and G2019S (*Lrrk2*^{tm1.1Hlme}) homozygous knockin mice were obtained from The Jackson Laboratory⁷⁶.

Euthanize mice via CO₂ or isoflurane inhalation and cervical dislocation.

2

Collect femurs and flush cavities with  5 mL ice-cold PBS.

3

Then, collect bone marrow cells by centrifugation  0.3 rcf, 00:05:00 .

5m



4

Resuspend the resulting pellet and differentiate for 6 days in culture media containing: DMEMF12 supplemented with 20% FBS, 20% L929 conditioned media, 1% penicillin-streptomycin and 1% GlutaMAXTM.