


Nov 16, 2018

Version 3

Bounce PCR V.3

DOI

dx.doi.org/10.17504/protocols.io.vhge33w



Sam Mugford¹, Saskia Hogenhout¹

¹John Innes Centre

OpenPlant Project



Sam Mugford

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.vhge33w>

Protocol Citation: Sam Mugford, Saskia Hogenhout 2018. Bounce PCR. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.vhge33w>



Manuscript citation:

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 09, 2018

Last Modified: November 16, 2018

Protocol Integer ID: 17672

Keywords: PCR, Cloning, Synthetic Biology, optimisation free pcr for synthetic biology cloning application, free pcr method suitable for the rapid cloning, bounce pcr bounce pcr, bounce pcr, version of touchdown pcr, touchdown pcr, free pcr method, gene amplification, rapid cloning, molecular cloning, synthetic biology cloning application, cloning of multiple dna component, spurious priming during gene amplification, pcr, expensive alternative to traditional cloning, gene synthesis, advantage of primer, free pcr, multiple dna component, traditional cloning, biomolecular technique, synthetic biology project, cloning, most dna fragments without time, nucleic acids research, most dna fragment, increased primer, gateway recombination site, primer, target amplicon, amplification efficiency in later cycle, journal of biomolecular technique, gibson cloning, dna, early amplification cycle, amplification efficiency, vector overlaps for gibson cloning

Abstract

Bounce PCR: optimisation free PCR for synthetic biology cloning applications.

Synthetic biology projects require the cloning of multiple DNA components, and increasingly this is done through in-vitro gene synthesis. This remains an expensive alternative to traditional cloning by PCR. However, PCR often requires time-consuming optimisation. Bounce PCR is a largely optimisation-free PCR method suitable for the rapid cloning of multiple targets. Bounce PCR is a modified version of Touchdown PCR that takes advantage of primer-extension sequences -commonly used in molecular cloning – to successfully amplify most DNA fragments without time consuming optimisation. Touchdown (and related methods) employ a sequential lowering of the annealing temperature over successive cycles to minimise non-specific primer-binding during the early amplification cycles, but maximise amplification efficiency in later cycles when the target amplicon is more abundant (Don et al. 1991; Rowther et al. 2012). However, the lower annealing temperature in later cycles may still allow for non-specific amplification. Bounce PCR takes advantage of the increased primer melting temperature caused by the addition to the template of primer-extension sequences (for example: restriction enzyme recognition sites, Gateway recombination sites, or vector overlaps for Gibson cloning) to maintain the specificity in later cycles- without sacrificing efficiency- by increasing the annealing temperature again after an initial round of touchdown.

Don, R. H., Cox, P. T., Wainwright, B. J., Baker, K., & Mattick, J. S. (1991). 'Touchdown' PCR to circumvent spurious priming during gene amplification. *Nucleic acids research*, 19(14), 4008.

Rowther, F. B., Kardooni, H., & Warr, T. (2012). TOUCH-UP gradient amplification method. *Journal of biomolecular techniques : JBT*, 23(1), 1-3.

Guidelines

This protocol describes only the thermal cycling parameters of the PCR, and should work with your preferred DNA polymerase and protocol for setting up the reaction mixture. The advantage of Bounce PCR relies on the presence of primer-extension sequences, if your primers are an exact match to your template, without 5' extension sequences, then it is unlikely to offer an advantage over standard Touchdown PCR.


Materials

PCR mixture in suitable reaction tube.

Thermal cycler.

Troubleshooting

Initial denaturation

1  94 °C for  00:02:00

Denaturation

2  94 °C for  00:00:30

Annealing

3

	Temperat ure	Duration
	60°C	40 seconds
	59°C	40 seconds
	58°C	40 seconds
	57°C	40 seconds
	56°C	40 seconds
	55°C	40 seconds
	54°C	40 seconds
	53°C	40 seconds
	52°C	40 seconds
	51°C	40 seconds
	50°C	40 seconds





49°C	40 seconds
48°C	40 seconds
47°C	40 seconds
46°C	40 seconds
45°C	40 seconds
45.6°C	40 seconds
46.2°C	40 seconds
46.8°C	40 seconds
47.4°C	40 seconds
48°C	40 seconds
48.6°C	40 seconds
49.2°C	40 seconds
49.8°C	40 seconds
50.4°C	40 seconds
51°C	40 seconds
51.6°C	40 seconds
52.2°C	40 seconds
52.8°C	40 seconds
53.4°C	40 seconds
54°C	40 seconds



54.6°C	40 seconds
55.2°C	40 seconds
55.8°C	40 seconds
56.4°C	40 seconds
57°C	40 seconds
57.6°C	40 seconds
58.2°C	40 seconds
58.8°C	40 seconds
59.4°C	40 seconds
60°C	40 seconds

60 degree C for 40 seconds (decreasing by 1 degree C per cycle for the first 15 cycles ("Touchdown"), then increasing by 0.6 degree C per cycle for the next 25 cycles ("Bounce"))

Extension


4  72 °C for  00:01:00 per kb of target length

Repeating

5 Repeat steps 2-4 for 40 cycles, varying the temperature of step 3 as described

 go to step #2 undefined

Final extension

6  72 °C for  00:06:00

