Feb 12, 2019 Version 1

# O Bovine satellite cell Pax7 ICC V.1

Forked from Bovine satellite cell Pax7 ICC

DOI

dx.doi.org/10.17504/protocols.io.xzyfp7w


Andrew Stout<sup>1</sup>, Natalie Rubio<sup>1</sup>

<sup>1</sup>Tufts University



Andrew Stout

Tufts University



# 



#### DOI: <u>dx.doi.org/10.17504/protocols.io.xzyfp7w</u>

Protocol Citation: Andrew Stout, Natalie Rubio 2019. Bovine satellite cell Pax7 ICC. protocols.io <u>https://dx.doi.org/10.17504/protocols.io.xzyfp7w</u>

License: This is an open access protocol distributed under the terms of the <u>Creative Commons Attribution License</u>, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's working

Created: February 12, 2019

Last Modified: February 12, 2019

Protocol Integer ID: 20248

Keywords: Staining, Immunocytochemistry

#### Abstract

Staining primary bovine satellite cells for Pax7, a common marker of satellite cells and myogenic potential. Protocol developed for <u>https://www.thermofisher.com/antibody/product/PAX7-Antibody-Polyclonal/PA5-68506</u> (Thermo Fisher Pax7 antibody PA5-68506; rabbit IgG anti-Pax7)

# Guidelines

For reference general volumes for given well formats are:

- 96-well = 100 uL
- 48-well = 150 uL
- 24-well = 300 uL
- 12-well = 500 uL
- 6-well = 750 uL

\*\* recommended to use a PAP-pen to select a smaller region of 6-well plates after initial fixing / washing, as this will save antibody

# Materials

#### MATERIALS

X 4% paraformaldehyde/1XPBS solution

Soat-anti-rabbit-Alexafluor 488 Thermo Fisher Scientific Catalog #A11008

🔀 PBS

X VECTASHIELD<sup>®</sup> Hardset<sup>™</sup> Antifade Mounting Medium Catalog #H-1400

X PAX7 Polyclonal Antibody Thermo Fisher Scientific Catalog #PA5-68506

🔀 Wash buffer (PBS / 5% goat serum / 0.05% NaAzide)

X Permeabilization solution (PBS / 0.5% Triton X-100)

**X** PBST (PBS 1:1000 Tween-20)

X Phalloidin 594 Thermo Fisher Scientific Catalog #A12381

## Fixation and Permeabilization (1 hour)

- 1 Aspirate media from cells
- 2 add cold 4% PFA to cells (enough to cover cells or scaffolds)

3 Incubate at room temperature for 30 minutes 🚫 00:30:00

- 4 Wash 3x with room temperature PBS
  - NOTE: at this point, can parafilm and leave in the fridge overnight (or up to 1 week) before staining
- 5 Aspirate PBS and add cold Permeabilization solution for 15 minutes 👀 00:15:00
- 6 Wash 3x with cold PBST

### Primary Stain (1 hour, overnight incubation)

7 Aspirate PBST and add cold Wash buffer for 45 minutes 🕚 00:45:00

Note

During soak, can move to step 8

- 8 Dilute primary antibodies in wash buffer and keep on ice (protected from light). For given antibody, use the following dilutions:
  - anti-Pax7 (1:500)
  - Phalloidin-594 (1:100)

note\* prepare enough antibody solution for all conditions (a little extra is usually good to make sure there is enough)

- 9 After step 7 incubation, wash 3x with cold PBST
- 10 Add primary antibody solutions and incubate overnight at 4C (parafilm to avoid evaporation)

#### Secondary Stain (1.5 hours)

- 11 Wash 3x with cold PBST
- Aspirate PBST and add cold Wash buffer for 15 minutes 🕥 00:15:00

Note	
during soak, can move to step 13	

- 13 Dilute secondary antibodies in wash buffer and keep on ice (protected from light). For given antibody, use the following dilutions:
  - 488 goat-anti-rabbit (1:500)

note\* prepare enough antibody solution for all conditions (a little extra is usually good to make sure there is enough)

STEP CASE

for 3D, when not using DAPI mounting media 3 steps

In the case where you're not planning to use a dapi mounting media (ie 3D constructs), prepare a DAPI solution in a suitable blocking buffer (ie Wash Buffer or a BSA-containing buffer), and use that to prepare antibody solutions, instead of plain wash buffer

- 14 After step 12 incubation, aspirate Wash buffer from cells, and add secondary antibody solutions. Incubate in the dark at room temperate for 60 minutes 01:00:00
- 15 Wash cells 3X with cold PBST, leaving the last wash to soak for 5 minutes 👀 00:05:00

- 16 Aspirate PBST, and add DAPI mounting media. Cover with cover-slip, and image after 10 minutes 00:10:00
  - Pax7 = green
  - Actin cytoskeleton = red
  - nuclei = blue