



Oct 14, 2019

## Bradford protein assay in 96-Well plate

DOI

<https://dx.doi.org/10.17504/protocols.io.76xhrfn>

iGEM Dusseldorf<sup>1</sup>

<sup>1</sup>Heinrich-Heine Universität Düsseldorf



Igem Dusseldorf

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.76xhrfn>

**Protocol Citation:** iGEM Dusseldorf 2019. Bradford protein assay in 96-Well plate. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.76xhrfn>

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 14, 2019



**Last Modified:** October 14, 2019

**Protocol Integer ID:** 28599

**Keywords:** well plate assay for quantification,  $\mu\text{l}$  of the diluted bradford, bradford protein, measurements protein concentration measurement, bsa stock solution to  $200\mu\text{g}$ , diluted bradford, well plate assay, bradford reagent, bovine serum albumin,  $\mu\text{l}$  of solution buffer, solution buffer needed 002002020180404016060601408080120100100100, reagent to each well, serum albumine standard, solution buffer, comparing measured absorbance, dilute the roti, amount of protein, bsa stock solution, measured absorbance, pipette  $50\mu\text{l}$  of the standard, expected protein level, pipette  $50\mu\text{l}$ , od595 for each well, dilution, ml roti, pipette  $50\mu\text{l}$ , concentration, ml of water,  $200\mu\text{g}$ , ml in the same buffer, standard concentration

## Abstract

Assay for quantification of protein by comparing measured Absorbance at 595 nm to bovine serum albumine standard

- Bovine serum albumin (concentrated >200 µg/ml)
- Roti-Quant 5x (ROTH)
- 96 Well Plate
- Plate Reader

### Standard measurements

Dilute your BSA stock solution to 200µg/ml in the same buffer used for your solution of interest.

Dilute the Roti-Quant 5x reagent at a rate of 2:7,5. For 15 ml solution, this means 4 ml Roti-Quant added to 11 ml of Water.

Create 200 µl each of standard concentrations according to the following table:

|  | Con<br>cent<br>ratio<br>n<br>[µg/<br>ml] | Volu<br>me<br>of<br>stoc<br>k<br>nee<br>ded | Volu<br>me<br>of<br>solut<br>ion<br>buff<br>er<br>nee<br>ded |
|--|------------------------------------------|---------------------------------------------|--------------------------------------------------------------|
|  | 0                                        | 0                                           | 200                                                          |
|  | 20                                       | 20                                          | 180                                                          |
|  | 40                                       | 40                                          | 160                                                          |
|  | 60                                       | 60                                          | 140                                                          |
|  | 80                                       | 80                                          | 120                                                          |
|  | 100                                      | 100                                         | 100                                                          |

Work in Triplicates. Pipette 50µl of the standards into the 96-Well plate

Add 200 µl of the diluted Bradford reagent to each well.

Incubate for 5 minutes at room temperature

Measure the OD<sub>595</sub> for each well. Create a regression line from the averages of each concentrations' measurements

### Protein concentration measurements

Dilute your solution of interest 1:20 or 1:40, depending on expected protein levels.

Work in triplicates and use blanks. Pipette 50 µl of each dilution into a well on the 96-Well plate.



When working with solutions containing pigments, prepare additional wells to which NO Bradford reagent will be added

Add 200  $\mu$ l of the diluted Bradford reagent to each well. For the pigment controls, add 200  $\mu$ l of solution buffer instead. Optimally, use a multi-pipette to reduce the time needed for adding the reagent.

Incubate for 5 minutes at room temperature

Measure the OD<sub>595</sub> for each well.

When processing the data, correct the measurements for the blanks. When pigments are present in the solution, correct the OD<sub>595</sub> by subtracting the OD<sub>595</sub> of the pigment controls from the appropriate measurements.

Determine the amount of protein by putting the achieved OD<sub>595</sub> of your solution of interest into the regression formula and multiplying by the dilution factor.

