

Aug 07, 2024

## Brain Clarification for Neuromelanin visualisation by OPT

 In 2 collections

DOI

[dx.doi.org/10.17504/protocols.io.n92ld8qdov5b/v1](https://dx.doi.org/10.17504/protocols.io.n92ld8qdov5b/v1)

Ariadna Laguna<sup>1</sup>, Miquel Vila<sup>1</sup>

<sup>1</sup>VHIR-CIBERNED-ASAP

Vilalab Public

Nuria



Miquel Vila

VHIR-CIBERNED-ASAP

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ld8qdov5b/v1>

**Protocol Citation:** Ariadna Laguna, Miquel Vila 2024. Brain Clarification for Neuromelanin visualisation by OPT. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n92ld8qdov5b/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 06, 2024

**Last Modified:** August 13, 2024

**Protocol Integer ID:** 104759

**Keywords:** neuromelanin visualisation, whole mouse brain

## Disclaimer

The **protocols.io** team notes that research involving animals and humans must be conducted according to internationally-accepted standards and should always have prior approval from an Institutional Ethics Committee or Board.

## Abstract

Protocol for the clarification of whole mouse brains to allow neuromelanin visualisation by Optical Projection Tomography

## Troubleshooting

## Brain processing

- 1 Mice are transcardially perfused in 4% PFA and postfixed in the same fixative overnight
- 2 After rinsing several times in PBS, brains are embedded in 1% low melting point agarose in water
- 3 Brains are dehydrated in 100% methanol covered in aluminum foil to protect the specimen from light; change the methanol three times during the next 24 h or until the specimen is completely dehydrated
- 4 Once dehydrated the brains are incubated during 4 hours shaking in 66% dichloromethane (DCM) /33% methanol at room temperature
- 5 The brains are incubated with shaking in 100% DCM 30 minutes twice to complete the delipidation process
- 6 The brains are chemically cleared with BABB (a combination of 1 part of benzyl alcohol and 2 parts of benzyl benzoate)
- 7 The brains are subsequently immersed in a BABB-filled chamber for Optical projection tomography (OPT) imaging

## Brain imaging

- 8 The brains are scanned with transmission light with three different filters (a cyan fluorescent protein -CFP-filter: emission 460-500nm, a green fluorescent protein -GFP-filter: emission 500-550nm and a DSR filter: emission >610nm) using an OPT imaging system mounted on a Leica MZ 16 FA microscope
- 9 The brain images are visualised using the open source software Fiji