

Apr 05, 2024

Brain Homogenization and MSD Protocol for Mouse Brain and Serum

DOI

<https://dx.doi.org/10.17504/protocols.io.kxygxprqkl8j/v1>

Scott Vermilyea¹

¹University of Minnesota

Metabolomics Protocols ...

Team Lee



Jane Balster

ASAP - Team Lee

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygxprqkl8j/v1>

Protocol Citation: Scott Vermilyea 2024. Brain Homogenization and MSD Protocol for Mouse Brain and Serum. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.kxygxprqkl8j/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 12, 2021

Last Modified: May 31, 2024

Protocol Integer ID: 55080

Keywords: Homogenization, MSD Protocol, Mouse Brain and Serum, ASAPCRN, msd protocol of the mouse brain, msd protocol for mouse brain, brain homogenization, mouse brain, msd protocol, serum this protocol detail, homogenization, mouse

Disclaimer

The protocols.io team notes that research involving animals and humans must be conducted according to internationally-accepted standards and should always have prior approval from an Institutional Ethics Committee or Board.

Abstract

This protocol details about the homogenization and MSD protocol of the mouse brain using various panel kits.

Attachments



[320-681.pdf](#)

54KB

Materials

Reagent, Kits and Equipment List

1.

Equipment	
gentleMACS™ C Tubes	NAME
C Tubes	TYPE
gentleMACS™	BRAND
130-093-237	SKU
https://www.miltenyibiotec.com/US-en/products/gentlemacs-c-tubes.html#gref	LINK

2.  RIPA Lysis and Extraction Buffer **Thermo Fisher Catalog #89901**

3. Halt PIC (ThermoFisher 1861279)

4.

Equipment	
Direct Detect® Assay-free Cards	NAME
Assay-free Cards	TYPE
Direct Detect®	BRAND
DDAC00010-GR	SKU
https://www.emdmillipore.com/US/en/product/Direct-Detect-Assay-free-Cards,MM_NF-DDAC00010-GR?bd=1	LINK
Membrane cards used for protein quantitation in the Direct Detect® Infrared Spectrometer	SPECIFICATIONS

5. HCl (Millipore Sigma H1758-500ML)

- NaOH (Millipore Sigma 415413-500ML)
-  HEPES Merck MilliporeSigma (Sigma-Aldrich) Catalog #83264-500ML-F
-  Proinflammatory Panel 1 mouse MESO SCALE DIAGNOSTICS, LLC. Catalog #MSD K15048D
-  Cytokine Panel 1 Mouse Kit MESO SCALE DIAGNOSTICS, LLC. Catalog #K15245D
-  U-PLEX TGF- β Combo mouse MESO SCALE DIAGNOSTICS, LLC. Catalog #MSD K15242K-2
-

Equipment	
Octo Dissociator with Heaters	NAME
Octo Dissociator	TYPE
gentleMACS™	BRAND
130-096-427	SKU
https://www.miltenyibiotec.com/US-en/products/gentlemacs-octo-dissociator-with-heaters.html	LINK

12.

Equipment	
Direct Detect Spectrometer	NAME
Spectrometer	TYPE
Direct Detect®	BRAND
DDHW00010-WW	SKU
https://www.emdmillipore.com/US/en/product/Direct-Detect-Spectrometer,MM_NF-DDHW00010-WW?bd=1	LINK

13. Meso Sector S 600 (MSD 1201)



Brain Homogenization Using Miltenyi gentleMACS Octodissociator

43m

- 1 Weight out brain using pre-chilled C-tube.
- 2 Use  2 mL RIPA (1x Halt PIC) per whole brain,  1 mL RIPA (1x Halt PIC) per half brain.
- 3 Homogenize brain using "protein_01_01" routine on Miltenyi gentle MACS Octodissociator.
- 3.1 Repeat homogenization routine if not fully homogenized.
- 4 Vortex briefly and incubate  00:30:00  On ice . 
- 5 Centrifuge at  3500 rpm, 4°C, 00:03:00 . 
- 6 Transfer homogenate to two  1.5 mL eppendorf tubes. 
- 7 Vortex briefly and centrifuge at  14000 x g, 4°C, 00:10:00 . 
- 8 Transfer supernatant to new  1.5 mL eppendorf tubes 
- 9 Perform 1:10 serial dilution. 
- 9.1 Dilute using  10 μ L sample to  90 μ L RIPA (1x Halt PIC) to achieve 1:100 dilution. (1/2)



- 9.2 Dilute using  10 μL sample to  90 μL RIPA (1x Halt PIC) to achieve 1:100 dilution. (2/2)
- 10 Blot  2 μL diluted sample onto Direct Detect assay card (run duplicates).
- 11 Use  2 μL RIPA (1x Halt PIC) as blank.
- 12 Read assay card on Direct Detect Spectrometer.
- 13 Detectable range is  0.2 mg/mL to  5.0 mg/mL , dilute or concentrate as needed.
- 14 Average all readings for each sample to determine concentration.

Brain MSD Protocol Proinflammatory Panel 1 mouse and Cytokine Panel 1 mouse

- 15 Using concentration determined, calculate volume needed for  200 μg protein.
- 16 Use RIPA(1x PIC) to dilute  200 μg protein to  100 μL ( 2 $\mu\text{g}/\mu\text{L}$). 
- 17 Combine  100 μL sample and  100 μL MSD diluent in U-bottom plate. 
- 18 Load  50 μL diluted sample to MSD plate per duplicate. 
- 19 Proceed with MSD assay protocol from manufacturer.

Brain MSD Protocol U-PLEX TGF- β Combo mouse

10m



20 Using concentration determined, calculate volume needed for $\text{200 } \mu\text{g}$ protein.

21 Use RIPA(1x PIC) to dilute $\text{200 } \mu\text{g}$ protein to $\text{100 } \mu\text{L}$ ($\text{2 } \mu\text{g}/\mu\text{L}$).

22 Load $\text{100 } \mu\text{L}$ sample into U-bottom plate.

23 Add $\text{20 } \mu\text{L}$ 1M 1 Mass Percent HCl into U-bottom plate and shake for 00:10:00 at $\text{25 } ^\circ\text{C}$.

24 Neutralize sample with $\text{14 } \mu\text{L}$ 1M 1.2 Mass Percent NaOH in 0.5 Mass Percent HEPES.

25 Combine $\text{100 } \mu\text{L}$ treated sample and $\text{100 } \mu\text{L}$ MSD diluent in new U-bottom plate.

26 Load $\text{50 } \mu\text{L}$ diluted sample to MSD plate per duplicate.

27 Proceed with MSD assay protocol from manufacturer.

Serum MSD Protocol Proinflammatory Panel 1 mouse

28 Add $\text{25 } \mu\text{L}$ MSD Diluent to MSD Plate per duplicate.

29 Add $\text{25 } \mu\text{L}$ serum sample to MSD Plate per duplicate.

30 Proceed with MSD assay protocol from manufacturer.

Serum MSD Protocol Cytokine Panel 1 mouse



31 Add  37.5 μL MSD Diluent to MSD Plate per duplicate.



32 Add  12.5 μL serum sample to MSD Plate per duplicate.



33 Proceed with MSD assay protocol from manufacturer.

Serum MSD Protocol U-PLEX TGF- β Combo mouse

10m

34 Add  50 μL serum sample to U-bottom plate.



35 Add  10 μL [IM] 1 Mass Percent HCl into U-bottom plate and shake for  00:10:00 at  25 $^{\circ}\text{C}$.

10m



36 Neutralize sample with  7 μL [IM] 1.2 Mass Percent NaOH in [IM] 0.5 Mass Percent HEPES.



37 Add  25 μL MSD Diluent to MSD Plate per duplicate.



38 Add  25 μL treated serum sample to MSD Plate per duplicate.



39 Proceed with MSD assay protocol from manufacturer.