



May 02, 2024

🌐 BrainSaw 50 mM pH 7.4 PB slicing buffer

DOI

dx.doi.org/10.17504/protocols.io.q26g714r3gwz/v1

Rob Campbell¹

¹Sainsbury Wellcome Centre, UCL



Rob Campbell

Sainsbury Wellcome Centre, UCL

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g714r3gwz/v1>

Protocol Citation: Rob Campbell 2024. BrainSaw 50 mM pH 7.4 PB slicing buffer. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.q26g714r3gwz/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 02, 2024

Last Modified: May 02, 2024

Protocol Integer ID: 99107

Keywords: brainsaw 50 mm ph, phosphate buffer, pb slicing buffer, pb slicing buffer this protocol, 2l of ph, brainsaw, mm ph, ph, photon imaging, imaging

Abstract

This protocol makes 2L of pH 7.4 phosphate buffer. This solution is used for serial section 2-photon imaging.

Materials

2 L GL45 bottle

Sodium phosphate monobasic monohydrate **S23120-1000.0 1 kg(CAS 10049-21-5)**

Sodium phosphate dibasic heptahydrate **5 kg**

magnetic stirrer

weighing scales

Protocol materials

⊗ Sodium phosphate monobasic monohydrate **Catalog #S23120-1000.0**

⊗ Sodium phosphate dibasic heptahydrate **Melford Catalog #S23175-5000.0**

⊗ Sodium phosphate monobasic monohydrate **Catalog #S23120-1000.0**









⊗ Sodium phosphate dibasic heptahydrate **Melford Catalog #S23175-5000.0**

Troubleshooting

Safety warnings

- ⚠ Ensure the first solute is completely dissolved before adding the second. Not doing this will result in a small quantity of fine precipitate being formed.



- 1 Fill 2L flask with MQ water to the 2L notch. Place on magnetic stirrer and set to 1000 RPM.
- 2 Weigh  3.1 g of  Sodium phosphate monobasic monohydrate **Catalog #S23120-1000.0** and pour into the flask.
- 3 Weigh out  20.8 g of  Sodium phosphate dibasic heptahydrate **Melford Catalog #S23175-5000.0** but do not pour into flask until existing solute is dissolved. 
- 4 Wait until  Sodium phosphate monobasic monohydrate **Catalog #S23120-1000.0** has dissolved. Skipping this step will result in a small quantity of fine precipitate being formed.  5m
- 5 Pour  Sodium phosphate dibasic heptahydrate **Melford Catalog #S23175-5000.0** weighed out above into flask. Screw on lid and shake once to ensure no solute is stuck around the rim of the flask. Return to stirring plate at 1000 RPM. 10m
- 6 Use PB once all solute is dissolved. Do not refrigerate: cold PB will form bubbles as it warms and this interferes with imaging.