

Oct 17, 2019

Version 1

Breast tumours dissociation V.1

DOI

dx.doi.org/10.17504/protocols.io.7m9hk96



Samah El Ghamrasni¹

¹Princess Margaret Cancer Research Centre

Human Cell Atlas Metho...

CZI START Project



Samah El Ghamrasni

Princess Margaret Cancer Centre Research

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.7m9hk96>

Protocol Citation: Samah El Ghamrasni 2019. Breast tumours dissociation . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.7m9hk96>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 25, 2019

Last Modified: October 17, 2019

Protocol Integer ID: 28065

Keywords: Breast, single-cell RNAseq, tissue dissociation , cell rnaseq, cryopreserved tissue, fresh breast tissue, dissociation

Abstract

A protocol designed to dissociate fresh breast tissues (surgical specimens and biopsies) for single-cell RNAseq.

The protocol has been demonstrated to work successfully with fresh and cryopreserved tissues.

Guidelines

Requires access to a flow sorter

Materials


MATERIALS

 BSA


 TrypLE

 Gibco Penicillin-Streptomycin (10,000 U/mL) (Pen/Strep) **Fisher Scientific Catalog # 15-140-122**

 FBS **Invitrogen - Thermo Fisher**

 Liberase TL **Roche Catalog #05 401 020 001**

 DMEM **Gibco - Thermo Fisher Scientific Catalog #11885**

 PBS

 SYTOX[®]; Blue Dead Cell Stain, for flow cytometry **Thermo Fisher Catalog #S34857**

Base Media: DMEM + Penstrep + 10%FBS

Resuspension buffer: PBS+0.01%BSA

Troubleshooting



Tissue Dissociation

- 1 Transfer the tissue onto a 10 cm petri dish
- 2 Rinse 1x briefly with ice cold PBS and aspirate PBS off.
- 3 Use a blade to carefully cut the sample into small pieces, approximately 3-4 mm in diameter.
- 4 Transfer pieces into 50ml tube and Resuspend in 5ml of Base media + Liberase (200ug/ml)
- 5 Incubate 2 hours at 37C
- 6 Mix with a 5 ml serological pipet 5 times to break up the pieces.
- 7 Let the pellet settle at the bottom of the tube, and transfer supernatant to a new falcon tube (Filter supernatant using 40um mesh)
- 8 Resuspend remaining tissue in 2mL of TrypLE, incubate for 10 min at 37C
- 9 Repeat Step 6 and 7
- 10 Spin cells down (300g for 10min at 4C) and Resuspend cells in resuspension media
- 11 Count cells

Flow Sorting

- 12 Stain cells with SytoxBlue (1ul/ml)



- 13 Incubate at room temperature for 15 min
- 14 Transfer to ice and proceed to flow sorting