

Mar 21, 2020

Brooks Lab Western Blotting Protocol

Genome Biology

DOI

[dx.doi.org/10.17504/protocols.io.bcsmiwc6](https://doi.org/10.17504/protocols.io.bcsmiwc6)


Brooks Lab University of California¹, Eva Robinson¹, Alison Tang¹

¹University of California, Santa Cruz

BrooksLabUCSC



Alison Tang



OPEN ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.bcsmiwc6](https://doi.org/10.17504/protocols.io.bcsmiwc6)

External link: <https://doi.org/10.1186/s13059-024-03301-y>

Protocol Citation: Brooks Lab University of California, Eva Robinson, Alison Tang 2020. Brooks Lab Western Blotting Protocol. [protocols.io](https://doi.org/10.17504/protocols.io.bcsmiwc6) <https://dx.doi.org/10.17504/protocols.io.bcsmiwc6>

Manuscript citation:

Tang AD, Felton C, Hrabeta-Robinson E, Volden R, Vollmers C, Brooks AN (2024) Detecting haplotype-specific transcript variation in long reads with FLAIR2. Genome Biology 25(). doi: [10.1186/s13059-024-03301-y](https://doi.org/10.1186/s13059-024-03301-y)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 21, 2020

Last Modified: March 21, 2020

Protocol Integer ID: 33325

Keywords: western blot, western blotting, protein expression, proteins,

Abstract

This is a Western Blotting Protocol for Brooks Lab, Department of Biomolecular Engineering, University of California, Santa Cruz.

Attachments



Brooks_Lab_Western_P..



109KB



Materials

MATERIALS

- ✕ Pierce BCA Protein Assay Kit **Thermo Fisher Scientific Catalog #23225**
- ✕ 4–15% Mini-PROTEAN® TGX™ Precast Protein Gels, 10-well, 30 µl **Bio-Rad Laboratories Catalog #4561083**
- ✕ Pierce RIPA Lysis and Extraction Buffer **Thermo Fisher Scientific Catalog #P189900**
- ✕ cOmplete™ Mini Protease Inhibitor Cocktail (Roche) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #04693124001**
- ✕ 4x Laemmli Sample Buffer **Bio-Rad Laboratories Catalog #1610747**
- ✕ 2-Mercaptoethanol **Bio-Rad Laboratories Catalog #1610710**
- ✕ Precision Plus Protein™ Dual Color Standards 500 µl **Bio-Rad Laboratories Catalog #1610374**
- ✕ Precision Plus Protein™ Dual Color Standards 500 µl **Bio-Rad Laboratories Catalog #1610374**
- ✕ Western Blot Box (black size 8.9 cm × 6.5 cm × 2.5 cm) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #Z742099-5EA**
- ✕ β-Actin Antibody (C4) **Santa Cruz Biotechnology Catalog #sc-47778**
- ✕ SuperSignal™ West Pico PLUS Chemiluminescent Substrate **Thermo Fisher Scientific Catalog #34577**

Note: You can use Pierce BCA Protein Assay Kit #23225 or #23227

Additional Reagents required:

- 100 mM Tris-Cl, pH 8.0
- 300 mM NaCl
- 10% NP-40 in ddH₂O
- 10% Na-deoxycholate in ddH₂O (light sensitive!)
- 10% SDS
- 10x PBS
- Tween-20
- 190 proof ethanol
- 10x TGX running buffer
- Safeway brand non-fat milk pouch
- PBST

- BioRad Transblot kit


- LoBind microcentrifuge tubes Protein, 1.5 ml

Additional Equipment





- Sonicator or Bioruptor
- Microplate reader
- 200 µl pipettor
- Varioskan
- BioRad Gel Doc EZ Gel Documentation System
- BioRad ChemiDoc XRS+

Safety warnings



 Please refer to the Safety Data Sheets (SDS) for safety and environmental hazards.



RIPA Lysis Buffer

- 1 Make  100 mL RIPA w/o protease inhibitor , cover, and freeze in *10 ml aliquots* at  -20 °C .

OR

Use pre-made Pierce RIPA buffer, 100 ml, which is aliquoted in  10 mL aliquots and stored in  -20 °C

- 2 Before use, thaw and add **1 tablet** of protease inhibitor (PI) **per 10 ml aliquot**. (Roche cOmplete™, Mini Protease Inhibitor Cocktail).



Harvesting Cells and Preparing Lysate



3

Note



Always keep everything on ice, unless otherwise indicated.

Wash confluent 10 cm plate of cells **2X** in *ice cold* PBS.


- 4 Add  1 mL cold RIPA with PI to cells, and scrape the cells to remove them from the dish. Transfer to a pre-chilled 2 ml tube  On ice .

- 5 Incubate  On ice  00:00:20 , with periodic vortexing.



- 6 Pellet the insoluble material by spinning at max speed in refrigerated ( 4 °C) microcentrifuge  00:00:10 .








- 7 Transfer supernatant to a clean LoBind Protein 1.5 ml tube as lysate. To avoid multiple freeze-thaw cycles, make aliquots (generally  200 µL each).

Note

Additional Sonicator or BioEruptor sonication step often needed here for more complete nucleic acid removal. See below.

- 8 Store lysates aliquots at  -80 °C .

Sonication of Lysate to break up nucleic acids

- 9 Place samples in a 1.5 ml LoBind Protein microfuge tube and prepare a volume balance tube.
- 10 Bring a timer, ear cuffs,  On ice samples, and balance to second floor to Kamakaka Lab to get the plastic tube adaptor (in back of second drawer in 2nd bay from the back of the lab.)
- 11 Sonicator is in cold room down the hall. Put the outflow tube up on the shelf and fill glass reservoir with water to top of white line.
- 12 Sonicator settings: timer on hold, use max setting (10), constant, then flip on.
- 13 Hold tubes in the adaptor in the sonicating water bath for  00:00:30 .
- 14 Place tubes  On ice to cool for  00:01:00 .
- 15 Repeat the 30' sonication (= total of two rounds on max at  00:00:30 .)

BCA Assay for Protein Concentration Determination

- 16 **Pierce BCA Protein Assay Kit** #23225 or #23227, Thermo Sci with 2 mg/ml BSA standard.
- Use instructions from kit.



Equilibrate reagents, samples, and standards to Room temperature . Use the *Microplate reader instructions*.

- 17 Prepare a dilution series of BSA in working range of – from 1 glass vial of stock 2 mg/ml BSA in the kit for the standard curve:

Vial	RIPA buffer, µl	2 mg/ml BSA or standard, µl	Final concentration, µg/ml
A	0	300 of stock	2000
B	125	375 of stock	1500
C	325	325 of stock	1000
D	175	175 of B	750
E	325	325 of C	500
F	325	325 of E	250
G	325	325 of F	125
H	400	100 of G	25
I	400	0	0

- 18 Prepare fresh working reagent 50:1 of A:B, enough for all standards, unknowns, and replicates of the unknowns (should have n=3).
- 19 Add 25 µL of unknown or standard to 96-well plate well followed by 200 µL working reagent using 200 µl pipettor and *gently mix avoiding spillover*.
- 20 Cover plate with foil and incubate 37 °C for 00:30:00 in shaking incubator, gently shaking.
- 21 Cool to Room temperature .





- 22 Measure Absorbance on the VarioSkan at 562 nm.
- 23 Determine protein concentration of your unknown from the standard curve (the Skanit software will plot the curve and the unknowns on it if you edit the standards in the plate layout to include the concentrations above in the table.)
- 24 Export report as an excel doc.



SDS PAGE

- 25 Prepare the following solutions for SDS-PAGE, Transfer, and Antibody Incubations:

Runni ng Buffe r (Tris- Glyci ne)	
100 ml	10x TGX runni ng buffer
900 ml	dd water

Trans fer Buffe r	
200m l	5x Biora d Trans fer Buffe r
200m l	190 proof ethan ol
600m l	dd water








PBST	
------	--



50 ml	10x PBS
0.5 ml	Tween-20
bring vol to 500 ml with dd water	

5 % Milk Block (prep are immediately before use!)	
5 g	Safe way brand non-fat milk pouch
100 ml	PBST

SDS PAGE: *Denaturing samples in Laemmli reducing buffer*

- 26 In a fume hood, Add  100 μ L 2-mercaptoethanol to  900 μ L 4x Laemmli sample buffer to make MLB.
- 27 Dilute samples with MLB at a ratio of *3 parts sample to 1 part MLB*. So for  30 μ L final volume, add  7.5 μ L MLB to  22.5 μ L sample in RIPA in 1.5 ml LoBind safety lock tube.
- 28 Denature at  95 $^{\circ}$ C for  00:05:00 .
- 29 Transfer to ice.

SDS PAGE: *Loading and Running Gel*



- 30 Remove gel from package, take off green strip at the bottom, and rinse the wells three times with about 1 mL running buffer (use 1 ml pipettor).
- 31 Place gel, tall plate facing out in outer side of a holder and the reservoir block at the other to make a running buffer reservoir. Put the other gel holder in place to take up space.
- 32 Fill the chamber to the 2-gel mark. Make sure there are no leaks before testing the circuit. Test the circuit by checking for bubbles after turning on power to **70 volts**.
- 33 Using a 20 ul conte-tipped pipettor, load 30 μ L of the samples in MLB onto gel, as well as 10 μ L of Biorad Precision Plus protein ladder Dual Color.
- 34 Run at 100 – 120V, depending on desired resolution, until adequate separation of ladder lanes in sizes regions of interest.
- 35 Photograph the gel after the run with your phone to help keep track of orientation.

Semi-dry Transfer to PVDF Membrane

- 36 Pre-soak blotting stacks in 1X transfer buffer made according to directions on the bottle.
- 37 Pre-wet PVDF membrane from *Biorad Transblot kit* in 190 proof ethanol, then in transfer buffer.
- 38 Layer into Transblot drawer (stack, membrane, gel, second stack) before rolling out bubbles gently with conical tube
- 39 Run the transfer on Biorad setting, standard Mininigel TGX, 25mA, 25V, 00:03:00 .
- 40 While running, pour 5% milk block in 1X PBST directly into cleaned small black Western blot box.
- 41 Incubate transferred blot for 01:00:00 at Room temperature in 5% milk block in 1X PBST.
- 42 Photograph the blot after the transfer to help keep track of orientation.



Antibody Incubations and ECL visualization of target protein bands: *Primary and secondary antibody binding*

43 Primary and secondary antibody binding



- β -Actin (C4): sc-47778
- ECL SuperSignal West Pico Chemiluminescent Substrate

Use 1:500 Dilutions for primary antibodies, in 5% milk/PBST



Note




Note: Primary antibody (**unconjugated with HRP**) can be re-used 3x if kept frozen.



44 Drain milk block from blots and add primary Ab (in milk block).



45 Incubate in cold room  Overnight covered and shaking on orbital shaker, covered with plastic. 

46 **The next day**, wash blot **3x**/  00:15:00 in PBST. 

47 Freshly prepare secondary antibody at 1:10,000 dilution ( 2 μ L to  20 mL 5% milk/PBST).

48 Incubate in secondary antibody for  01:00:00 ,  Room temperature . 

49 Discard secondary and wash blot **3x**/  00:05:00 in 1xPBST. 

50 Discard secondary and in the same box or moving the blot to a new box, add ECL working solution ( 1 mL of both reagents mixed in a  15 mL conical and applied to blot.)

51 Incubate in ECL for  00:05:00 , shaking occasionally by hand, covered. 




- 52 Drain blot and place between plastic sheets for imaging using Biorad GelDoc, as below. After imaging, the blot can be probed for actin as a protein loading control without stripping the blot, using anti-actin-HRP as described below.

Antibody Incubations and ECL visualization of target protein bands: *Actin staining with anti-actin-HRP primary*

- 53 *Actin staining with anti-actin-HRP primary*



Wash blot 3x with 1X PBST < 00:10:00 .

- 54 Incubate 00:20:00 ,  Room temperature in the primary actin antibody 1:500 freshly made dilution in 5% Milk block.



- 55 Wash 2 – 3x in 1x PBST 00:05:00 each.



- 56 Incubate in ECL as above (step 51) and image.



Antibody Incubations and ECL visualization of target protein bands: *Imaging Blot on BioRad ChemiDoc XRS+ (Vollmers lab)*

- 57 *Imaging Blot on BioRad ChemiDoc XRS+ (Vollmers lab)*

Place blot in between sheet protector plastic and transport to imager in cassette to keep dark.

- 58 Place blot on white screen.

- 59 Place screen into chamber drawer.

- 60 Select new protocol, then blot, then chem hi sensitivity and position blot under live focus setting.

- 61 Under Applications, select chemidoc hi sensitivity to photograph target bands. Note the image size, the gain (2X), and bin (2×2) to make sure the image size is the same when the protein standards are photographed in order to be able to merge the images.



- 62 Under Live Acquire, select acquisition settings: 1, 600, for range and total =100 acquisitions (the 100 can be lowered for longer exposure), starting at 0.25 start time for high expression target.
- 63 Freeze and save the image of your choice as it comes up.
- 64 **Without moving the blot**, take a picture to visualize the protein standards by selecting Custom under Applications, then create an Epi illumination protocol with the same gain and binning as used for photographing the bands (double check the image size prediction under the options...)
- 65 Merge the target bands and the protein standards pictures in Image Lab software using Image tools.

Note

Note: Imaged blots can be stored at **4°C in PBST** for stripping and re-probing. Some folks also store them frozen in **-20°C** flat in plastic bag.