

Jan 12, 2026

Bulk Astrocyte RNA Seq Extraction

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g77y29gwz/v1>

Mackenzie Ferguson¹

¹Drexel University College of Medicine

MF Space



Mackenzie Ferguson

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g77y29gwz/v1>

Protocol Citation: Mackenzie Ferguson 2026. Bulk Astrocyte RNA Seq Extraction. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.q26g77y29gwz/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: January 10, 2026

Last Modified: January 12, 2026

Protocol Integer ID: 238407

Keywords: seq, rna, extraction, bulk

Disclaimer

Only draft form

Abstract

Bulk Astrocyte RNA Seq preparation

Cell Preparation

- 1 Assume microglia and OPCs have been removed prior to this procedure, if not refer to proper protocols to have pure astrocyte culture
- 2 Remove media and discard
- 3 Rinse twice with PBS, aspirate the PBS, add 5 ml trypsin-EDTA (.05%) and incubate in the CO₂ incubator at 37 °C. Check detachment of astrocytes every 5 min and enforce detachment of astrocytes by hitting the flask against the palm of your hand (2-3 times)
- 4 After astrocytes are detached from the culture flask, add 5 ml of astrocyte culture media to inactivate trypsin
- 5 Cell Counting
 - 5.1 Take a previously prepared 1.5 mL microtube (with  10 µL of trypan blue) and combine with  10 µL of cell suspension
 - 5.2 Load a hemacytometer with 10 µL of the cells and trypan blue solution and examine immediately under a microscope at low magnification.
 - 5.3 Count the number of viable cells. Those will light up white and are not invaded with the blue dye
 - 5.4 To calculate the number of viable cells per mL of culture, use the formula: Number of viable cells × 10⁴ = cells/mL culture. Remember to correct for the dilution factor.
 - 5.5 For Plasmidsaurus e cell couFor Plasmidsaurus bulk RNA seq, recommended cell count per well is 2*10⁵
- 6 Transfer proper amount of cell suspension to microtubes. Spin cells at 180 x g for 5 min, aspirate supernatant and add  50 µL of Zymo 1x DNA/RNA shield and resuspend
- 6.1 If the lysed cell solution appears highly viscous, vortex it vigorously for about 30s to shear chromosomal DNA.



- 7 Transfer prepared supernatant to properly labeled  200 μ L strip tubes or 96 well plate according to <https://plasmidsaurus.com/sample-prep/rna#section1> instructions