

Feb 18, 2026

Bulk ATAC-seq

DOI

<https://dx.doi.org/10.17504/protocols.io.dm6gp12kdgzp/v1>

Lexi Bounds¹, Alejandro Barrera¹, Maria ter Weele¹, Siyan Liu¹, Euphy Wu², Shengyu Li¹, Revathy Venukuttan¹, Ruhi Rai¹, Wancen Mu², Nahid Iglesias¹, Paola Giusti-Rodriguez³, Timothy Reddy¹, Yun Li², Raluca Gordan², Andrew Allen², Michael Love², Patrick Sullivan², Greg Crawford¹, Charles Gersbach¹

¹Duke University; ²University of North Carolina; ³University of Florida

Gersbach Lab



Andrea R Daniel

Duke University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gp12kdgzp/v1>

Protocol Citation: Lexi Bounds, Alejandro Barrera, Maria ter Weele, Siyan Liu, Euphy Wu, Shengyu Li, Revathy Venukuttan, Ruhi Rai, Wancen Mu, Nahid Iglesias, Paola Giusti-Rodriguez, Timothy Reddy, Yun Li, Raluca Gordan, Andrew Allen, Michael Love, Patrick Sullivan, Greg Crawford, Charles Gersbach 2026. Bulk ATAC-seq. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.dm6gp12kdgzp/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 17, 2026

Last Modified: February 18, 2026

Protocol Integer ID: 243537

Keywords: bulk atac, using neural progenitor cell, neural progenitor cell, npc, seq, seq method, seq this protocol

Funders Acknowledgements:

NIH

Grant ID: HG012053

Disclaimer

This protocol was adapted from the work of Lexi Bounds in the Gersbach lab at Duke University.

Abstract

This protocol describes bulk ATAC-seq methods and analysis using Neural Progenitor Cells (NPCs).

Materials

KAPA Library Quantification Kit (Roche #7960255001)

High Sensitivity DNA ScreenTape Analysis (Agilent #5067)

Bulk ATAC-seq

- 1 100,000 NPCs are harvested using Accutase and resuspended in 1 mL of cold ATAC-seq resuspension buffer and the Omni-ATAC seq protocol was used to generate ATAC-sequencing libraries (Ref. 1).
- 2 Following pre-amplification, 3 additional cycles of PCR were performed prior to final amplification and cleanup.
- 3 Library concentrations and sizes were obtained using the KAPA Library Quantification Kit (Roche #7960255001) and High Sensitivity DNA ScreenTape Analysis (Agilent #5067).
- 4 Perform ATAC-seq for two biological replicates. The final libraries are pooled to 3 nM and sequenced on a HiSeq 4000 with 1×50bp configuration.

Bulk ATAC-seq analysis

- 5 The pipeline for processing and analysis of the bulk ATAC-seq data can be found here: https://github.com/Duke-GCB/GGR-cwl/blob/master/v1.0/ATAC-seq_pipeline/pipeline-se-blacklist-removal.cwl.
- 6 Quality control of FASTQ files is performed using FASTQC (Ref.2).
- 7 Adapter reads are then trimmed using Trimmomatic.
- 8 Reads are then aligned to the reference genome using Bowtie (Ref. 3).
- 9 Duplicate reads are removed using Picard MarkDuplicates and ENCODE hg38 blacklist reads were removed using BEDtools2 (Ref. 4).
- 10 Peak calling is performed using MACS2 (Ref. 5) with narrowPeak settings.
- 11 Sequencing-depth normalized ATAC bigWig files were generated using deeptools (Ref. 6) bamCoverage.



Protocol references

1. Corces, M. R. *et al.* An improved ATAC-seq protocol reduces background and enables interrogation of frozen tissues. *Nat. Methods* **14**, 959–962 (2017).
2. Babraham Bioinformatics - FastQC A Quality Control tool for High Throughput Sequence Data.
<https://www.bioinformatics.babraham.ac.uk/projects/fastqc/>.
3. Langmead, B. Aligning short sequencing reads with Bowtie. *Curr. Protoc. Bioinforma.* **Chapter 11**, Unit 11.7 (2010).
4. Quinlan, A. R. BEDTools: The Swiss-Army Tool for Genome Feature Analysis. *Curr. Protoc. Bioinforma.* **47**, 11.12.1–34 (2014).
5. Zhang, Y. *et al.* Model-based analysis of ChIP-Seq (MACS). *Genome Biol.* **9**, R137 (2008).
6. Ramírez, F., Dündar, F., Diehl, S., Grüning, B. A. & Manke, T. deepTools: a flexible platform for exploring deep-sequencing data. *Nucleic Acids Res.* **42**, W187–191 (2014).