

Jul 16, 2019

Version 1

C. elegans bleaching solution preparation V.1

DOI

dx.doi.org/10.17504/protocols.io.3kkgkuw

Cristian Riccio¹

¹University of Cambridge



Cristian Riccio

University of Cambridge

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3kkgkuw>

Protocol Citation: Cristian Riccio 2019. C. elegans bleaching solution preparation. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.3kkgkuw>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 01, 2019



Last Modified: July 16, 2019

Protocol Integer ID: 23916

Keywords: nematode, bleaching, *C. elegans*, *Caenorhabditis elegans*, worm, synchronisation, synchronization, decontamination, bleaching solution preparation, bleaching solution, resistant to the bleaching, worms at different developmental stage, fungal cell, hatched worm, bacteria, most fungi, larvae

Abstract

The *C. elegans* bleaching solution is used for several purposes:

- destroying and removing bacterial or fungal cells from a contaminated *C. elegans* population. This is possible because *C. elegans* eggs are resistant to the bleaching used in this protocol but most fungi and bacteria are not.
- synchronising a *C. elegans* population consisting of worms at different developmental stages. This is possible because eggs are resistant to the bleaching used in this protocol but hatched worms (larvae and adults) are not.

Materials

MATERIALS

✕ Disposable gloves, nitrile

✕ sodium hydroxide 10 M

✕ DEPC water

✕ Sodium hypochlorite solution **Honeywell Fluka Catalog #239305-500ML**

✕ 50 ml Falcon tube

Troubleshooting

Safety warnings

- ! Sodium hypochlorite and sodium hydroxide are highly toxic and should be handled with care, wearing a lab coat and goggles. Read relevant information.

Before start

Prepare the bleaching solution.

Get a plate of gravid (= pregnant) worms or a plate with lots of unhatched eggs if you want to carry out an egg prep.



- 1 Add 42 ml DEPC water to a 50 ml Falcon tube. 3m
- 2 Add 3 ml 10 M NaOH 3m
- 3 Add 5 ml sodium hypochlorite solution 3m
- 4 See relevant protocol for details but 6 ml of bleaching solution for 5 minutes on a vortex in 14 ml tubes should be sufficient to dissolve worms and spare eggs.