



Jan 02, 2025

C. elegans culture and maintenance

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlk8y6wl5r/v1>

Justin Donnelly¹

¹Arcadia Science

Arcadia Science

c elegans



Arcadia Science

Arcadia Science

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlk8y6wl5r/v1>

Protocol Citation: Justin Donnelly 2025. *C. elegans* culture and maintenance. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.j8nlk8y6wl5r/v1>

Manuscript citation:

Avasthi P, Borges AL, Cheveralls K, Donnelly J, Mets DG, Reiter T. (2025). An experimental and computational workflow to characterize nematode motility behavior. <https://doi.org/10.57844/arcadia-b89a-7e76>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 23, 2024

Last Modified: March 13, 2025

Protocol Integer ID: 106301

Keywords: C. elegans, cryopreservation, procedure for cryopreservation, cryopreservation, revival of this organism, homozygous strain, maintenance this protocol, maintenance

Abstract

This protocol describes the maintenance of *C. elegans* WT and homozygous strains, as well as a procedure for cryopreservation and revival of this organism.

Materials

6 cm Petri dishes with NGM(OP50)

Spatula

Platinum wire

1X trehalose freezing buffer

Cryo-tubes

 Ethanol 200 Proof **Decon Labs Catalog #2716**

Protocol materials

 Ethanol 200 Proof **Decon Labs Catalog #2716**



C. elegans culture protocol (WT and homozygous mutants)

1 Grow worms in Petri dishes on nematode growth medium (NGM) with OP50 *E. coli* at  20 °C .

2 For optimal health, subcultivate worms every 2–4 days.

2.1 Subcultivate via chunking

Dip a spatula in 100%  Ethanol 200 Proof **Decon Labs Catalog #2716** . Flame-sterilize. Once the flame is extinguished, cut a chunk from the plate with an edge length of ~  0.5 cm . Select a chunk with several larvae. Transfer this chunk to a new plate under flame, with the worm side facing down.

Note

If worms are dauer, you can use this same procedure to revive them from the dauer plate.

2.2 Subcultivate via picking

Flame-sterilize a platinum wire. Pick up one or more adult worms from plate, taking care to prevent puncturing the NGM. Transfer worms by gently touching the platinum wire to the new plate under flame. Hold the wire in this position for several seconds to allow worms to crawl off. Flame-sterilize the wire before picking additional worms. Repeat 1–4 times, depending on the desired density.

3 Return worms to incubator.

Freezing & thawing *C. elegans*

3m

4 Prepare synchronized L1s via bleach synchronization. See:



Protocol

NAME

Bleach life-stage synchronization of *C. elegans*

CREATED BY

Arcadia Science

[Preview](#)

- 5 Centrifuge on clinical centrifuge at speed 5 for 00:01:00 – 00:02:00 . 3m
- 6 Resuspend worms in 15 mL 1× trehalose freezing buffer. Centrifuge again on clinical centrifuge at speed 5 for 00:01:00 – 00:02:00 . 3m
- 7 Resuspend worms in 3.6 mL 1× trehalose freezing buffer. Mix to ensure worms are dispersed.
- 8 Quickly after resuspension, pipette 1.2 mL into each of three pre-labeled cryotubes.
- 9 Freeze the worms in a Mr. Frosty (or similar) at -80 °C . Keep the container on its side so that worms don't all sink to the bottom before freezing.
- 10 After several days, retrieve the worms and run a test thaw by scraping some amount of frozen worm and buffer, and transfer to a fresh NGM (OP50) plate.
- 11 Over the following several days, look to see if any worms have revived. If so, you may remove cryotubes from the Mr. Frosty and transfer to long-term storage at -80 °C .
- 12 Once you've confirmed successful revival, you can thaw worms as in step 10.



Protocol references

Trehalose freezing buffer recipe:

<http://wbg.wormbook.org/2017/07/11/an-improved-method-for-freezing-nematodes-using-a-trehalose-dmso-cryoprotectant/>