



Feb 28, 2023

Version 2

Calcium fluorimetry with the FLIPR Calcium 6 kit on FlexStation 3 V.2

 [PLOS One](#)

DOI

dx.doi.org/10.17504/protocols.io.3byl4jdmjlo5/v2

Angus Li¹, Samuel Liu¹, Rennica Huang², Seungkirl Ahn¹, Robert J Lefkowitz^{1,2,3}

¹Department of Medicine, Duke University Medical Center, Durham, North Carolina, United States of America;

²Department of Biochemistry, Duke University Medical Center, Durham, North Carolina, United States of America;

³Howard Hughes Medical Institute, Duke University Medical Center, Durham, North Carolina, United States of America



Angus Li

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4jdmjlo5/v2>



Protocol Citation: Angus Li, Samuel Liu, Rennica Huang, Seungkirl Ahn, Robert J Lefkowitz 2023. Calcium fluorimetry with the FLIPR Calcium 6 kit on FlexStation 3. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.3byl4jdmjlo5/v2> Version created by **Angus Li**

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 27, 2023

Last Modified: February 28, 2023

Protocol Integer ID: 77743

Keywords: calcium fluorimetry with the flipr calcium, calcium fluorimetry, flipr calcium, receptor, coupled receptor

Funders Acknowledgements:

Howard Hughes Medical Institute

National Institutes of Health

Grant ID: R01 HL16037

Abstract

This protocol details an experimental procedure used to generate results described in the manuscript Li, A., Liu, S., Huang, R., Ahn, S., & Lefkowitz, R. J. (2023). Loss of biased signaling at a G protein-coupled receptor in overexpressed systems.

Troubleshooting

- 1 Plate U2OS-TetOn-AT1R in microplates (Black with clear bottom, lysine- coated Corning 3842) at 15000 cells/well 2 days prior to assay
- 2 Add doxycycline and optionally PTX to wells 14 hours before replacement with loading buffer
- 3 Remove one vial of Calcium 6 Assay Reagent (Component A, good for 2 plates) from the freezer and equilibrate to room temperature.
- 4 Dissolve contents of one Component A vial by adding 10 ml of 1X HBSS Buffer plus 20 mM HEPES. Mix by vortexing (~1-2 min) until contents of vial are dissolved. It is important that contents are completely dissolved to ensure reproducibility between experiments.
- 5 Remove cell plates from the incubator, Calcium 6 Kit - add an equal volume of loading buffer to each well in the dark (i.e. 120ul HBSS+HEPES + 40ul of Dye/well (4X) for a 96-well plate in the original protocol). Use 11.8ml of HBSS+HEPES +4.8ml of Dye (4X) for one plate. Save the dissolved dye solution in -20.
- 6 Return plates to the incubator and incubate two hours at 37 C, 5% CO₂
- 7 Optionally, add YM-254890 1 hour prior to reading
- 8 30 minutes before experiment (1.5 hours into incubation), prepare the reagents (5X, 40ul/well+20ul volume). The ligand needs to be plated in a 96 well with round bottom (The FlexStation will add ligand, shake, and measure sequentially).
- 9 Turn on FlexStation 30 minutes before experiment and warm up to 37.
- 10 Place tips and ligand plate
- 11 After the two-hour incubation period, transfer the assay plate directly to the FlexStation instrument assay plate carriage 10 minutes prior to reading to re-equilibrate temperature, and then run the assay.



- 12 In an individual well or column of wells, the calcium flux peak(s) should be complete within 1 to 3 minutes after addition. For an entire plate however, the protocol will not complete until all chosen columns are finished. 40-45 min for entire 96 well plate.