



Nov 11, 2022

## CALCIUM STAINING PROTOCOL

DOI

[dx.doi.org/10.17504/protocols.io.kxygx9zewg8j/v1](https://dx.doi.org/10.17504/protocols.io.kxygx9zewg8j/v1)

Machlusil Husna<sup>1</sup>, Kusworini Handono<sup>2</sup>, Hidayat Sujuti<sup>3</sup>, Aulanni'am<sup>4</sup>, Ettie Rukmigarsari<sup>5</sup>

<sup>1</sup>Department of Neurology Faculty of Medicine, Universitas Brawijaya / dr. Saiful Anwar General Hospital, Malang, Indonesia;

<sup>2</sup>Department of Clinical Pathology Faculty of Medicine, Universitas Brawijaya / dr. Saiful Anwar General Hospital, Malang, Indonesia;

<sup>3</sup>Department of Ophthalmology Faculty of Medicine, Universitas Brawijaya / dr. Saiful Anwar General Hospital, Malang, Indonesia;

<sup>4</sup>Department of Chemistry, Faculty of Science Universitas Brawijaya, Malang, Indonesia;

<sup>5</sup>Faculty of Teacher Training and Education, Islamic University of Malang, Indonesia



Machlusil Husna

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygx9zewg8j/v1>



**Protocol Citation:** Machlusil Husna, Kusworini Handono, Hidayat Sujuti, Aulanni'am, Ettie Rukmigarsari 2022. CALCIUM STAINING PROTOCOL. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kxygx9zewg8j/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 18, 2022

**Last Modified:** November 11, 2022

**Protocol Integer ID:** 71480

**Keywords:** hippocampal excitotoxicity, hippocampal excitotoxicity process, immediate effect after kainic acid exposure, kainic acid exposure, protocol neurodegeneration due to neurotoxicity, organotypic hippocampal slice culture, neurotoxicity, staining protocol neurodegeneration, neurodegeneration process, temporal lobe epilepsy, phenomena in temporal lobe epilepsy, neuronal death, increased calcium, astrocyte reactivity

## Abstract

Neurodegeneration due to neurotoxicity is one of the phenomena in temporal lobe epilepsy. Experimentally, hippocampal excitotoxicity process can occur due to kainic acid exposure, especially in the CA3 area. Neuronal death, astrocyte reactivity and increased calcium also occur in hippocampal excitotoxicity, but few studies have investigated immediate effect after kainic acid exposure. The organotypic hippocampal slice culture (OHSC) is a useful model for studying the neurodegeneration process, but there are still many protocol differences. In this study, minor modifications were made in the OHSC protocol.

## Troubleshooting



## Calcium assay kit

- 1 We use Fluo-4 assay kit (calcium ab 228555) that provides a homogenous fluorescence-based assay for detecting the intracellular calcium mobilization

## Material supplied in kit

- 2 Fluo-4 AM 1 vial  
10x F127 Plus 1 bottle (10 mL)  
HHBS (Hank's Buffer with 20 mM Hepes) 1 bottle (100 mL)

## Reagent Preparation

- 3 Briefly centrifuge small vials at low speed prior to opening  
Thaw all the kit components at room temperature before use  
Fluo-4 AM stock solution : add 200  $\mu$ L of DMSO into the vial of Fluo-4 AM and mix well  
1x Assay buffer : make 1x assay buffer by adding 1 mL of 10x F127 plus into 9 mL of HHBS buffer and mix them well  
Fluo-4 AM dye-loading solution : add 20  $\mu$ L of Fluo-4 AM stock solution into 10 mL of 1x assay buffer and mix them well. This working solution is stable for at least 2 hours at room temperature

## Assay procedure

- 4 Equilibrate all materials and prepared reagents to room temperature just prior to use and gently agitate  
Add 20  $\mu$ L Fluo-4 AM dye-loading solution into the cell plate  
Incubate the dye-loading plate in incubator with temperature 37°C for 1 h  
Prepare the plate with HHBS  
Run the calcium flux assay by monitoring the fluorescence intensity at Ex/Em = 490/525 nm with CLSM