


Sep 23, 2021

Version 2

🌐 Calibration of glass electrode half-cells V.2

DOI

dx.doi.org/10.17504/protocols.io.eq2lpyrmlx9/v2



Agnes Heering¹, Ivo Leito¹

¹University of Tartu



Agnes Heering

University of Tartu

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2lpyrmlx9/v2>

External link: <https://analytical.chem.ut.ee/>

Protocol Citation: Agnes Heering, Ivo Leito 2021. Calibration of glass electrode half-cells. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.eq2lpyrmlx9/v2> Version created by **Agnes Heering**

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: September 23, 2021

Last Modified: September 23, 2021

Protocol Integer ID: 53498

Keywords: electrode calibration procedure with aqueous buffer, electrode calibration procedure, calibration of glass, calibration, cells the purpose, glass, cell

Abstract

The purpose of this document is give glass-electrode calibration procedure with aqueous buffers.

Troubleshooting

Before start

Turn on the thermostat and electrometer at least half an hour before starting the work.

Software

- 1 Start Quick IV Measurement Software.

Note


Computer cannot go to sleep during measurements or the communication between computer and instruments is lost and data collection stops.

- 2 On the left hand pick Function "Source & Sampling".



Position of the Function choice "Source & Sampling".

- 3 Right hand side click on tab "Channel Communication". Click on "Search Channels", which opens a new window. Select USB interface and search for the channel. After the search click on the channel and close the window.




Position of the tabs “Channel Communication” and “Search Channels”.

Note

This must be done every time the program is started.

- 4 Settings can be saved as QIVM file, which can be used for measurements. The saved file can be opened under tab “File”. The settings will remain the same between the measurements if the program is not closed.



Program settings.

Note

Range is "Fixed 20V" in case of EST-0601 calibration.

Filling the cell

8m

5



🌡 25.0 °C

Equipment

| | |
|--------------|-------|
| VistaShield | NAME |
| Faraday cage | TYPE |
| Gamry | BRAND |
| VistaShield | SKU |

Note

The side ports must be open.

**Note**

A beaker can also be used if temperature is controlled.

6 Fill capillary with standard aqueous buffer with pH 7.00.  50 mL

2m

7 Insert glass electrode and a reference electrode.

1m

Equipment

| | |
|-----------------------------|-------|
| K401 | NAME |
| Calomel reference electrode | TYPE |
| Radiometer | BRAND |
| K401 | SKU |

7.1 Electrode is rinsed with water, gently dried, and then rinsed with solution to be measured.

1m

8 Connect electrodes to the instrument.

Equipment

B2987A Electrometer / High Resistance Meter NAME

Electrometer TYPE

Keysight BRAND

B2987A SKU

<https://www.keysight.com/zz/en/home.html> LINK

- 8.1 Crocodile clips are as following:
red is signal and connects to the glass electrode,
black is shield and is connected to reference electrode and
green is guard. Green stays unused.

Note

Cable goes through the port in Faraday cage. Electrometer is not in the Faraday cage.

Equipment

N1415A Triax to Alligator Cable NAME

Cable TYPE

Keysight Technologies BRAND

N1415A SKU

<https://www.keysight.com/zz/en/product/N1415A/triax-alligator-cable-200-v-1-5-m.html> LINK



Measurement

1h

9 Start the measurement by clicking on the "Measure" button on the upper panel.

10 Data collection. Point is taken at 10 s interval.

1h

Note

Calibration duration must be the same as later used in pH measurements.

11 Save data. Files are named as Date_GE_vs_ref_pHx.

Note


SAVE THE FILE BEFORE NEW MEASUREMENT! Otherwise, data is lost. Data is given in volts and seconds.


11.1 Go to table view, right click on the table, choose "Export as CSV..." and save the results. The data can be automatically imported into excel file while measuring or after (both under "Excel Navigation"), but this causes constant shift in the comma position.

Washing

12 Remove the electrodes.

13 Remove the solution.

14 Rinse three times with plenty of water  20 mL each

15 Rinse once with acetone to ease the drying.  2 mL

- 16 Dry cell with compressed air.

Measurement


- 17 Repeat measurements with standard buffers pH 4 and pH 10.

Data analysis

- 18 Points from 30 min to 60 min are used for analysis.

- 19 Take the average of the chosen data points.

- 20 Plot ΔE vs pH.



Calibration graph of a glass electrode half-cell.

- 21 Obtain the slope and intercept of the glass electrode-half cell.