


May 06, 2019

## Calibration Protocol - Particle Standard Curve with Microspheres

DOI

[dx.doi.org/10.17504/protocols.io.zgjf3un](https://dx.doi.org/10.17504/protocols.io.zgjf3un)



Richard Tennant<sup>1</sup>, Paul Rutten<sup>1</sup>

<sup>1</sup>iGEM Measurement Committee

iGEM Measurement

Tech. support email: [pauljrutten@gmail.com](mailto:pauljrutten@gmail.com)



Paul Rutten

The University of Oxford

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.zgjf3un>

External link: <https://2019.igem.org/Measurement>

**Protocol Citation:** Richard Tennant, Paul Rutten 2019. Calibration Protocol - Particle Standard Curve with Microspheres. protocols.io <https://dx.doi.org/10.17504/protocols.io.zgjf3un>

## Manuscript citation:

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**Created:** March 24, 2019

**Last Modified:** August 07, 2019

**Protocol Integer ID:** 21739

**Keywords:** particle standard curve with microsphere, nm absorbance measurement, optical characteristics of these microsphere, standard curve of particle concentration, microsphere, calibration protocol, calibration, plate reader, particle concentration, amount of particle, particle standard curve, measurement, particle, optical characteristic, measure

## Abstract

You will prepare a dilution series of monodisperse silica microspheres and measure the Abs<sub>600</sub> in your plate reader.

The size and optical characteristics of these microspheres are similar to cells, and there is a known amount of particles per volume. This measurement will allow you to construct a standard curve of particle concentration which can be used to convert 600 nm absorbance measurements into an estimated equivalent number of cells.

## Attachments




iGEM Data Analysis T...


25KB

## Materials


### MATERIALS


 96 well plate

 double distilled water (ddH<sub>2</sub>O)

 300µl Silica beads

### STEP MATERIALS

 ddH<sub>2</sub>O

 300µl Silica beads

300 µL Silica beads are provided in the iGEM Measurement Kit. The 96-well plate should preferably be black with a clear flat bottom.



## Protocol materials

⊗ ddH2O

⊗ 300μl Silica beads

⊗ 96 well plate

⊗ double distilled water (ddH2O)

⊗ 300μl Silica beads

⊗ 300μl Silica beads

⊗ ddH2O


## Troubleshooting

## Before start

Read through this entire protocol carefully before you start your experiment and prepare any materials you may need.


## Prepare the Microsphere Stock Solution

- 1 Obtain the tube labeled "Silica Beads" from the Measurement Kit and vortex vigorously for 30 seconds.

 300µl Silica beads

### Note

**Microspheres should NOT be stored at 0°C or below**, as freezing affects the properties of the microspheres. If you believe your microspheres may have been frozen, please contact the iGEM Measurement Committee for a replacement (measurement@igem.org).


- 2 Immediately pipet 100 µL microspheres into a 1.5 mL eppendorf tube
- 3 Add 900 µL of ddH<sub>2</sub>O to the microspheres  
 ddH<sub>2</sub>O
- 4 Vortex well. This is your Microsphere Stock Solution

## Prepare the serial dilution of microspheres

- 5 Accurate pipetting is essential. Serial dilutions will be performed across columns 1-11. **Column 12 must contain ddH<sub>2</sub>O only.**

Initially you will setup the plate with the microsphere stock solution in column 1 and an equal volume of 1x ddH<sub>2</sub>O in columns 2 to 12.

You will perform a serial dilution by consecutively transferring 100 µl from column to column with good mixing.



- 6 Add 100  $\mu$ l of ddH<sub>2</sub>O into wells A2, B2, C2, D2....A12, B12, C12, D12
- 7 Vortex the tube containing the stock solution of microspheres vigorously for 10 seconds
- 8 Immediately add 200  $\mu$ l of microspheres stock solution into A1
- 9 Transfer 100  $\mu$ l of microsphere stock solution from A1 into A2
- 10 Mix A2 by pipetting up and down 3x and transfer 100  $\mu$ l into A3
- 11 Mix A3 by pipetting up and down 3x and transfer 100  $\mu$ l into A4
- 12 Mix A4 by pipetting up and down 3x and transfer 100  $\mu$ l into A5
- 13 Mix A5 by pipetting up and down 3x and transfer 100  $\mu$ l into A6
- 14 Mix A6 by pipetting up and down 3x and transfer 100  $\mu$ l into A7



- 15 Mix A7 by pipetting up and down 3x and transfer 100  $\mu$ l into A8
- 16 Mix A8 by pipetting up and down 3x and transfer 100  $\mu$ l into A9
- 17 Mix A9 by pipetting up and down 3x and transfer 100  $\mu$ l into A10
- 18 Mix A10 by pipetting up and down 3x and transfer 100  $\mu$ l into A11
- 19 Mix A11 by pipetting up and down 3x and transfer 100  $\mu$ l into liquid waste

**Note**

**Take care not to continue serial dilution into column 12**

- 20 Repeat dilution series for rows B, C, D

**21 IMPORTANT!**

Re-Mix (pipette up and down) each row of your plate **immediately before** putting in the plate reader! (This is important because the beads begin to settle to the bottom of the wells within about 10 minutes, which will affect the measurements.)

**Note**

**Take care to mix gently and avoid creating bubbles on the surface of the liquid**

**Measure OD**

- 22 Measure OD<sub>600</sub> of all samples in instrument
- 23 Record the data in your notebook



24 Import data into this Excel sheet:



iGEM Data Analysis Template - Part...

## **Congratulations!**

---

25 You have now completed this calibration protocol