


Aug 23, 2023

CAMbank: cfDNA BCT Field Processing v1

DOI

dx.doi.org/10.17504/protocols.io.6qpvr33yovmk/v1



Elijah G Overbey^{1,2}, Krista A Ryon¹, JangKeun Kim¹, Christopher E Mason^{1,2}

¹Weill Cornell Medicine; ²BioAstra

Cornell Aerospace Medi...



Elijah G Overbey

University of Austin (UATX), Weill Cornell Medicine, BioAstr...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvr33yovmk/v1>

Protocol Citation: Elijah G Overbey, Krista A Ryon, JangKeun Kim, Christopher E Mason 2023. CAMbank: cfDNA BCT Field Processing v1. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.6qpvr33yovmk/v1>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 22, 2023

Last Modified: June 27, 2024

Protocol Integer ID: 86810

Keywords: CAMbank, astronaut, SOMA, cfDNA BCT, Streck, plasma, cfdna bct field processing v1 field processing, cfdna bct, cornell aerospace medicine biobank, cambank, rbc pellet

Funders Acknowledgements:

WorldQuant Foundation

NASA

Grant ID: NNX14AH50G

NASA

Grant ID: NNX17AB26G

NASA

Grant ID: 80NSSC22K0254

NASA

Grant ID: NNH18ZTT001N-FG2

NASA

Grant ID: 80NSSC22K0254

NASA

Grant ID: NNX16AO69A

NASA

Grant ID: 80NSSC21K0316

Abstract

Field processing of cfDNA BCTs for the Cornell Aerospace Medicine Biobank (CAMbank).

Instructions for preserving: plasma and RBC Pellets.

Materials

Tube Type: Streck Cell-Free DNA BCT[®] RUO (Streck: #230470)

Troubleshooting



Perform Venipuncture

5m

- 1 After venipuncture, invert the tubes gently 8 to 10 times to fully mix tube anticoagulant with blood sample.

5m

Store the tube upright at room temperature until centrifugation.

Note: The tube stabilizes cell-free DNA for up to 14 days at 6 °C to 37 °C and CTCs for up to 7 days at 15 °C to 30 °C.

Centrifuge Settings: Plasma Separation

25m

- 2 Note: A **swing bucket** centrifuge is required.

3m

Set centrifuge:

- acceleration: 9
- deceleration: 9
- temperature: RT
- duration: 20 minutes
- speed: 300xg

- 3 Place the cfDNA BCTs in the centrifuge.

22m

Place a protective cover over the swing buckets in case of tube leakage.

Start the centrifuge.

Stand by the centrifuge until the centrifuge reaches max speed. Listen for signs of imbalance or compromised tube integrity.

Separate Plasma

7m

- 4 Carefully remove cfDNA BCTs from the centrifuge and inspect for separation of red blood vs plasma layers.

2m

Transfer cfDNA BCTs to a sanitized laminar flow hood.

- 5 Remove the upper plasma layer and transfer to a new 15mL conical tube. Cap the conical tubes and return to the centrifuge.

5m

Centrifuge Settings: Plasma Purification

15m

- 6 Note: A **swing bucket** centrifuge is required.

3m

Set centrifuge:



- acceleration: 9
- deceleration: 9
- temperature: RT
- duration: 10 minutes
- speed: 5000xg

- 7 Place the conical tubes in the centrifuge. 12m
- Place a protective cover over the swing buckets in case of tube leakage.
- Start the centrifuge.
- Stand by the centrifuge until the centrifuge reaches max speed. Listen for signs of imbalance or compromised tube integrity.

Freeze Red Blood Cell Pellets

5m

- 8 Recap and freeze the cfDNA tubes to preserve DNA in the red blood cell pellet 5m

Aliquot Plasma

17m

- 9 Carefully remove conical tubes from centrifuge. 2m
- Transfer conical tubes to a sanitized laminar flow hood.
- 10 Using a P1000, carefully aliquot plasma into 2D barcoded tubes at **500µL each**. 15m
- Return 2D tubes with plasma to the rack and place in the -80C freezer.