

May 06, 2019

🌐 Carbon and nitrogen stable isotope analysis in harbor seal pup fur at CICESE

📖 [PLOS One](#)

DOI


dx.doi.org/10.17504/protocols.io.2m3gc8n

Juan Carlos Herguera-García¹, Gisela Heckel¹

¹Centro de Investigación Científica y de Educación Superior de Ensenada, Mexico (CICESE)



Gisela Heckel



OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.2m3gc8n

External link: <https://doi.org/10.1371/journal.pone.0225889>

Protocol Citation: Juan Carlos Herguera-García, Gisela Heckel 2019. Carbon and nitrogen stable isotope analysis in harbor seal pup fur at CICESE. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.2m3gc8n>

Manuscript citation:

Juárez-Rodríguez, M., Heckel, G., Herguera-García, J.C., Elorriaga-Verplancken, F.R., Herzka, S.Z., Schramm, Y. Trophic ecology of Pacific harbor seals in Mexican colonies based on $\delta^{13}\text{C}$ and $\delta^{15}\text{N}$ stable isotopes. Submitted to PLOS ONE.

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 06, 2019

Last Modified: May 06, 2019

Protocol Integer ID: 22939

Keywords: Stable isopes, Pacific harbor seal, *Phoca vitulina richardii*, pup fur,

Abstract

In the Stable Isotope Laboratory at CICESE (Ensenada Center for Scientific Research and Graduate Education, Mexico) animal tissues (pup fur from Pacific harbor seals, *Phoca vitulina richardii*) were ground for isotopic determinations. Samples were stored in a desiccator until they were loaded into a Costech[®] Carrousel. Samples were combusted at high temperature (1000°C) in a pure oxygen atmosphere in an elemental analyzer coupled to a Finnigan MAT Delta V Advantage continuous flow stable isotope mass spectrometer.

Stable carbon and nitrogen ratios are expressed as $\delta^{13}\text{C}$ or $\delta^{15}\text{N}$ according to the following equation:

$$\delta^{13}\text{C} \text{ or } \delta^{15}\text{N} (\text{‰}) = [\text{R}_{\text{sample}}/\text{R}_{\text{standard}} - 1] * 1000 \text{ (‰)}$$

where R is $^{13}\text{C}/^{12}\text{C}$ or $^{15}\text{N}/^{14}\text{N}$.

Isotopic values of carbon and nitrogen are reported relative to Pee Dee Belemnite and atmospheric nitrogen standards, respectively. The accuracy of isotopic measurements was verified using secondary standard reference materials (Glutamic acid, a pure carbonate Merck and the laboratory internal references Lanugo for N, and CH94 for C).

Attachments



CICESE_Stable_Isotop...

13KB

