

Aug 27, 2025

Carver et al., Aged Brain Spatial Profiling - Immunofluorescence Imaging

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vzx6drgx1/v1>

Chase Carver¹

¹Mayo Clinic

Cellular Senescence Net...



Chase Carver

Mayo Clinic

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzx6drgx1/v1>

Protocol Citation: Chase Carver 2025. Carver et al., Aged Brain Spatial Profiling - Immunofluorescence Imaging. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.rm7vzx6drgx1/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: September 29, 2023

Last Modified: August 27, 2025

Protocol Integer ID: 88602

Keywords: modifiable features of aged brain white matter, aged brain white matter, old brain accumulation of lipofuscin granule, aged brain spatial profiling, immunofluorescence imaging, associated microglia, immunofluorescence imaging this protocol, floating mouse brain section, autofluorescent noise, old brain accumulation, autofluoresce, mouse brain section, fluorophore, immunohistochemistry, mouse brain, traditional immunohistochemistry, brain slice, immunostaining process, lipofuscin granule

Funders Acknowledgements:

NIH Cellular Senescence Network

Grant ID: UG3 CA275669-01

Abstract

This protocol provides the preparation and staining process to perform traditional immunohistochemistry on free-floating mouse brain section used in in Carver et al., "Senescent- and disease-associated microglia are modifiable features of aged brain white matter". The immunostaining process involves unconjugated- or fluorophore-conjugated- primary antibodies and secondary antibodies conjugated to fluorophores. Due to old brain accumulation of lipofuscin granules that autofluoresce, we note that in order to achieve maximal signal, a quenching step is performed to extinguish the autofluorescent noise.

Sections were cut at 30 μ m with a Leica CM3050 S cryostat

Brain slices were incubated in free-floating solutions in 12-well or 24-well plates

Imaging was performed with a Nikon Ti2 Eclipse inverted microscope and Orca Fusion BT sCMOS camera.



Block and permeabilization

2h 25m

- 1 Transfer 30 μ m free-floating brain section to well containing 8% donkey serum, 0.1% Triton-X, 0.1% Tween-20 in PBS. Incubate for 2 hours at 25C. 2h
- 2 Remove blocking solution and wash sections thoroughly in cold PBS 5 times for 5 minutes each on shaker 25m

Quench autofluorescence

35m

- 3 Dilute stock 40X True Black Plus aqueous (Biotium, catalogue #23014) to 1X in PBS. Incubate sections in 1X True Black Plus for 10 min. Note: any residual detergents may interfere with quenching efficacy 10m
- 4 Remove quenching solution and wash sections thoroughly in cold PBS 5 times for 5 minutes each on shaker 25m

Primary antibody stain

12h

- 5 Stain sections in primary antibody diluted in solution of 8% donkey serum in PBS overnight at 4C on a shaker rotating at 50 rpm. Antibody combinations consist of:
rabbit polyclonal anti-IBA1 (Fujifilm Wako, #019-19741, RRID:AB_839504)
dilution: 1:500
goat polyclonal anti-IBA1 (Abcam, #ab5076, RRID:AB_2224402)
dilution: 1:200
rat IgG_{2A} monoclonal anti-galectin-3 clone eBioM3/38 (Invitrogen, #14-5301-85, RRID:AB_837133)
dilution: 1:200
goat polyclonal anti-GFAP (Abcam, #ab53554, RRID:AB_880202)
dilution: 1:2000
rabbit polyclonal anti-apolipoprotein E clone EPR19392 (Abcam, #ab183597, RRID:AB_2832971)
dilution: 1:4000 12h
- 6 Remove antibody solution and wash sections thoroughly in cold PBS 3 times for 5 minutes each on shaker 15m

Secondary antibody stain

2h 15m



- 7 Stain sections with donkey-host secondary antibodies conjugated to fluorophores AF488, AF594, or AF647 (Jackson ImmunoResearch) for 2 hours at 25C in dark room environment. Antibodies are diluted 1:250 in PBS. 2h
- 8 Remove antibody solution Wash sections thoroughly in cold PBS 3 times for 5 minutes each on shaker 15m

Mount sections onto slides

10m

- 9 Mount sections onto Superfrost Plus microscope slides, let dry for 10 min. in dark 10m
- 10 Add Vectashield with DAPI(Vector Laboratories, #H-1200) to slide and cover with a 1.5 glass coverslip

Microscopy

- 11 Image slides with Nikon Ti2 Eclipse Inverted microscope using 10X, 20X, and 40X Plan Apo objectives
Detect each fluorophore channel with 8-channel Spectra III light engine and Orca Fusion BT sCMOS camera in Nikon Elements AR software.
Illuminate with laser lines of 405, 488, 594, and 647 nm with an image exposure time of 200 ms in 16-bit readout mode.
Process fluorescence intensity, colocalization, and morphometry data with Fuji + ImageJ
Calculate total corrected cellular fluorescence intensity for each cell from the formula:
mean integrated density - (mean local background X cell area)