

Jun 19, 2024

Cas9 enrichment for DiMeLo-sequencing

DOI

<https://dx.doi.org/10.17504/protocols.io.5jyl8p8n7g2w/v1>

Pragya Sidhwani¹, Aaron Straight²

¹Stanford University; ²Stanford University, Department of Biochemistry



Pragya Sidhwani

Stanford University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl8p8n7g2w/v1>

Protocol Citation: Pragya Sidhwani, Aaron Straight 2024. Cas9 enrichment for DiMeLo-sequencing. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.5jyl8p8n7g2w/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 13, 2023



Last Modified: June 19, 2024

Protocol Integer ID: 92204

Keywords: cas9 enrichment protocol for nanopore, cas9 enrichment for dimelo, cas9 enrichment protocol, cas9 enrichment, enrichment of dna, read sequencing, sgrna, cas9, genomic loci of interest, sequencing, genomic loci, dna interaction, nanopore, seq protocol, dna

Abstract

This protocol builds upon two previously published protocols: [this](#) DiMeLo-seq protocol to map protein-DNA interactions at a single molecule level using long-read sequencing and [this](#) Cas9 enrichment protocol for Nanopore sequencing to enrich for genomic loci of interest for DiMeLo-sequencing. Some unique features of our protocol include:

- using sgRNAs instead of duplexing crRNA and and tracrRNA for Cas9 enrichment
- Cas9-enrichment of DNA containing methyl-A marks
- multiplexing of Cas9-enriched samples
- Potential to use >10 sgRNAs (while we have not tested this, we think our edits to step 4 enable this)

Introduction

- 1 This protocol builds upon two previously published protocols: [this](#) DiMeLo-seq protocol to map protein-DNA interactions at a single molecule level using long-read sequencing and [this](#) Cas9 enrichment protocol for Nanopore sequencing to enrich for genomic loci of interest. Some unique features of our protocol include:
 - using sgRNAs instead of duplexing crRNA and and tracrRNA for Cas9 enrichment
 - Cas9-enrichment of DNA containing methyl-A marks
 - multiplexing of Cas9-enriched samples
 - Potential to use >10 sgRNAs (while we have not tested this, we think our edits to step 4 enable this)
- 2 Designing sgRNAs for Cas9 enrichment
 - 2.1 We use [CRISPOR](#) to design sgRNAs to target the loci of interest. We have tested this protocol by cutting with 2 sgRNAs each on either end (so 4 sgRNAs total) of a 10 kb centromeric region.
 - 2.2 Instead of duplexing a crRNA and tracrRNA, we order 1.5 nmol sgRNAs from [Synthego](#), which are resuspended in 15 uL 1x TE for 100 uM.

Isolating pA-Hia5-treated DNA from human K562 cells

1h

- 3 Extract DNA that has been methylated *in situ* by pA-Hia5 at sites of protein-DNA interactions
 - 3.1 Use the NEB Monarch HMW DNA Extraction Kit. Follow protocol for genomic DNA isolation using cell lysis buffer. Include RNase A. Perform lysis with 2000 rpm agitation. We have validated 2000 rpm gives N50 ~50-70 kb but if longer reads are desired we expect 300 rpm would work.

Apart from using a different kit, all of the steps for the long fragment DNA extraction are the same as the general protocol. To reiterate, make the following changes to the protocol:
 - 3.2 If fixation was performed, be sure to do the 56°C incubation for lysis for 1 hour (not just 10 minutes) to reverse crosslinks. Agitate for 10 minutes and then keep at 56°C without agitation for 50 minutes.



- 3.3 Quantify DNA yield by Qubit dsDNA BR Assay Kit. We need 3-10 ug of high quality DNA in less than 24 uL of 1x TE.

Preparing RNPs for Cas9 enrichment

45m

- 4 Preparing the ribonucleoprotein (RNP) complex for Cas9 enrichment
- 4.1 For each sgRNA, assemble the RNP complex as follows:
0.4 uL IDT HiFi Cas9 (62 uM)
0.5 uL Synthego sgRNA (100 uM)
Assemble each RNP for 30 minutes at room temperature, keep on ice till ready to use
- 4.2 Prior to adding to dephosphorylated DNA from step 5, pool and dilute the RNPs as follows:
3.6 uL Pooled RNPs (for 4 sgRNAs)
5 uL 10x rCutSmart
41.4 uL water
This results in 50 uL of a 100x solution enough for 4-5 samples (10-12 uL/sample)

Cas9-based enrichment of pA-Hia5 methylated DNA

2h 40m

- 5 Dephosphorylate DNA isolated in step 3
- 5.1 In a PCR tube, mix the following:
3 uL 10x rCutSmart buffer
3-10 ug of DNA
3 uL of NEB QuickCIP
water up to 27 uL

Mix by resuspending with a cut pipet tip
- 5.2 Incubate at
37 C for 60 mins (dephosphorylation)
80 C for 2 mins (enzyme inactivation)
hold at 20 C

1h



6 RNP-based cleavage and A-tailing

6.1 To the dephosphorylated DNA, add the following:

12 uL of pooled 100x RNP complex

1 uL of 10 mM dATP

1 uL Taq polymerase

Mix by resuspending with a cut pipet tip

6.2 Incubate at:

37 C for 60 minutes (Cas9 cleavage)

72 C for 30 minutes (monoadenylation of cut ends)

Hold at 12 C

1h 30m

7 Do a cleanup using the Zymogen genomic DNA clean and concentrator following the manufacturer's protocol. Elute DNA in 30 uL EB that has been pre-heated.

10m

Barcoding and library prep

10m

8 To 500 ng of each sample, add barcodes from the SQK-LSK 109 kit and 2x TA/ligase mix

38 uL DNA (500 ng)

4 uL barcode

42 uL 2x TA/ligase mix

10m

9 Pool samples and do a cleanup using the Zymogen genomic DNA clean and concentrator kit into 65 uL of EB. Proceed with adaptor ligation following ONT protocol.