

May 09, 2019

Case - Intraperitoneal Insulin Tolerance Test

DOI

dx.doi.org/10.17504/protocols.io.yefftbn



Henri Brunengraber¹

¹Case Western Reserve University

Mouse Metabolic Phenotyping Centers
Tech. support email: info@mmpc.org



Lili Liang

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.yefftbn

External link: <https://mmpc.org/shared/document.aspx?id=272&docType=Protocol>

Protocol Citation: Henri Brunengraber 2019. Case - Intraperitoneal Insulin Tolerance Test. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.yefftbn>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 21, 2019

Last Modified: May 09, 2019

Protocol Integer ID: 20647

Keywords: insulin tolerance testing

Abstract

Summary:

This is the standard protocol for routine insulin tolerance testing. It is performed on awake mice, fasted for 18 hr (overnight) or for 6 hours.

Materials

MATERIALS

⊗ Insulin Eli Lilly

⊗ 1cc syringe 26g Becton Dickinson (BD)

⊗ Glucometer Alphasat Abbott Catalog #cfigs185-m0285

⊗ Alpha trak2 Glucose strips Abbott Catalog #Alpha trak2 Glucose strips

⊗ Clotisol Amazon Catalog #T212002

Reagents and Materials:

Reagent/Material	Quantity Required	Vendor	Stock Number
Insulin	0.75U-1.0U per mouse	Eli Lilly	
1cc syringe 26g	1/mouse	BD	
Glucometer Alphasat	1	Abbott	cfigs185-m0285
Alpha trak2 Glucose strips	5 per mouse	Abbott	
Clotisol	100ul	Amazon	T212002

Note:

BD Biosciences, RRID:SCR_013311

Abbott, RRID:SCR_010477

- 1 Fast mice for overnight or for 6 hours. Remove mouse from cage and put into a clean cage with water and no food.
- 2 Label the mice using a marker on their tail to quickly tell them apart
- 3 Weigh each mouse and record weight. Calculate the amount of insulin needed for 0.75U to 1.0 U /kg body weight. Record the volume for each mouse.
- 4 Dilute the stock of insulin (100U/ml) to 0.1u/ml
- 5 Take a glucose reading and record it as time zero. Dip tail into clotisol to clot blood.
- 6 Inject each mouse 3' apart. Repeat #4 at 15min, 30 min, 60 min after IP injection of insulin.
- 7 Return mice to cage with food.