



Feb 10, 2025

Caspase-3/7 activity

 [Nature Cell Biology](#)

 In 2 collections

DOI

<https://dx.doi.org/10.17504/protocols.io.8epv52ow4v1b/v1>

Ching-Chieh Chou¹, Judith Frydman¹

¹Stanford University



Ching-Chieh Chou

Stanford University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv52ow4v1b/v1>

External link: <https://doi.org/10.1038/s41556-025-01623-y>

Protocol Citation: Ching-Chieh Chou, Judith Frydman 2025. Caspase-3/7 activity. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.8epv52ow4v1b/v1>



Manuscript citation:

Chou, CC., Vest, R., Prado, M.A. *et al.* Proteostasis and lysosomal repair deficits in transdifferentiated neurons of Alzheimer's disease. *Nat Cell Biol* **27**, 619–632 (2025). <https://doi.org/10.1038/s41556-025-01623-y>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 07, 2025

Last Modified: February 10, 2025

Protocol Integer ID: 119805

Keywords: ASAPCRN, detecting cell apoptosi, cell apoptosi

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000282

Aligning Science Across Parkinson's

Grant ID: ASAP-024268

Abstract

Detecting cell apoptosis by measuring Caspase-3/7 activation

Protocol for detecting Caspase-3/7 activation

- 1 Caspase-3/7 substrate FAM-DEVD-FMK (ImmunoChemistry) is reconstituted in DMSO and stored at -20°C. The FAM-DEVD-FMK reagent is one of the fluorochrome-labeled inhibitors of caspases that covalently and irreversibly binds to the active caspases. The green fluorescence intensity is a direct measurement of Caspase-3/7 activity.
- 2 To test toxic or rescuing effects on cell apoptosis, the cells are pre-treated with DMSO or any pharmacological perturbations.
- 3 When cells are available, first dilute the 150X stock concentrate 1:5 in PBS to form the staining solution at 30X.
- 4 Add the 30X staining solution at a dilution of 1:30 to cell culture medium to form 1X staining solution. The solution should be used immediately.
- 5 Prepare 1X Apoptosis Buffer by diluting 10X Apoptosis Buffer 1:10 in dH₂O.
- 6 Add 100-300 µL FAM-FLICA solution to each well and incubate for 30 min at 37°C.
- 7 Aspirate the FAM-FLICA staining solution and rinse 1X Apoptosis Buffer once.
- 8 The cells are fixed by 4% paraformaldehyde (PFA) for 15 min and stained with antibodies and 5 µg/mL Hoechst for confocal microscope imaging.