

Mar 21, 2024

Cassiopea xamachana Cellular Dissociation

Forked from a private protocol

DOI

dx.doi.org/10.17504/protocols.io.3byl4qnnovo5/v1

Anthony Bonacolta¹, Victoria Sharp², Marta Mammone³

¹University of Miami; ²Pennsylvania State University; ³Texas A&M University


Cassiopea protocols

Tech. support phone: +1 (814) 321-5684 email: momedinamunoz@gmail.com



Victoria Sharp

Pennsylvania State University



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4qnnovo5/v1>

Protocol Citation: Anthony Bonacolta, Victoria Sharp, Marta Mammone 2024. Cassiopea xamachana Cellular Dissociation. [protocols.io https://dx.doi.org/10.17504/protocols.io.3byl4qnnovo5/v1](https://dx.doi.org/10.17504/protocols.io.3byl4qnnovo5/v1)



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 09, 2024

Last Modified: March 21, 2024

Protocol Integer ID: 94965

Keywords: dissociation, cassiopea, cytometry, calcein, cnidaria, cnidarian, enzyme, fluorescence, cassiopea xamachana cellular dissociation this protocol, cassiopea xamachana cellular dissociation, cassiopea xamachana cells for cell, cassiopea xamachana cell, viable cell, cassiopea xamachana, cell, lyse cell, scrna

Abstract

This protocol is to optimized to dissociate and fix *Cassiopea xamachana* cells for cell sorting and scRNA-seq. The dissociation by itself results in 53-55% of viable cells. Cells cannot be sorted without fixative, unless your machine can sort a seawater solution. Any other solution will lyse cells.

Guidelines

Make sure to work in an RNase free-environment when able to. Use RNase-ZAP or work in a UV sterilized hood if possible.

The tissue should be less than 1 cm long



Materials

- Sterile razors
- Sterile forceps
- Ice
- Wide-bore pipette tips, or cut 1000 uL pipette tips
- 15 mL tube for digestion buffer
- 2 petri dishes to wash and incubate tissue in Ca^{2+} and Mg^{2+} free seawater
- 2 70-um filters
- 2 mL microcentrifuge tubes
- pipettes

Equipment:

microcentrifuge

rocking plate

Reagents:

NaCl

KCl

NaSO_4

NaHCO_3

Dispase II

Liberase

L-cysteine

PBS

BSA

DNase/RNase-free distilled water

methanol

glacial acetic acid

glycerol

RNAse Inhibitor

Troubleshooting



Before start

Treat reagents and materials with UV-light for ~15 mins before beginning protocol.

Set 15 mL and microcentrifuge to 4° C.

Prepare Reagents:

Ca²⁺ Mg²⁺ free seawater (Roger et al. 2021)

To 1 L Distilled Water add:

- 23 g NaCl
- 0.763 g KCl
- 3 g NaSO₄
- 0.25 g NaHCO₃

- Dissociation Mix:

To Ca²⁺ Mg²⁺ free ASW add:

- 3.6 mg/mL Dispase II
- 0.25 mg/mL of Liberase
- 4% L-cysteine

- 1x PBS 0.5% bovine serum albumin (BSA)

- Add 0.25 g to 50 mL 1x PBS.

- Fresh ACME Solution

- 13:3:3:2 ratio of DNase/RNase-free distilled water, methanol, glacial acetic acid, and glycerol
- Prep about 15 mL, FRESH, per sample each time

- RNase Inhibitor



Dissociation

1h

- 1 Cut the jellyfish tissue with a sterile razor to encourage permeability of reagents. 2m
- 2 Gently wash the jellyfish tissue in 10 mL Ca-Mg-Free SW for 00:01:00 then transfer to fresh 10 mL Ca-Mg-Free SW and let incubate at Room temperature for 00:02:30 . 3m 30s
- 3 Using sterilized forceps, place the jelly tissue into a clean 15 mL tube then add 1 mL dissociation mix on top of the jelly, or enough to submerge the tissue. 30s
- 4 Incubate the tissue on a rocker for 00:30:00 at room temperature. 30m
- 5 Pipette up and down using a wide-bore tip 10 times. 5m
- 6 Repeat steps 4 and 5. 35m
- 7 Add 80 μ L fetal bovine serum to the cell suspension to create a 8% FBS solution to halt enzyme digestion. 2m
- 8 After the incubation period, filter the sample through a 70 μ m filter Keep sample On ice moving forward. 2m
- 9 Resuspend in 1000 μ L Ca-Mg-free SW . Gently pipette up and down 10 times with a wide-bore tip to dissociate clumps. 2m

Staining and Fixation

- 10 Add 1 micromolar (μ M) Calcein Violet 450 AM to cells. Incubate in the dark for 00:30:00 . 30m



- 11 After the incubation period, filter the sample through a 70 μm filter 5m
- 12 Resuspend in 500 μL ACME . Incubate for 00:30:00 . 35m
- 13 Centrifuge 350 x g, 4°C, 00:07:00 then carefully discard the supernatant. 7m
- 14 Re-suspend the pellet in 800 μL 1x PBS w/ 1% BSA using gentle pipetting with a wide-bore pipette then add 1 μL RNA Inhibitor . 1m

FACS & Cryopreservation


- 15 Pre-chill chambers of FACS Machine to 4 °C 5m
- 16 Sort at the slowest rate (High-purity) with less than 50 PSI at 4 °C . Gate for Calcein Violet (450 nm) and chlorophyll autofluorescence (650–700 nm) for viable jelly cells and symbiont cells. 2h
- 17 After sorting, cryopreserve by adding 10 % volume DMSO and 1 μL RNA Inhibitor and immediately putting the sample at -80 °C . 1m

Thaw and Sample Submission

38m

- 18 Thaw the sample On ice 30m
- 19 Centrifuge 350 x g, 4°C, 00:07:00 then carefully remove the supernatant. 7m
- 20 Re-suspend the cells in 1 mL 3.3x PBS w/ 1% BSA then add 1 μL RNA Inhibitor . 1m



21 Submit to sequencing center  On ice .

Protocol references

https://github.com/thesyntheticcoral/SynCoral_Protocols/tree/master/Culturing

Roger LM, Reich HG, Lawrence E, Li S, Vizgaudis W, Brenner N, et al. (2021) Applying model approaches in non-model systems: A review and case study on coral cell culture. PLoS ONE 16(4): e0248953.

<https://doi.org/10.1371/journal.pone.0248953>