



Aug 23, 2022

## cDNA Library Preparation for scRNA-seq of Human Meniscus (10x Genomics)

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvkw32vmk/v1>

molmer<sup>1</sup>, Martin Lotz<sup>1</sup>, Tony Mondala<sup>1</sup>, Steven Head<sup>1</sup>

<sup>1</sup>Scripps Research

Human BioMolecular Atlas Program (HuBMAP) Method Development Community

Tech. support email: [Jeff.spraggins@vanderbilt.edu](mailto:Jeff.spraggins@vanderbilt.edu)



**Merissa Olmer**

Scripps Research

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.e6nvkw32vmk/v1>

**Protocol Citation:** molmer , Martin Lotz, Tony Mondala, Steven Head 2022. cDNA Library Preparation for scRNA-seq of Human Meniscus (10x Genomics) . **protocols.io** <https://dx.doi.org/10.17504/protocols.io.e6nvkw32vmk/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 10, 2022

**Last Modified:** August 23, 2022

**Protocol Integer ID:** 68486

**Keywords:** 10x genomic, chemistry kits of 10x genomic, cdna library preparation, scrna, barcoded antibody, cells from individual donor, cell hashing, cdna, barcoding chemistry kit, cell

## Abstract

Samples are processed using V2 barcoding chemistry kits of 10x Genomics. For each run, 10,000 cells from individual donors are labeled with distinct oligo-barcoded antibodies (cell hashing), enabling us to pool samples from both conditions on each 10X Genomics run and reduce batch effects. All samples from a given experiment are processed in parallel in the same thermal cycler. The 10x user guide can be found [https://assets.ctfassets.net/an68im79xiti/1C16trEdzy1Folq5xbOijE/7e6fb1f504e130bd561d898384da99d9/CG000315\\_ChromiumNextGEMSingleCell3-\\_GeneExpression\\_v3.1\\_DualIndex\\_RevB.pdf](https://assets.ctfassets.net/an68im79xiti/1C16trEdzy1Folq5xbOijE/7e6fb1f504e130bd561d898384da99d9/CG000315_ChromiumNextGEMSingleCell3-_GeneExpression_v3.1_DualIndex_RevB.pdf) and is also an attached document.

## Attachments



CG000315\_ChromiumNe

X...

4.5MB

## Materials

Chromium Next GEM Single Cell 3' GEM Kit v3.1, 16 rxns PN-1000123  
Library Construction Kit, 16 rxns PN-1000190  
Chromium Next GEM Single Cell 3' Gel Bead Kit v3.1, 16 rxns PN-1000122  
Dynabeads™ MyOne™ SILANE PN-2000048  
Chromium Next GEM Chip G Single Cell Kit, 16 rxns PN-1000127  
Dual Index Kit TT Set A, 96 rxns PN-1000215  
BD Hu Single Cell Sample Multiplexing Kit cat# 633781



- 1 Single cells are isolated from human knee meniscus using doi [dx.doi.org/10.17504/protocols.io.n2bvj6k7blk5/v1](https://doi.org/10.17504/protocols.io.n2bvj6k7blk5/v1).
- 2 Single cell suspensions are multiplexed and converted to barcoded scRNAseq libraries using the Chromium Single Cell 3' Library, Gel Bead and Chip Kit (10x Genomics), and the **BD™ Hu Single Cell Sample Multiplexing Kit**.

#### Note

Chromium Next GEM Single Cell 3' Reagent Kits v3.1 (Dual Index) User Guide  
[https://assets.ctfassets.net/an68im79xiti/1C16trEdzy1Folq5xbOijE/7e6fb1f504e130bd561d898384da99d9/CG000315\\_ChromiumNextGEMSingleCell3-\\_GeneExpression\\_v3.1\\_DualIndex\\_\\_RevB.pdf](https://assets.ctfassets.net/an68im79xiti/1C16trEdzy1Folq5xbOijE/7e6fb1f504e130bd561d898384da99d9/CG000315_ChromiumNextGEMSingleCell3-_GeneExpression_v3.1_DualIndex__RevB.pdf)

- 3 Libraries are sequenced on an Illumina NextSeq2000 sequencer targeting 20,000-25,000 reads per cell.