

Oct 20, 2019

cDNA synthesis

DOI

dx.doi.org/10.17504/protocols.io.8hhht36

Igem Dusseldorf¹

¹Heinrich-Heine Universität Düsseldorf



Igem Dusseldorf

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.8hhht36

Protocol Citation: Igem Dusseldorf 2019. cDNA synthesis. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.8hhht36>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 20, 2019

Last Modified: October 20, 2019

Protocol Integer ID: 28937

Materials

200 ng of RNA

preferred cDNA synthesis Kir



- 1 for cDNA synthesis use sythesis Kit for e.g. the BioRad iScript cDNA synthesis Kit (Cat.# 1708891)
- 2 use 200 ng of RNA for reaction
- 3 the reaction protocol reverse transcription was carried out at 42° C for 30 min
- 4 The reverse transcriptaseinactivation was performed at 85° C for 5 min
- 5 used directly or stored at -20° C