

Oct 02, 2019

## CDO expression into OnePot PURE

DOI

[dx.doi.org/10.17504/protocols.io.7nhhmb6](https://dx.doi.org/10.17504/protocols.io.7nhhmb6)

Laura Kvedarauskaite<sup>1</sup>

<sup>1</sup>EPFL - EPF Lausanne

iGEM EPFL



Laura Kvedarauskaite

---

OPEN  ACCESS



**DOI:** [dx.doi.org/10.17504/protocols.io.7nhhmb6](https://dx.doi.org/10.17504/protocols.io.7nhhmb6)

**Protocol Citation:** Laura Kvedarauskaite 2019. CDO expression into OnePot PURE. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.7nhhmb6>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development

**We are still developing and optimizing this protocol**

**Created:** September 25, 2019

**Last Modified:** October 02, 2019

**Protocol Integer ID:** 28073

**Keywords:** OnePot PURE, CDO, Catechol, Cell Free

## Abstract

This protocol explains the procedure for expressing the catechol degrading enzyme *Catechol-2,3-deoxygenase* (CDO) into a "homemade" OnePot PURE cell-free transcription/translation system.

Once expressed, this enzyme degrades catechol, a colorless substrate, into a yellow colored one, 2-hydroxymuconate semialdehyde (2-HMS), providing a colorimetric signal that can be easily implemented in another protocol.

## Guidelines

We used a CDO plasmid with A T7 promoter added to it.

## Materials

	<b>Materials μl</b>	<b>CDO expression (1)</b>	<b>Control without plasmid (2)</b>	<b>Control without catechol (3)</b>
	Energy solution	2	2	2
	Ribosomes	0.9	0.9	0.9
	Proteins	0.65	0.65	0.65
	CDO Plasmid	25 ng	-	25 ng
	Catechol 10mM	0.5	0.5	-
	Water	up to 5 μl	up to 5 μl	up to 5 μl
	<b>Total</b>	<b>5</b>	<b>5</b>	<b>5</b>

## Before start

Preheat the incubator at 37°C



- 1 Label 3 PCR tubes according to the reactions
- 2 In each tubes add the Energy solution, proteins and ribosomes.
- 3 Add the catechol, CDO and water as needed according to the materials chart.
- 4 Make a quick spin in the centrifuge to have all the liquid in the bottom.
- 5 Incubate one hour at 37°C