



Apr 12, 2023

Cell culture and Western blot

 [The Journal of Biological Chemistry](#)

DOI

dx.doi.org/10.17504/protocols.io.kxygx9bozg8j/v1

Andrea Dickey¹, rabrisch¹

¹UCSD Department of Cellular and Molecular Medicine



Yu Xuan Lin

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygx9bozg8j/v1>

External link: <https://dx.doi.org/10.1016/j.jbc.2024.107469>

Protocol Citation: Andrea Dickey, rabrisch 2023. Cell culture and Western blot. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.kxygx9bozg8j/v1>

Manuscript citation:

Dederer V, Murillo MS, Karasmanis EP, Hatch KS, Chatterjee D, Preuss F, Azeez KRA, Nguyen LV, Galicia C, Dreier B, Plückthun A, Versees W, Mathea S, Leschziner AE, Reck-Peterson SL, Knapp S (2024) A designed ankyrin-repeat protein that targets Parkinson's disease-associated LRRK2. The Journal of Biological Chemistry 300(7). doi: [10.1016/j.jbc.2024.107469](https://doi.org/10.1016/j.jbc.2024.107469)



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 11, 2023

Last Modified: April 12, 2023

Protocol Integer ID: 80356

Keywords: rab7a phosphorylation at s72, detecting rab7a phosphorylation, lrrk1 mutant, s72, lrrk1, cell culture

Funders Acknowledgements:

ASAP

Grant ID: ASAP-000519

Abstract

Protocol for detecting Rab7a phosphorylation at S72 by LRRK1 and LRRK1 mutants in HEK 293T cells.

Troubleshooting




Cell culture and transfection

1d 0h 15m

1 Day - 1

1d

Split 293T cells in 6-well dishes  24:00:00 before transfection at 200K/well in

 2 mL DMEM (with FBS and with P/S)

2 Day - 2

Transfection

2.1 Warm PEI and Opti-MEM at RT for 00:15:00

15m

Label a tube for each well.

2.2 Add DNA (500 ng of **GFP-Rab7** + 1 µg of **LRRK1 construct of interest**) to 150 µL of Opti-MEM and gently vortex for 5 sec

2.3 Add 3 µL room-temperature PEI to each tube containing Opti-MEM/DNA mix. Mix gently by pipetting and incubate DNA/DMEM/PEI mixture at room temperature in the hood for 00:15:00 .

15m

2.4 In the meantime...remove media from 6-well dish containing cells and replace with fresh 1 mL DMEM (with FBS and withoutP/S)...put back in 37 deg incubator until ready to add transfection mixture

2.5 After 15 min incubation, add DNA/Optimem/PEI mixture dropwise 150 µL onto each well Give a bit of a light swirl before putting back in incubator

Cell lysis








1d 12h

3 36:00:00 after transfection, begin cell lysis

1d 12h

3.1 Wash plate on ice with cold PBS 1x






- 3.2 Add  300 μL RIPA buffer (0.5% Triton, 50 mM Tris pH7.5, 150 mM NaCl, 0.1% SDS) with protease and phosphatase inhibitors (cOmplete mini EDTA free + PhosSTOP tablets)
- 3.3 Lift with cell lifters on ice
- 3.4 Pipet up lysate, put in eppendorf tube, and shake 15 mins in the cold room
- 3.5 Spin at MAX at 4 deg for 15 mins
- 3.6 Remove supernatant and make sample (boil at  95 $^{\circ}\text{C}$ for  00:10:00), store in  -80 $^{\circ}\text{C}$. *I like to store the lysate and take some out to make a sample –I use the 4x NuPage LDS sample buffer:*  65 μL lysate +  10 μL 10x Reducing Agent +  25 μL 4x LDS buffer

10m






Western blot

1h

- 4 SDS-PAGE with Bis-Tris gel and MOPS running buffer
- 4.1 Load a 4-12% Bis-Tris gel with  25 μL of prepared lysate in sample buffer and run at 180V for ~  00:50:00 or until dye front has reached the bottom of the gel.
- 4.2 Assemble gel with Immobilon-FL PVDF membrane for transfer according to instructions from your western blot transfer apparatus. When using fluorescence detection for Western blot it is important to use low fluorescence background membrane (Immobilon-FL or equivalent).
- We activate membrane with MeOH and rinse with water.
Transfer in Western transfer buffer with Tris/Glycine and 20% MeOH at 200 mA for 4 hr at  4 $^{\circ}\text{C}$
- 4.3 After transfer is complete, rinse membrane with water and allow to dry between sheets of Whatman paper.

50m



- 4.4 Block in 5% milk in TBS (no Tween 20)
- 4.5 Dilute primary antibody in 5% milk with TBST (with Tween 20)
-rabbit anti Rab7 phospho-S72 (MJF-38) at 1:1000
-mouse anti-GFP at 1:2500 (Santa Cruz) (for total Rab quantification)
-rabbit anti-LRRK1 (ab228666) at 1:500
-rabbit anti-GAPDH at 1:3000 (Cell signaling technologies)
- 4.6 Rock  Overnight at  4 °C
- 4.7 Rinse 3x with TBST for  00:05:00
- 4.8 Rinse 1x with 5% milk in TBST
- 4.9 Add secondary antibodies in 5% milk with TBST and incubate at  Room temperature
for 1 hr
-LiCor mouse and rabbit secondary IRdye antibodies at 1:5000
- 4.10 Rinse 3x with TBST for  00:05:00
- 4.11 Image on LiCor Odyssey CLx

5m

5m