

Apr 18, 2022

Cell culture, transfection, and imaging

DOI

dx.doi.org/10.17504/protocols.io.36wgq41j3vk5/v1

Will Hancock-Cerutti^{1,2,3}, Pietro De Camilli^{1,3}

¹Departments of Neuroscience and of Cell Biology, Howard Hughes Medical Institute, Program in Cellular Neuroscience, Neurodegeneration and Repair, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

²Interdisciplinary Neuroscience Program and MD-PhD Program, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

³Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815



William Hancock-Cerutti

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.36wgq41j3vk5/v1

Protocol Citation: Will Hancock-Cerutti, Pietro De Camilli 2022. Cell culture, transfection, and imaging. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.36wgq41j3vk5/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 08, 2021

Last Modified: May 31, 2024

Protocol Integer ID: 51390

Keywords: Cell culture, Transfection, Spinning-disk confocal, Fluorescent microscopy, ASAPCRN, culturing hela cell, cell culture, hela cell, transfection, transient transfection, andor dragonfly, cell, imaging, disk confocal system, spinning disk

Abstract

This protocol describes general procedures for culturing HeLa cells, transient transfection, and imaging using an Andor Dragonfly spinning disk confocal system.

Attachments



[dxwubkjz7.pdf](#)

142KB

Materials

DMEM Solution:

A	B
FBS	10%
Penicillin	100 U/ml
Streptomycin	100 mg/ml
L-glutamine	2 mM

* (all from Gibco)

General preparation

- 1 Culture the HeLa-M cells at  37 °C in 5% CO₂ and DMEM containing 10% FBS,  100 U/ml penicillin,  100 mg/mL streptomycin, and  2 millimolar (mM) L-glutamine (all from Gibco).

Note

Note: For general maintenance, when cells reached 80-90% confluency, they were deattached from the dish with Trypsin and diluted 1:20 in a new dish.

- 2 For live-cell imaging experiments, seed the cells on glass-bottomed dishes (MatTek) at a concentration of 35,000 cells per dish and transfect after  24:00:00 using FuGene HD (Promega) in Opti-MEM (Gibco).

1d

- 3 Image the cells  24:00:00 after transfection.

1d



- 4 Just before imaging, remove the growth medium and replace with prewarmed live-cell imaging solution (Life Technologies).

- 5 For lysotracker experiments, incubate the cells in  50 nanomolar (nM) LysoTracker Red DND-99 (ThermoFisher) in complete DMEM for  00:30:00, wash twice with media, then image in live-cell imaging solution.

30m



- 6 Perform all live-cell imaging at  37 °C and 5% CO₂.

- 7 Perform spinning-disk confocal microscopy using an Andor Dragonfly system equipped with a plan apochromat objective (63×, 1.4 NA, oil) and a Zyla scientific CMOS camera.



- 8 For any given experiment, use the same exposure time, laser power, and gain for image acquisition to allow for quantitative comparison.



Imaging of cells stably expressing STING-GFP

3d 14h

- 9 Generate the cells stably expressing STING-GFP as described elsewhere.
- 10 Culture the stable STING-GFP HeLa-M cells at  37 °C in 5% CO₂ and DMEM containing 10% FBS,  100 U/ml penicillin,  100 mg/mL streptomycin, and  2 millimolar (mM) Lglutamine (all from Gibco).
- 11 For experiments using siRNA, transfect 60 pmols of the indicated siRNA using  6 µL Lipofectamine RNAiMax (ThermoFisher) in Opti-MEM (Gibco) per dish according to manufacturer protocol. Image the cells  72:00:00 after siRNA transfection.
- 12 For experiments using cGAMP, transfect  50 µg/L of cGAMP using  18 µL Lipofectamine RNAiMax (ThermoFisher) in Opti-MEM (Gibco) per dish according to manufacturer protocol. Image the cells  14:00:00 after transfection.

3d



14h

