



Mar 19, 2024

Cell culture, transfection, immunocytochemistry, and imaging

DOI

dx.doi.org/10.17504/protocols.io.5qpvokx3bl4o/v1

Nisha Mohd Rafiq¹, Pietro De Camilli^{2,3,4,5}

¹University of Tuebingen;

²Department of Neuroscience, Yale University School of Medicine, New Haven, Connecticut 06510, USA.;

³Department of Cell biology, Yale University School of Medicine, New Haven, Connecticut 06510, USA.;

⁴Program in Cellular Neuroscience, Neurodegeneration and Repair.;

⁵Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815, USA.



Berrak Ugur

Yale University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpvokx3bl4o/v1>

Protocol Citation: Nisha Mohd Rafiq, Pietro De Camilli 2024. Cell culture, transfection, immunocytochemistry, and imaging. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5qpvokx3bl4o/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 12, 2024

Last Modified: May 31, 2024

Protocol Integer ID: 96916

Keywords: ASAPCRN, imaging of rpe1, cell culture, imaging, imaging of ipsc, transfection, immunocytochemistry, i3 neuron, cell, neuron, da neuron, rpe1, protocol

Funders Acknowledgements:

ASAP

Grant ID: ASAP-000580

Abstract

This protocol describes the maintenance, transfection, immunocytochemistry, and imaging of RPE1 and also transfection, immunocytochemistry, and imaging of iPSCs, i³ Neurons and DA neurons.

Attachments





Rafiq_transfection.i...

18KB

Materials


Reagents:

-  Lipofectamine™ 2000 Transfection Reagent Thermo Fisher Scientific Catalog #11668019
-  Lipofectamine™ Stem Transfection Reagent Thermo Fisher Scientific Catalog #STEM00008

Troubleshooting



General cell culture for RPE1








- 1 Grow hTERT-RPE1 cells in DMEM/F12 (Thermo Fisher Scientific) supplemented with 10% FBS (Thermo Fisher Scientific), 1% glutaMAX, 1% penicillin-streptomycin. Keep cells at  37 °C with 5% CO₂ in an enclosed incubator.

Note

For general maintenance, when cells reach 80-90% confluency, detach them from the dish with Trypsin and dilute 1:10-20 in a new dish.

Cell transfection for RPE1

5d





- 2 For live-cell imaging experiments, seed the cells on glass-bottom dishes (MatTek; 35mm) at the concentrations ranging from $1-2 \times 10^5$ cells.
- 3 For RPE1: allow the cells to adhere for  08:00:00 -  24:00:00 before being transiently transfected using  4 μ L Lipofectamine™ 2000 Transfection Reagent (Invitrogen) in Opti-MEM media, mix them with the respective plasmids ( 1 μ g -  2 μ g) and visualize after  48:00:00 .
- 4 For cilia generation, serum-starve the cells in DMEM/F12 media (without FBS) for  48:00:00 .

3d

2d

Cell transfection for iPSCs, i³Neurons and DA Neurons

2d

- 5 For live-cell imaging experiments, seed cells on glass-bottom dishes (MatTek; 35mm) at the concentrations ranging from $3-5 \times 10^5$ cells.
- 6 Use  2.4 μ L of Lipofectamine™ Stem Transfection Reagent (Invitrogen) in Opti-MEM media with respective plasmids ( 1 μ g -  3 μ g) and visualize at least  48:00:00 later.

2d

Immunocytochemistry of RPE1, iPSCs, i³Neurons and DA Neurons

2h 7m



7 For fixed-cell imaging experiments, seed cells on glass-bottom dishes (MatTek; 35mm) at the concentrations ranging from $1-3 \times 10^5$ cells.

20m

- For immunofluorescence visualization, fix the cells with 4% (v/v) paraformaldehyde (Electron Microscopy Sciences) in 1× phosphate-buffered saline (PBS) for

00:20:00 .

8 Wash cells thrice in PBS.



9 Perform cell extraction with 0.25-0.5% (v/v) Triton X-100 in PBS for 00:10:00 .

10m

10 Wash cells thrice in PBS.

11 For removal of free aldehyde groups, quench the cells with fresh 1 μ L sodium borohydride (Sigma-Aldrich) in PBS for 00:07:00 , and then wash thrice in PBS.

7m

12 Further block the cells for 00:30:00 in 5% bovine serum albumin (BSA, Sigma-Aldrich) in PBS and then incubate Overnight at 4 °C with the respective antibodies listed in Table S1.

30m



13 Wash the cells with PBS thrice the following day and incubate with Alexa Fluor-conjugated secondary antibodies (Thermo Fisher Scientific) for 01:00:00 at Room temperature .


1h

14 Wash the cells thrice in 1×PBS.

15 Use DAPI (Thermo Fisher Scientific) for nuclear staining, when necessary.

Imaging



- 16 For live imaging, maintain cells in a caged incubator with humidified atmosphere (5% CO₂) at  37 °C .

Note

The Yokogawa spinning disk field scanning confocal system with microlensing (CSU-W1 SoRa, Nikon) controlled by NIS elements (Nikon) software was used for imaging. Excitation wavelengths between 405-640 nm, CFI SR Plan ApoIR 60XC WI objective lens and SoRa lens-switched light path at 1x, 2.8x or 4x were used. SoRa images were deconvolved using the Batch Deconvolution (Nikon) software.