


May 13, 2020

## Cell-free 3PGA energy solution

DOI

[dx.doi.org/10.17504/protocols.io.8pshvne](https://dx.doi.org/10.17504/protocols.io.8pshvne)



Nadanai Laohakunakorn<sup>1</sup>

<sup>1</sup>University of Edinburgh



Nadanai Laohakunakorn

University of Edinburgh

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.8pshvne](https://dx.doi.org/10.17504/protocols.io.8pshvne)

**Protocol Citation:** Nadanai Laohakunakorn 2020. Cell-free 3PGA energy solution. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.8pshvne>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 24, 2019

**Last Modified:** May 13, 2020

**Protocol Integer ID:** 29138

**Keywords:** cell-free protein synthesis, synthetic biology, in vitro transcription translation



## Abstract

Energy solution for E. coli lysate based on 3PGA. Adapted from Sun 2013 and Cai 2015. Successfully implemented at the University of Edinburgh by Nadeen Laohakunakorn, LBNC-EPFL by Zoe Swank.

### CITATION

Sun ZZ, Hayes CA, Shin J, Caschera F, Murray RM, Noireaux V (2013). Protocols for implementing an Escherichia coli based TX-TL cell-free expression system for synthetic biology.. Journal of visualized experiments : JoVE.

LINK

<https://doi.org/10.3791/50762>

### CITATION

Cai Q, Hanson JA, Steiner AR, Tran C, Masikat MR, Chen R, Zawada JF, Sato AK, Hallam TJ, Yin G (2015). A simplified and robust protocol for immunoglobulin expression in Escherichia coli cell-free protein synthesis systems.. Biotechnology progress.

LINK

<https://doi.org/10.1002/btpr.2082>



## Materials

- Amino acids LAA21-1KT Sigma
  - Mg-glutamate 49605-250G Sigma
  - K-glutamate 49601-500G Sigma
  - DTT 10708984001 Sigma
  - NTP set R1481 ThermoFisher
  - tRNA 10109541001 Sigma
  - CoA C4282-10MG Sigma
  - NAD 10127981001 Sigma
  - cAMP A9501-1G Sigma
  - folinic acid PHR1541-1G Sigma
  - spermidine S2626-1G Sigma
  - PEG-8000 89510 Sigma
  - 3PGA P8877-1G Sigma
  - HEPES H3375-100G Sigma
  - Tris base T1503-100G Sigma
- 
- KOH
  - mass balance





## Amino acids stock solution

### 1 Make stock solution of amino acids, 1000 uL at 50 mM.

- 1.1 Weigh each amino acid (excluding tyrosine) on parafilm paper, record the weight, and carefully add together in one tube. **Alternative: carrying out at 10x quantities makes weighing much easier.**






Amino acid	weight (mg)	added (mg)
Alanine	4.5	
Arginine	8.7	
Asparagine	6.6	
Aspartate	6.7	
Cysteine	6.1	
Glutamate	7.3	
Glutamine	7.3	
Glycine	3.8	
Histidine	7.8	
Isoleucine	6.6	
Leucine	6.6	
Lysine	9.1	
Methionine	7.5	
Phenylalanine	8.3	
Proline	5.8	

<b>Serine</b>	5.3	
<b>Threonine</b>	6.0	
<b>Tryptophan</b>	10.2	
<b>Valine</b>	5.9	

- 1.2 Add  1000  $\mu\text{L}$  of deionized water to the tube to make a 50 mM stock solution, vortex to mix, and adjust pH with KOH to ~5.2 (approximately 80  $\mu\text{L}$  of 1M KOH and 920  $\mu\text{L}$  dH<sub>2</sub>O required). pH can be measured roughly by spotting 1-2  $\mu\text{L}$  of the solution on appropriate pH paper. If powder does not fully dissolve, add up to ~50  $\mu\text{L}$  more of 1M KOH. pH will be ~8. Keep solution  On ice

- 1.3 Weigh tyrosine and add to a separate tube

<b>Amino acid</b>	<b>weight (mg)</b>	<b>added (mg)</b>
<b>Tyrosine</b>	9.1	

- 1.4 Add  900  $\mu\text{L}$  of [M] 1 millimolar (mM) KOH, and dissolve as far as possible. The tyrosine powder will not be entirely soluble.
- 1.5 Add  50  $\mu\text{L}$  of 15% KOH, which should fully dissolve the powder. Then add  50  $\mu\text{L}$  of deionized water. Vortex well, and measure pH, which should be around ~11-12. Keep solution  On ice
- 1.6 Keep the tyrosine and the rest of the amino acids separate. If storage is required keep at  -80  $^{\circ}\text{C}$  (flash-freezing with liquid nitrogen is optional)

## Other components stock solution

### 2 Make stock solution of other components

## 2.1 Prepare 2M stock solution of tris base.

Com pone nt	Mass to add (g)	Wate r to add (mL)	Final conc entra tion
Tris base	60.57	250	2 M

## 2.2 Prepare 1M KOH stock, and 15% KOH stock.

## 2.3 Weigh out and make the following stock solutions. Four species require titration; their pH can be approximately measured by spotting 1 uL on appropriate pH paper.

Com pone nt	Mass to add (g)	Wate r to add (uL)	Tris to add (uL)	Final conc entra tion	Note s
L- gluta mic acid mono potas sium salt (K- gluta mate)	1.219	1000		6 M	
L- gluta mic acid hemi magn esium salt (Mg- gluta mate)	0.388 6	1000		1 M	
DTT	0.154 3	1000		1 M	
CoA	0.049 8	1000		65 mM	
NAD	0.1161	to 1000	~90	175 mM	pH 7.5-8, titrate with

					2M tris
Folinc acid	0.0160	1000		33.9 mM	
cAMP	0.2139	to 1000	~365	650 mM	pH 8, titrate with 2M tris
3-PGA	0.2604	to 1000	~540	1.4 M	pH 7.5, titrate with 2M tris
HEPES	0.4766	to 1000		2 M	pH 8, titrate with 1M KOH (around 25 uL required)
	<b>Volume to add (uL)</b>	<b>Water to add (uL)</b>		<b>Final concentration</b>	
spermidine	23.55	126.45		1 M	heat up stock solution in your hand
PEG-8000	50	50		50%	

## Energy solution

### 3 Preparation of final energy solution (4x)

This energy solution will form 25% of the final reaction volume.

#### 3.1

Add the components together to produce the final energy solution, in the following order (not critical), vortexing after the addition of each one, and keeping tube On ice

	Com pone nt	Stock (mM)	Final conc (mM)	Volu me to add (uL)
	HEPE S	2000	200	100
	Water			114.2
	ATP	100	6	<b>60</b>
	GTP	100	6	<b>60</b>
	CTP	100	3.6	<b>36</b>
	UTP	100	3.6	36
	tRNA (in mg/ml)	43.75	0.8	18.29
	CoA	65	1.04	16
	NAD	175	1.32	7.54
	cAMP	650	3	4.62
	Folini c acid	33.9	0.27	8.02
	Sper midin e	1000	4	4
	3- PGA	1400	120	85.7
	amin o acids	50	6	120
	tyrosi ne	50	3	60
	PEG- 8000	50%	8%	160
	<b>Mg- gluta mate</b>	<b>1000</b>	<b>42</b>	<b>42</b>
	<b>K- gluta mate</b>	<b>6000</b>	<b>400</b>	<b>66.6 7</b>
	DTT	1000	1	1
	<b>Total</b>			<b>1000</b>

- 3.2 Measure and record pH of final solution using pH paper (should be ~8). Aliquot into storage tubes (25 uL recommended) and (optionally) flash freeze in liquid nitrogen.





3.3 Store at -80 °C

3.4 It is possible to calibrate the energy solution for maximum yield, in which case the optimization proceeds sequentially by determining optimum concentrations for Mg-glutamate, then K-glutamate. More details are given in Sun 2013.

- Sun 2013 report final optimal concentrations of 4.5-10.5 mM Mg-glutamate, 40-160 mM K-glutamate
- Kwon and Jewett 2015 report 12 mM Mg-glutamate, 130 mM K-glutamate

This step is not necessary if all that is required is functional extract.

For optimisation of energy solution, add all components apart from PEG, Mg-glu, K-glu.

This makes a solution of 731.33 µL and can be aliquoted into 10 tubes of

73.13 µL .

## Citations

Sun ZZ, Hayes CA, Shin J, Caschera F, Murray RM, Noireaux V. Protocols for implementing an Escherichia coli based TX-TL cell-free expression system for synthetic biology.

<https://doi.org/10.3791/50762>

Cai Q, Hanson JA, Steiner AR, Tran C, Masikat MR, Chen R, Zawada JF, Sato AK, Hallam TJ, Yin G. A simplified and robust protocol for immunoglobulin expression in Escherichia coli cell-free protein synthesis systems.

<https://doi.org/10.1002/btpr.2082>