

May 23, 2023

Cell lysis and immunoblotting for protein and phospho-protein quantification

 [Nature Communications](#)

DOI

dx.doi.org/10.17504/protocols.io.5jyl8j5zrg2w/v1

Dan Dou^{1,2}, C. Alexander Boecker³, Erika L.F. Holzbaur^{1,2}

¹Department of Physiology, Perelman School of Medicine, University of Pennsylvania, Philadelphia, PA 19104, USA;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, USA;

³Department of Neurology, University Medical Center Goettingen, 37077 Goettingen, Germany



Dan Dou

University of Pennsylvania

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl8j5zrg2w/v1>

External link: <https://doi.org/10.1038/s41467-025-62379-5>



Protocol Citation: Dan Dou, C. Alexander Boecker, Erika L.F. Holzbaur 2023. Cell lysis and immunoblotting for protein and phospho-protein quantification. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5jyl8j5zrg2w/v1>

Manuscript citation:

Basak, B., Holzbaur, E.L.F. Mitochondrial damage triggers the concerted degradation of negative regulators of neuronal autophagy. *Nat Commun* **16**, 7367 (2025). <https://doi.org/10.1038/s41467-025-62379-5>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 13, 2022

Last Modified: May 31, 2024

Protocol Integer ID: 71304

Keywords: Lysate, Immunoblot, Western, Phospho, LRRK2, ASAPCRN, mouse embryonic fibroblast, embryonic fibroblast, derived neuron, levels of protein, protein quantification, protein, human ipsc, neuron, using western blot

Funders Acknowledgements:

ASAP

Grant ID: ASAP-000350

Abstract

Here, we describe the procedure by which human iPSC-derived neurons or mouse embryonic fibroblasts (MEFs) were lysed and probed for levels of proteins of interest using Western blot.

Attachments



[552-1148.pdf](#)

77KB






Materials

Reagents

▪ RIPA buffer

	A	B
	Tris-HCl	50 mM
	NaCl	150 mM
	Triton X-100	0.1%
	Deoxycholate	0.5%
	SDS	0.1%

-  HALT phosphatase and protease inhibitor cocktail (100x) **Thermo Fisher Scientific Catalog #78442**
-  Microcystin-LR Microcystis aeruginosa - CAS 101043-37-2 - Calbiochem **Merck MilliporeSigma (Sigma-Aldrich) Catalog #475815m**
-  Pierce BCA Protein Assay Kit **Thermo Fisher Scientific Catalog #23225**

▪ 4x Protein Loading Buffer

	A	B
	Tris-HCl, pH 6.8	125 mM
	Glycerol	50%
	SDS	4%
	Orange G	0.2%

- Acrylamide
- **4x Running buffer**






	A	B
	Trizma base	48 g
	Glycine	230.4 g
	NaN ₃	20 mL
	ddH ₂ O	Diluted to 4 L

▪ Running buffer

	A	B
	4x running buffer	250 mL
	ddH ₂ O	750 mL
	10% SDS	10 mL

▪ Transfer buffer

	A	B
	4x running buffer	125 mL
	ddH ₂ O	875 mL
	10% SDS	500 µL
	For RABs add 20% Methanol	


-  Immobilon®-FL PVDF Membrane **Merck MilliporeSigma (Sigma-Aldrich) Catalog #Ipfl00010**
-  Chameleon® Duo Pre-stained Protein Ladder **LI-COR Catalog #928-60000**
-  Revert™ 700 Total Protein Stain for Western Blot Normalization (250 ml) **LI-COR Catalog #926-11021**

▪ Revert Wash Solution

	A	B
	Acetic acid	6.7%
	Methanol	30%
	in ddH ₂ O	

▪ Revert Reversal Solution

	A	B
	NaOH	0.1 M
	Methanol	30%
	in ddH ₂ O	

-  EveryBlot Blocking Buffer 500 ml **Bio-Rad Laboratories Catalog #12010020**
- Primary antibodies (see Materials and Methods for specific antibodies used)
- Secondary antibodies (see Materials and Methods for specific antibodies used)

Equipment

- ODYSSEY CLx Imaging System (LI-COR)

Equipment

Mini-PROTEAN Tetra Vertical Electrophoresis Cell

NAME

Electrophoresis Cell

TYPE

Bio-Rad

BRAND

1658004

SKU

<https://www.bio-rad.com/en-in/product/mini-protean-tetra-vertical-electrophoresis-cell?ID=N3F2UD4VY>

LINK



Equipment

Mini Trans-Blot Electrophoretic Transfer Cell

NAME

Electrophoretic Transfer Cell

TYPE

Bio-Rad

BRAND

1703930

SKU

<https://www.bio-rad.com/en-in/sku/1703930-mini-trans-blot-electrophoretic-transfer-cell?ID=1703930>^{LINK}

Troubleshooting




Safety warnings






- Microcystin-LR is an extremely potent hepatotoxin and should be handled with great care.
- Acrylamide is a neurotoxin and should be handled with care.
- Methanol-containing reagents should be handled carefully, as methanol can penetrate single-layer laboratory gloves.



Preparation of cell lysates


- 1 Quickly wash cells twice with ice-cold PBS. After the second wash, tilt the dish and completely aspirate all residual PBS. 
- 2 Immediately add ice-cold lysis buffer, ensuring that the entire surface is covered by lysis buffer. Place cells  On ice . 

Note

The amount of lysis buffer to use depends on cell confluency / cell number, cell type, and cell culture dish. In most cases, using  100 μL -  150 μL lysis buffer per well of a 6-well plate should result in a protein concentration >  1 $\mu\text{g}/\mu\text{L}$.

Note

Halt protease and phosphatase inhibitor cocktail and microcystin-LR should be added fresh on the day of use.

- 3 Scrape cells off the dish using a cell lifter.
- 3.1 Transfer the cell lysate to an Eppendorf tube  On ice .




3.2

20m

Leave lysates  On ice for  00:20:00 to allow for efficient lysis.

Note

Cell lysates can be snap frozen in liquid nitrogen and stored at  -80 °C for future use.

4

Centrifuge for  00:10:00 at  17000 x g and  4 °C .

10m



4.1

Discard pellet and use clarified supernatant to determine protein concentration by BCA assay following the manufacturer's instructions, performing all measurements in triplicates.



5

Add  100 µL β-mercaptoethanol to  900 µL of 4x Protein Loading Buffer and mix well.



5.1

Add complete 4x Protein Loading Buffer to cell lysates, mix well, and boil for

 00:05:00 at  95 °C .

5m

**Note**

Do not store Protein Loading Buffer with BME for more than two weeks.

SDS-polyacrylamide gel electrophoresis

6

Load samples onto 8% (for LRRK2 protein) to 15% (for PPM1H and Rab proteins) acrylamide gels alongside Chameleon Duo pre-stained protein ladder (LI-COR).





**Note**

Carefully rinse wells with running buffer before loading cell lysates.

- 7 Start electrophoresis at 80 V for  00:20:00 , then increase to 120 V and electrophorese until orange dye runs out.

20m

Protein transfer

- 8 Activate Immobilon-FL PVDF membrane by submerging in methanol for  00:00:30 -  00:01:00 .

1m 30s

- 8.1 Wash in ddH₂O and equilibrate in transfer buffer.



- 8.2 Soak sponges in methanol, wash in ddH₂O and equilibrate in transfer buffer.




- 8.3 Equilibrate filter paper in transfer buffer.

- 8.4 Assemble blotting sandwich.

- 8.5 Carefully remove any air bubbles between layers using a roller.

- 9 Fill transfer tank with ice-cold transfer buffer.



- 9.1 Place transfer system  On ice .

- 9.2 Transfer proteins from gel onto PVDV membrane at 100 V for  01:15:00 .

1h 15m






Total protein stain, membrane blocking, and antibody incubation


1h 41m



- 10 After wet-tank transfer, let membrane dry completely for at least 01:00:00 at Room temperature .
- 11 Rehydrate membrane for 00:01:00 in 100% methanol, then wash 00:05:00 in 1x TBS.
- 12 Incubate for 00:05:00 in Revert total protein stain (LI-COR) while gently shaking at Room temperature .
- 12.1 Wash twice with Revert wash solution, then rinse in ddH₂O and image membrane on ODYSSEY CLx imaging system.
- 13 Remove Revert total protein stain by incubating membrane in Revert reversal solution for 00:10:00 at Room temperature while gently shaking.
- 13.1 Rinse membrane with ddH₂O.
- 14 Block membrane for 00:05:00 in Everyblot Blocking Buffer (Bio-Rad) at Room temperature .
- 15 Dilute primary antibodies in Everyblot Blocking Buffer and incubate at 4 °C Overnight .
- 16 Wash membrane in 1x TBS + 0.1% Tween-20 (TBS-T) at Room temperature (4 washes, 00:05:00 each).
- 17 Dilute secondary antibodies 1:20,000 in Everyblot Blocking Buffer and 0.02% SDS.



Note

Incubate membrane in secondary antibodies for  01:00:00 at

 Room temperature .

18 Wash membrane in TBS-T at  Room temperature (4 washes,  00:05:00 each).

5m



19 Rinse membrane with TBS (no detergent), then image on ODYSSEY CLx imaging system. Quantify signal intensity using Image Studio Software.

