

Nov 02, 2019

## Cell lysis, detergent-free



Forked from [Cell lysis, detergent-free](#)

DOI

[dx.doi.org/10.17504/protocols.io.y4ffyt9](https://dx.doi.org/10.17504/protocols.io.y4ffyt9)

Teesha C Luehr<sup>1</sup>

<sup>1</sup>University of British Columbia

Leonard Foster's Lab



Teesha C Luehr

University of British Columbia

---

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.y4ffyt9](https://dx.doi.org/10.17504/protocols.io.y4ffyt9)

**Protocol Citation:** Teesha C Luehr 2019. Cell lysis, detergent-free. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.y4ffyt9>

**Manuscript citation:**  
<https://doi.org/10.1074/mcp.M114.047407>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** March 13, 2019

**Last Modified:** November 02, 2019

**Protocol Integer ID:** 21351



## Abstract

Detergents are generally not compatible with mass spectrometers, so this is a detergent-free method of cell lysis that is compatible with mass spectrometry. Since this protocol does not have a precipitation step, it saves time and minimizes sample loss as well.



## Growing Cells

- 1 Culture HeLa cells in Dulbecco's modified Eagle's medium (DMEM) containing 10% FBS, 20mM glutamine, and 1% PenStrep.
- 2 Original protocol uses  $5 \times 10^7$  cells per sample. Use  $2 \times 10^7$  cells. Use one 15cm Petri dish or 2-3 10 cm dishes. Each 10cm dish gives  $8.8 \times 10^6$  cells.

## Preparing Solutions and Materials Needed

- 3 Ice Bucket

- 4 [M] 2 Molarity (M) tris at pH 8.5

Add 1.211 g tris-base and 243  $\mu$ L concentrated HCl (37%) up to a total of 5 mL in Milli-Q water

### Note

Calculations done with <https://www.cytographica.com/lab/HHTris.html>.

### Note

Tris-base is in NCE 438 chemical room on the top right doubles shelf

### Note

A stock conical tube of [M] 2 Molarity (M) tris is stored on the bench.

- 5 [M] 100 millimolar (mM) TCEP

Add 250.187 mg tris(2-carboethyl)phosphine to 10 mL Milli-Q water

**Note**

TCEP is located in NCE 436 -20C freezer **door** (currently the top shelf)

**Note**

🧪 405  $\mu$ L aliquots of [M] 100 millimolar (mM) TCEP have been made and are stored in Teesha's 🧊 -20 °C storage box

**6**

[M] 400 Molarity (M) CAA

Add 🧪 374.04 mg 2-chloroacetamide to 🧪 10 mL Milli-Q water

**Note**

CAA is located in NCE 438 chemical room on the "C" shelf

**Note**

🧪 405  $\mu$ L aliquots of [M] 400 millimolar (mM) CAA have been made and are stored in Teesha's 🧊 -20 °C storage box

**7**

[M] 50 millimolar (mM)  $\text{NH}_4\text{HCO}_3$

Add 🧪 197.4 mg ammonium bicarbonate to 🧪 50 mL Milli-Q water

**Note**



Upscale or downscale volume needed. Up to 15 mL will be used for each replicate.

**8**


🧊 95 °C heating block




## Pellet Cells

- 9 Using a  15 mL conical tube pellet cells for  00:05:00 at 300 g

### Note

Pellet in a tube a minimum size of  15 mL as this tube will be used through till the end of digestion.


- 10 Wash cells with  10 mL cold PBS

- 11 Pellet cells for  00:05:00 at 300 g


- 12 Carefully discard supernatant

- 13 Store cell pellet on ice


### Note

If not performing cell lysis immediately, the pellet can be stored at  -80 °C until further use.


## Cell Lysis

- 14 Resuspend the cell pellet(s) in  1.5 mL ice cold Milli-Q water

### Note

Perform the lysis and digestion in the  15 mL conical tube . The lysis volumes, sonification, and digestion require the larger tube volume.



15 Add  1.5 mL trifluoroethanol

#### Safety information

This step should be done in the fume hood.

#### Note

1:1 water-TFE acts as a hypotonic aqueous buffer to lyse cells, eliminating the need for detergent. TFE helps protein solubility and denaturation; it readily evaporates, so removing it is easy.

#### Note


TFE is located in the NCE 435 flammable cabinet.

#### Note

TFE evaporates fast, so work quickly.

16 Cool for  00:10:00 on ice

17 Mix the sample for  00:01:00 with a vortex












18 Sonicate the sample for  00:02:00 with the Branson Digital Sonifier 250 at 30% amplitude in pulse mode (pulse ON for 0.2s and pulse OFF for 0.8s) with the tapered micro tip probe.

#### Note

The Branson Sonifier is the cell disrupter of choice, however, the bench top water bath sonicator can be used for 10 minutes.



## Reduction & Alkylation

- 19 Add  200  $\mu$ L 2M tris for a final concentration of  100 millimolar (mM)
- 20 Mix the sample for  00:00:05 with a vortex
- 21 Add  400  $\mu$ L 100 mM TCEP for a final concentration of  10 millimolar (mM)
- 22 Mix the sample for  00:00:05 with a vortex
- 23 Add  400  $\mu$ L 400 mM CAA for a final concentration of  40 millimolar (mM)
- 24 Mix the sample for  00:00:05 with a vortex
- 25 Incubate in the heating block for  00:10:00 at  95 °C

## LysC Digestion

- 26 Dilute the sample to a total of  15 mL 50mM NaHCO<sub>3</sub>

### Note

This is performed to dilute the TFE.

- 27 Measure the protein concentration with a NanoDrop (using SCOPES A205 Protein)
- 28 Calculate the total amount of protein that is desired to carry forward with the experiment. Keeping in mind that the final peptide concentration will be approximately 10-50% of the



protein concentration at this step. Transfer this volume to a new tube. If small enough, transfer to a 2 mL lo-bind tube .

29

Calculate how much LysC is needed for a 1 µg LysC : 100 µg protein

30

Add the calculated amount of LysC to the sample

#### Note

LysC is in the -80 °C freezer

31

Incubate for 02:00:00 at 37 °C in the digestion incubator

## Trypsin Digestion

32

Calculate how much trypsin is needed for a 1 µg trypsin : 100 µg protein

33

Add the calculated amount of trypsin to the sample

#### Note

Trypsin is in NCE 435 -20 freezer in the door, bottom shelf

34

Incubate for a minimum of 16:00:00 at 37 °C in the digestion incubator

## Second Trypsin Digestion

35

Calculate how much trypsin is needed for a 1 µg trypsin : 100 µg protein



#### Note

A second trypsin digestion is optional. User digression is advised.





36 Add the trypsin to the sample

37 Incubate for  05:00:00 at  37 °C in the digestion incubator

## STAGE Tip

---

38 Proceed to the STAGE tip protocol