



Cell Surface Immunofluorescent Staining of Whole Blood V.3



DOI

dx.doi.org/10.17504/protocols.io.tjwekpe

Kelsey Knight¹

¹BioLegend

BioLegend

Tech. support email: tech@biolegend.com



Sam Li

BioLegend

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account



DOI: <https://dx.doi.org/10.17504/protocols.io.tjwekpe>

External link:

http://www.biolegend.com/media_assets/support_protocol/BioLegend_Surface_Staining_Flow_Protocol_060215.pdf

Protocol Citation: Kelsey Knight . Cell Surface Immunofluorescent Staining of Whole Blood. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.tjwekpe>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Created: September 13, 2018

Last Modified: August 29, 2019

Protocol Integer ID: 15702

Keywords: cell surface immunofluorescent, whole blood, staining, cell, blood

Guidelines

Reagent List:

- Cell Staining Buffer (BioLegend Cat. No. 420201)
- Red Cell Lysis Buffer (BioLegend Cat. No. 420301)
- 7-AAD Viability Staining Solution (BioLegend Cat. No. 420403)
- TruStain FcX™ (anti-CD16/32, BioLegend Cat. No. 101319)
- Human TruStain FcX™ (Fc Receptor Blocking Solution, BioLegend Cat. No. 422301)

References:

Current Protocols in Cytometry (John Wiley & Sons, New York), Unit 6 Phenotypic Analysis.

Materials

MATERIALS

⊗ Cell Staining Buffer **BioLegend Catalog #420201**

⊗ Red Cell Lysis Buffer **BioLegend Catalog #420301**

⊗ 7-AAD Viability Staining Solution **BioLegend Catalog #420403**

⊗ TruStain FcX™ **BioLegend Catalog #101319**

⊗ Human TruStain FcX™ **BioLegend Catalog #422301**

STEP MATERIALS

⊗ Cell Staining Buffer **BioLegend Catalog #420201**

⊗ Cell Staining Buffer **BioLegend Catalog #420201**



Protocol materials

⊗ Cell Staining Buffer **BioLegend Catalog #420201**

⊗ Cell Staining Buffer **BioLegend Catalog #420201**

⊗ Red Cell Lysis Buffer **BioLegend Catalog #420301**

⊗ 7-AAD Viability Staining Solution **BioLegend Catalog #420403**

⊗ TruStain FcX™ **BioLegend Catalog #101319**

⊗ Human TruStain FcX™ **BioLegend Catalog #422301**








⊗ Cell Staining Buffer **BioLegend Catalog #420201**


⊗ Cell Staining Buffer **BioLegend Catalog #420201**

⊗ Cell Staining Buffer **BioLegend Catalog #420201**

Troubleshooting



- 1 Add predetermined optimum concentrations of desired fluorochrome conjugated, biotinylated, or purified primary antibodies to 100 µl of anti-coagulated whole blood.
- 2 Incubate at room temperature for 15-20 minutes in the dark.
 00:20:00
- 3 Dilute 10X Red Blood Cell (RBC) Lysis Buffer (BioLegend Cat. No.**420301**) to 1X working concentration with DI water. Warm to room temperature prior to use. Add 2ml of 1X RBC lysis solution to whole blood/antibody mixture. Incubate at room temperature for 10 minutes.  00:10:00
- 4 Centrifuge at 350 X g for 5 minutes, discard the supernatant.
 00:05:00
- 5 Wash 1X with at least 2 ml of Cell Staining Buffer by centrifugation at 350 x g for 5 minutes.
 00:05:00
 Cell Staining Buffer **BioLegend Catalog #420201**
- 6 If using a purified primary antibody, resuspend pellet in residual buffer and add a previously determined optimum concentration of anti-species immunoglobulin fluorochrome conjugated secondary antibody(e.g. FITC anti-mouse Ig) and incubate in the dark for 15-20 minutes.
 00:20:00
- 7 If using a biotinylated primary antibody, resuspend cell pellet in residual buffer and add a previously determined optimum concentration of fluorochrome conjugated Streptavidin (SAv) reagent (e.g. SAv-PE, BioLegend Cat. No.**405204**) and incubate for 15-20 minutes in the dark.  00:20:00
- 8 Repeat step 5.
- 9 Resuspend cells in 0.5 ml Cell Staining Buffer or 0.5 ml 2% paraformaldehyde-PBS fixation buffer.

Tip: For gentler fixation (particularly with tandem fluors), FluoroFix™ Buffer (Cat. No.**422101**) may be used.
 Cell Staining Buffer **BioLegend Catalog #420201**
- 10 Analyze with a Flow Cytometer.



Note: If you are unable to immediately read your samples on a cytometer, keep them shielded from light and in a refrigerator set at 4-8°C. The samples should be resuspended in Cell Staining Buffer. Note that samples should not remain in a fixation buffer for extended periods of time as this can affect fluor conformation and fluorescence.