

Apr 06, 2023

Cell viability assessment

DOI

<https://dx.doi.org/10.17504/protocols.io.eq2ly7rpelx9/v1>

Qing Wang^{1,2}

¹Massachusetts General Hospital; ²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network



Jackson Schumacher

Massachusetts General Hospital

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2ly7rpelx9/v1>

Protocol Citation: Qing Wang 2023. Cell viability assessment. protocols.io
<https://dx.doi.org/10.17504/protocols.io.eq2ly7rpelx9/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 06, 2023

Last Modified: December 12, 2024

Protocol Integer ID: 80147

Keywords: ASAPCRN, mts cell proliferation assay kit, assessing cell viability, cell viability assessment, cell viability, assay

Funders Acknowledgements:

From cancer associations to altered immunity in the pathogenesis of Parkinson's disease

Grant ID: ASAP-000312

Abstract

This is the protocol for assessing cell viability using a MTS cell proliferation assay kit.

- 1 1. Cells are seeded and differentiated for 4 days; assess cell viability using an MTS cell proliferation assay kit (ab197010, Abcam).
- 2 2. Treat the differentiated cells with synthetic DOPA pheomelanin or synthetic DOPA eumelanin or vehicle PBS to determine the time and dose responses of the cells to the treatment.
- 3 3. Remove culture medium and add a 110 μ L solution containing 10 μ L MTS and 100 μ L medium to each well, incubate in humidified air containing 5% CO₂ at 37°C for 1 hour.
- 4 4. Measure absorbance at 490 nm using a microplate reader.

Protocol references

Su, W., Zhu, S., Chen, K., Yang, H., Tian, M., Fu, Q., ... & Yang, C. (2021). Overexpressed WDR3 induces the activation of Hippo pathway by interacting with GATA4 in pancreatic cancer. *Journal of Experimental & Clinical Cancer Research*, 40(1), 1-15.

Ye, Q., Srivastava, P., Al-Kuwari, N., & Chen, X. (2022). Oncogenic BRAFV600E induces microglial proliferation through extracellular signal-regulated kinase and neuronal death through c-Jun N-terminal kinase. *Neural Regeneration Research*.