



Mar 28, 2025

Cell Viability Protocol using CellTiter-Glo 3D

 [Methods and Protocols](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.5qpvo055xv4o/v1>

Marine Simonneau¹

¹LabTAU



magali perier

INSERM U1032

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpvo055xv4o/v1>

External link: <https://doi.org/10.3390/mps8040075>

Protocol Citation: Marine Simonneau 2025. Cell Viability Protocol using CellTiter-Glo 3D. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.5qpvo055xv4o/v1>

Manuscript citation:

Perier M, Wang L, Simonneau M, Ngo-Reymond J, Guillermet-Guibert J, Lafond M, Lafon C (2025) Formation and Characterization of Two Magnetic Three-Dimensional Spheroid Models of Murine Pancreatic Adenocarcinoma. *Methods and Protocols* 8(4). doi: [10.3390/mps8040075](https://doi.org/10.3390/mps8040075)



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 28, 2025

Last Modified: March 28, 2025

Protocol Integer ID: 125616

Keywords: viability measurement for 3d culture, cell viability protocol, cell viability, glo 3d this protocol, using celltiter, 3d culture, cell, glo 3d, 3d, viability measurement

Abstract

This protocol provides a viability measurement for 3D culture.

Materials

ATP (Mw = 551.14 g/mol) (A7699-1G, Sigma-Aldrich)

spheroids

Incomplete DMEM/F-12 culture medium without FBS (1%L-Glutamine, 1% Penicillin/Streptomycin) (Gibco 21331-020)

CellTiter-Glo® 3D Cell Viability Assay (Promega G968A)

Aluminium

Shaker

Opaque 96-well plate for reading (Costar 96 Flat Bottom White Polystyrene)

Luminescence reader (Tecan Infinite M Plex)



Preparation of the CellTiter-Glo 3D Solution

- Thaw the CellTiter-Glo reagent in the fridge at 4°C the day before the experiment.
- Before use, let the kit reach room temperature for 30 minutes.
- Gently mix the solution.

Preparation of ATP Solutions at Different Concentrations for the Calibration Curve

- Dilution of the ATP Solution

	ATP	0	0.5	1	2	3
	Concentration (μM)					
	ATP at 10 μM (μL)	0	10	20	40	60
	DMEM without FBS (μL)	200	190	180	160	140

- Transfer the ATP solutions into the opaque 96-well plate for luminescence

Preparation of Spheroids in the Opaque Plate

- Transfer the spheroids from the incubation plate to the opaque 96-well plate for reading.
- Place the plate on the 28-magnet



plate (Plate Holding) and aspirate the medium.

- 8 ●
Add 50 μ L of DMEM without FBS per well.

Addition of CellTiter-Glo and Incubation of the Opaque Plate

- 9 ●
Add 50 μ L of CellTiter-Glo per well
(into wells containing 50 μ L of DMEM without FBS).
- 10 ●
Cover the plate with aluminum foil to
protect the samples from light, as CellTiter-Glo is light-sensitive.
- 11 ●
Incubate the plate at room temperature
for 25 minutes.
- 12 ●
During this incubation period, shake
the plate on a vortex for 5 minutes at 2 rpm.

Luminescence Reading

- 13 ●
Set
up the device:
☐ Shaking
mode: Orbital shaking
☐ Shaking
time: 3 s
☐ Attenuation:
None
☐ Integration
time: 1000 ms
☐ Settle
time: 0 ms



○

Temperature:

22°C

14



Place the plate without the lid on the reader platform and position the plate.

15



Start the reading and save the Excel file.