



Mar 28, 2025

Cell viability using PrestoBlue HS

 [Methods and Protocols](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl4z7k2vo5/v1>

Jacqueline Ngo-Reymond¹

¹INSERM U1032



magali perier

INSERM U1032

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4z7k2vo5/v1>

External link: <https://doi.org/10.3390/mps8040075>

Protocol Citation: Jacqueline Ngo-Reymond 2025. Cell viability using PrestoBlue HS. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.3byl4z7k2vo5/v1>

Manuscript citation:

Perier M, Wang L, Simonneau M, Ngo-Reymond J, Guillermet-Guibert J, Lafond M, Lafon C (2025) Formation and Characterization of Two Magnetic Three-Dimensional Spheroid Models of Murine Pancreatic Adenocarcinoma. *Methods and Protocols* 8(4). doi: [10.3390/mps8040075](https://doi.org/10.3390/mps8040075)



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 28, 2025

Last Modified: March 28, 2025

Protocol Integer ID: 125641

Keywords: cell viability, viability of cell, prestobblue hs this protocol, using prestobblue h, cell

Abstract

This protocol measures the viability of cells cultured in 2D.

Materials

KPC A219

iMEF

DMEM/F-12 culture medium (Gibco 21331-020) with FBS , 1%

L-Glutamine and 1% Penicillin/Streptomycin

DPBS

Trypsin EDTA 0,05%

96-well plate for culture bottom F

PrestoBlue HS Cell Viability Reagent (Invitrogen P50200)

Aluminium

Opaque 96-well plate for reading (Perkin Elmer Optiplate

96-F black Polystyrene 6005270)

Fluorescence reader (Horiba Fluoromax-4)

1. Day one : Cell seeding

- 1 Cells should not be confluent the day of test. They should be in exponential phase of growth. At confluency, cells may lower their metabolic rate. So before experiment, the appropriate cell density should be determined.
- 2 - Warm the complete DMEM/F12 medium at 37°C
- 3 - Under biosafety hood, aspirate the medium in 75cm² flask
- 4 -Add 5 ml of DPBS to rinse the cell and aspirate the PBS
- 5 - Add 2 ml of trypsin EDTA 0,05% and incubate at 37°C for 5 minutes. Verify that cells were detached from the flask wall. Add 10 ml of medium and flush 3 times to separate cells
- 6 -Take 20 µl of cells of each cell line and count
- 7 -Prepare cell suspension to have 4000 cells in 100 µl for drug incubation for 24 hours and 2000 cells in 100µl for 72 hours
- 8 - Dispense 100 µl in each well of 96-well plate
- 9 - Place the 96-well plates in incubator at 37°C, 5% CO₂

2. Day 2 : Drug treatment

- 10 - prepare solution to have 0-50-100-200 nM of Gemcitabine
- 11 - Aspirate the medium in 96-well plate



12 - Dispense 100µl of different solutions

13 -Incubate for 24 to 72 hours at 37°C, 5% CO₂

3. Day 3-5 : PrestoBlue Test

14 Warm PrestoBlue HS Cell Viability Reagent to room temperature

15 Aspirate the medium in each well

16 Add 90 µl of fresh medium in each well

17 Add 10 µl per well of PrestoBlue HS Cell Viability Reagent

18 Place the 96-well plate in incubator at 37°C, 5% CO₂ for 1 hour

19 Transfer medium of each well in a opaque black 96-well plate

20 Record fluorescence on plate reader (Enter value of emission (560nm) and excitation (590nm) wavelength)