

Jul 28, 2023

Cellular lipid uptake with flow cytometry readout

DOI

<https://dx.doi.org/10.17504/protocols.io.14egn2pnpg5d/v1>

Igor Beletchi¹, rosanne.wouters¹, Peter Vangheluwe¹

¹KU Leuven



rosanne.wouters

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn2pnpg5d/v1>

Protocol Citation: Igor Beletchi, rosanne.wouters , Peter Vangheluwe 2023. Cellular lipid uptake with flow cytometry readout. protocols.io <https://dx.doi.org/10.17504/protocols.io.14egn2pnpg5d/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 31, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 82703

Keywords: ASAPCRN, cellular lipid uptake with flow cytometry, cellular lipid uptake, flow cytometry, lipid

Funders Acknowledgements:

ASAP (Aligning Science Across Parkinson's)

Abstract

Cellular lipid uptake with flow cytometry readout

Protocol materials

- ✕ TrypLE Express or preferred cell dissociation reagent **Thermo Fisher Scientific**
- ✕ PBS without Ca² or Mg² **Gibco - Thermo Fisher Scientific Catalog #10010-031**
- ✕ HBSS **Gibco - Thermo Fisher Scientific Catalog #14170-112**
- ✕ NBD-SM (C6 NBD sphingomyelin) **Avanti Polar Lipids, Inc. Catalog #810218C**
- ✕ NBD-GluCer (C6-NBD Glucosyl Ceramide) **Avanti Polar Lipids, Inc. Catalog #810222C**
- ✕ NBD-PC (18:1-06:0 NBD PC) **Avanti Polar Lipids, Inc. Catalog #810132C**
- ✕ NBD-lyso-PC (12:0 lyso NBD PC) **Avanti Polar Lipids, Inc. Catalog #810128C**
- ✕ NBD-PS (18:1-06:0 NBD PS) **Avanti Polar Lipids, Inc. Catalog #810194C**
- ✕ NBD-PE (18:1-06:0 NBD PE) **Avanti Polar Lipids, Inc. Catalog #8105155C**
- ✕ Albumin Bovine Serum Fraction V Fatty Acid-Free **Merck Catalog #126575**



harvest cells

10m

1 harvest cells by detachment with

 TrypLE Express or preferred cell dissociation reagent **Thermo Fisher Scientific**

1.1 cells should be grown to 70% confluency, not higher



1.2 wash cells with

 TrypLE Express or preferred cell dissociation reagent **Thermo Fisher Scientific**

1.3 incubate cells with

 TrypLE Express or preferred cell dissociation reagent **Thermo Fisher Scientific**

(3ml/T175 flask) at room temperature

1.4 when cells detach resuspend with 7ml

 PBS without Ca²⁺ or Mg²⁺ **Gibco - Thermo Fisher Scientific Catalog #10010-031**

and keep on ice in 15 mL falcon

2 collect cells by centrifugation  400 x g, 4°C, 00:05:00

5m



3 wash pellet with with

 PBS without Ca²⁺ or Mg²⁺ **Gibco - Thermo Fisher Scientific Catalog #10010-031**

4 collect cells by centrifugation  400 x g, 4°C, 00:05:00

5m

5 resuspend pellet in  HBSS **Gibco - Thermo Fisher Scientific Catalog #14170-112**

6 count cells

7 prepare cell suspension with 1*10⁶ cells in 500µl

 HBSS **Gibco - Thermo Fisher Scientific Catalog #14170-112**



prepare lipid suspension

- 8 work under fume hood equipped with N2 flow
- 9 wash glass syringe 3x with 100% ethanol and 3x with chloroform, let dry
- 10 with glass syringe pipette needed volume of NBD-lipid into glass vial
 -  NBD-PC (18:1-06:0 NBD PC) **Avanti Polar Lipids, Inc. Catalog #810132C** ,
 -  NBD-lyso-PC (12:0 lyso NBD PC) **Avanti Polar Lipids, Inc. Catalog #810128C** ,
 -  NBD-PS (18:1-06:0 NBD PS) **Avanti Polar Lipids, Inc. Catalog #810194C** ,
 -  NBD-PE (18:1-06:0 NBD PE) **Avanti Polar Lipids, Inc. Catalog #8105155C** ,
 -  NBD-SM (C6 NBD sphingomyelin) **Avanti Polar Lipids, Inc. Catalog #810218C** ,
 -  NBD-GluCer (C6-NBD Glucosyl Ceramide) **Avanti Polar Lipids, Inc. Catalog #810222C**
- 10.1 keep glass bottle with NBD-labeled lipids in chloroform solution cold, by working on ice or in freezer block
- 10.2 for each sample 500µl lipid suspension with 2µM final concentration is needed
- 11 dry lipids by evaporating chloroform under N2 flow
- 12 resuspend lipids in  **HBSS Gibco - Thermo Fisher Scientific Catalog #14170-112** to final concentration of 2µM with magnetic stirbar until fully dissolved, keep lipids at 4°C while dissolving

lipid uptake

- 13 equilibrate cell suspension  0 rpm, 37°C, 00:15:00



- 14 equilibrate lipid suspension 🔗 0 rpm, 37°C, 00:15:00
- 15 prepare collection tubes with 200µl
 - ⊗ HBSS **Gibco - Thermo Fisher Scientific Catalog #14170-112** containing 5%
 - ⊗ Albumin Bovine Serum Fraction V Fatty Acid-Free **Merck Catalog #126575** (per sample 1 collection tube is needed for each timepoint) and keep on ice
- 16 for timepoint 0, take 200µl of equilibrated cell suspension and add to collection tube T=0 (containing 200µl of ice-cold
 - ⊗ HBSS **Gibco - Thermo Fisher Scientific Catalog #14170-112** containing 5%
 - ⊗ Albumin Bovine Serum Fraction V Fatty Acid-Free **Merck Catalog #126575**) on ice
- 17 add 500µl of lipid suspension to remaining 300µl of cell suspension
- 18 incubate at 37°C with thermoshaker with interval 1 min shaking every 5 min (700 rpm).
- 18.1 keep samples away from light during incubation
- 19 after 30min, take 200µl of lipid+cell suspension and add to collection tubes T=30 (containing 200µl of ice-cold
 - ⊗ HBSS **Gibco - Thermo Fisher Scientific Catalog #14170-112** containing 5%
 - ⊗ Albumin Bovine Serum Fraction V Fatty Acid-Free **Merck Catalog #126575**), keep on ice
- 20 after 60min, take 200µl of lipid+cell suspension and add to collection tubes T=60 (containing 200µl of ice-cold
 - ⊗ HBSS **Gibco - Thermo Fisher Scientific Catalog #14170-112** containing 5%
 - ⊗ Albumin Bovine Serum Fraction V Fatty Acid-Free **Merck Catalog #126575**), keep on ice
- 21 add sodium dithionite to all samples (1/100 dilution of freshly prepared 1M stock in tris-buffer, pH10) and vortex samples (make sure to mix well with quencher)
- 22 measure total internal fluorescence with flow cytometer

Protocol references

The lipid uptake assay was adapted from (Naito *et al.*, 2015).