

Apr 15, 2019

Centrifuge the whole blood to separate red blood cells

DOI

dx.doi.org/10.17504/protocols.io.zzwf77e

George Lykotrafitis¹

¹University of Connecticut



George Lykotrafitis

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.zzwf77e>

Protocol Citation: George Lykotrafitis 2019. Centrifuge the whole blood to separate red blood cells. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.zzwf77e>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol in our group and it is working.

Created: April 11, 2019

Last Modified: April 15, 2019

Protocol Integer ID: 22294

Keywords: red blood cell, whole blood in order, whole blood, cell

Abstract

Centrifuge the whole blood in order to obtain the red blood cells.

Materials

Human whole blood, centrifuge, centrifuge tubes.

Troubleshooting

- 1 Centrifuge the whole blood using the parameters: 500 ×g, 4 °C, and 5 min.
- 2 Remove the plasma and buffy coat.
- 3 Add the same volume of Alsever's solution as the remaining blood components.
- 4 Centrifuge again using the parameters: 500 ×g, 4 °C, and 5 min.
- 5 Remove the top layer.
- 6 Add the same volume of Alsever's solution as the remaining blood components.
- 7 Centrifuge a third time using the parameters: 500 ×g, 4 °C, and 5 min.
- 8 Remove the top layer.
- 9 Add the same volume of Alsever's solution as the remaining blood component - red blood cells, and store at 4 °C for up to 7 days.