



Jun 12, 2019

CGAP Human Oesophagus Epithelium Dissociation - Tissue Stability

 Forked from [CGAP Human Oesophagus Epithelium Dissociation](#)

DOI

dx.doi.org/10.17504/protocols.io.34fgqtn

Anna Wilbrey-Clark¹, Adam Hunter²

¹Wellcome Trust Sanger Institute, HCA; ²CGAP

Human Cell Atlas Metho...



Anna Wilbrey-Clark

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.34fgqtn>

Protocol Citation: Anna Wilbrey-Clark, Adam Hunter 2019. CGAP Human Oesophagus Epithelium Dissociation - Tissue Stability. [protocols.io](#) <https://dx.doi.org/10.17504/protocols.io.34fgqtn>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: June 12, 2019

Last Modified: June 12, 2019

Protocol Integer ID: 24423

Keywords: tissue stability

Troubleshooting



1

	Mat erial	Qua ntity	Sup plier Info
	100 mm Petri Dish	1	Corn ing (430 591)
	Cold PBS	50ml	GIB CO (141 90- 144)
	Forc eps	2	Ther moFi sher UK Ltd (152 322 90)
	Scal pels	2	Swa nn- Mort on Ltd (050 7)
	50ml Falc on Tub es	2	Falc on (352 098)
	15ml Falc on Tub es	1	Falc on (352 097)
	Diss ociat ion Age nt (Try psin - EDT A 0.25 %)	40ml	GIB CO (252 00- 056)
	DNa se I (Sto ck	400 ul	Sig ma (112 849)



	solution 10mg/ml)		320 01)
	RPMI + 20% FBS	25ml	Gibco (424 0104 2) + Sigma (F75 24- 50ML)
	70um Cell Strainer	2	ThermoFisher UK Ltd (153 708 01)
	2.0 ml Syringe	1	ThermoFisher UK Ltd. (106 7355 5)
	1.5ml DNA LoBind Eppendorf Tubes	1	Eppendorf (003 0 108. 051)
	BSA	267 µl 7.5% BSA in 50ml PBS -/-	Sigma (A84 12)

0.5ml DNA LoBind Eppendorf Tubes	2	Eppendorf (0030108.035)
Red Cell Lysis Buffer	5ml	Life Technologies Ltd. (00-4333-57)
Trypan Blue	40ul	Fisher Scientific (11414815)
C-Chips	2	Cambridge Bioscience (DH C-N01-50)

- 2 Receive oesophagus sample from mid-region in hypothermasol FRS solution (Sigma H4416).
- 3 Wash the samples with 10ml cold PBS to remove any residual contamination, stomach content and loose mucus.
- 4 Pour oesophagus onto 100mm glass petri dish and add another 10ml fresh cold PBS.
- 5 Open the samples longitudinally.
 - Epithelium/mucosa should be a relatively loose, yellowish layer on the lumen side.
- 6 Using two forceps (scissors if necessary) separate the mucosa/epithelium layer from the stroma. Place them onto separate 100mm petri dishes each with 10ml PBS (this is a

relatively crude separation protocol and is not 'pure' but avoids long overnight incubations in dispase).

- From this stage on only process the mucosa/epithelium part for dissociation.

- 7 Aspirate PBS.
- 8 In a few drops of PBS (~200µl), finely mince the epithelium using two scalpels simultaneously.
 - Too much PBS in the dish will make it more difficult to mince.
- 9 Add 19ml Trypsin-EDTA 0.25% containing 100µg/ml DNase I to the dish and transfer the tissue to a 50ml falcon tube. (DNase is critical for this protocol). Note later tests indicated that this volume can be reduced to 5ml total.
- 10 Wash the scalpel and the dish with a further 1ml of Trypsin-EDTA 0.25% and transfer it into the Falcon with the tissue.
- 11 Incubate the biopsies for 30 min at 37°C on a rocker.
- 12 Centrifuge sample 500g, 5 min at 4°C.
- 13 Add 20ml fresh 0.25% Trypsin-EDTA with 100µg/ml DNase and incubate at 37°C on rocker for a further 15min. (Note it was later determined that this second incubation step is not necessary for oesophagus, only for stomach, which causes a change in pH of the dissociation agent).
- 14 Pass through a 70µM cell strainer into a 50ml falcon tube.
- 15 Add 20ml of RPMI + 20% FBS to inactivate trypsin, mash the remaining tissue in the strainer with a syringe plunger and wash again with 5ml RPMI + 20% FBS. There is usually some undigested material remaining).
- 16 Centrifuge at 500g for 5 minutes, 4°C.
- 17 Wash cells with 5ml cold PBS-/-.
- 18 Centrifuge at 500g for 5 minutes, 4°C.



- 19 Add 5ml Red Cell Lysis buffer to the pellet and resuspend.
- 20 Place on ice 5 min with periodic agitation.
- 21 Add 10ml of Cold PBS-/-.
- 22 Filter through a 70µm cell strainer. (Often, small white filaments of undissociated material are still visible at this stage which must be removed before 10x loading).
- 23 Centrifuge at 500g for 5 minutes, 4°C. Remove supernatant.
- 24 Resuspend in 300-1000µl cold PBS with 0.04% BSA.
- 25 Count cells and viability manually using C-chip / trypan blue.
- 26 Proceed to 10X Preparation for scRNA sequencing.
- 27 Ensure all unused tissue, equipment and tubes that have been in contact with primary tissue are placed into Virkon in sweetie jar for a minimum of 1 hour. After this time aspirate and disposing in relevant sharps or waste routes.