


Jun 24, 2023

Characterizing spatial and temporal properties of ER-phagy receptors

 [Nature Cell Biology](#)

 In 3 collections

DOI

dx.doi.org/10.17504/protocols.io.6qpvr3en3vmk/v1

Melissa Hoyer¹, Harper JW¹

¹Harvard Medical School

Harper

ASAP



Melissa Hoyer

Harvard Medical School

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvr3en3vmk/v1>

External link: <https://doi.org/10.1038/s41556-024-01356-4>

Protocol Citation: Melissa Hoyer, Harper JW 2023. Characterizing spatial and temporal properties of ER-phagy receptors. [protocols.io https://dx.doi.org/10.17504/protocols.io.6qpvr3en3vmk/v1](https://dx.doi.org/10.17504/protocols.io.6qpvr3en3vmk/v1)

**Manuscript citation:**

Melissa J. Hoyer, Cristina Capitanio, Ian R. Smith, Julia C. Paoli, Anna Bieber, Yizhi Jiang, Joao A. Paulo, Miguel A. Gonzalez-Lozano, Wolfgang Baumeister, Florian Wilfling, Brenda A. Schulman, J. Wade Harper (2024) Combinatorial selective ER-phagy remodels the ER during neurogenesis. *Nature Cell Biology* doi: [10.1038/s41556-024-01356-4](https://doi.org/10.1038/s41556-024-01356-4)

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 24, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 83982

Keywords: ASAPCRN, axonal er accumulation within synaptic er bouton, phagy receptors the endoplasmic reticulum, quantitative landscape of er proteome remodeling, axonal er accumulation, specific autophagic capture of er, synaptic er bouton, embedded receptor, phagy receptor fam134b, endoplasmic reticulum, er remodeling by autophagy, extensive er remodeling during differentiation, er to degradative autophagy machinery, phagy receptor, deficient neurons in *vivodisplay*, deficient neuron, er proteome remodeling, selective autophagy, underlying er remodeling, neuron, receptor, refined tubular er network, vast proteomic landscape, monitoring extensive er remodeling, autophagy, degradative autophagy machinery, er proteome, tractable induced neuron, specific autophagic capture, temporal properties of er, receptor selectivity, selective er, proteome, axon, protein, linking er

Funders Acknowledgements:

ASAP

Grant ID: Grant ID: ASAP-000282

Abstract

The endoplasmic reticulum (ER) has a vast proteomic landscape to preform many diverse functions including protein and lipid synthesis, calcium ion flux, and inter-organelle communication. The ER proteome is remodeled in part through membrane-embedded receptors linking ER to degradative autophagy machinery (selective ER-phagy). A refined tubular ER network is formed in neurons within highly polarized dendrites and axons. Autophagy-deficient neurons *in vivodisplay* axonal ER accumulation within synaptic ER boutons, and the ER-phagy receptor FAM134B has been genetically linked with human sensory and autonomic neuropathy. However, mechanisms and receptor selectivity underlying ER remodeling by autophagy in neurons is limited. Here, we combine a genetically tractable induced neuron (iNeuron) system for monitoring extensive ER remodeling during differentiation with proteomic and computational tools to create a quantitative landscape of ER proteome remodeling via selective autophagy. Using spatial sensors or ER-phagy we demonstrate receptor-specific autophagic capture of ER in axons.

Materials

A	B	C
REAGENT or RESOURCE	SOURCE	IDENTIFIER
Cells lines		
NGN2 inducible embryonic stem cells	See protocol dx.doi.org/10.17504/protocols.io.br9em93e	CVCL_9773 (modified from this source line)
HEK293T cells	ATCC CRL-1573	CVCL_0045
Chemicals		
Dulbecco's MEM (DMEM), F12		
Phosphate Buffered Saline 1X	Corning	21-031-CV
E8 components	See protocol dx.doi.org/10.17504/protocols.io.bsacnaaw	
ND1 and ND2 components	See protocol dx.doi.org/10.17504/protocols.io.br9em93e	
Hygromycin B	Life Technologies	10687-010
pAC150- FAM134C-GFP	This paper	Addgene 201932
pAC150- TEX264-GFP	This paper	Addgene 201931
pAC150- TEX264(deltaLIR, F273A)-GFP	This paper	Addgene 201930
pCMV-hyPBase hyperactive piggyBac vector	Yusa et al 2011	Available upon request at the Sanger Institute Achives
pHAGE-TEX264-GFP	An et al 2019	Addgene 201925
pHAGE-TEX264(deltaLIR,F273A)-GFP	An et al 2019	Addgene 201926
pHAGE-mCherry-LC3B	An et al 2019	Addgene 201924
Software		
Nikon Imaging Software Elements	5.21.3 (Build 1489)	https://cellprofiler.org/ SCR_014329



	A	B	C
	Cell Profiler	CellProfiler v4.0.6	https://cellprofiler.org/ SCR_007358
	Fiji	ImageJ V.2.0.0	https://imagej.net/software/fiji/ SCR_002285
	Hardware		
	Thermo Neon™ Transfection System	Thermo Fisher Scientific	MPK5000
	Yokogawa CSU-X1 spinning disk confocal on a Nikon Ti-E inverted microscope		

Troubleshooting

Generation of Stable Cell ES H9 Ngn2 line expressing ER receptors

- 1 Electroporation of PB vectors. Use ThermoFisher kit and ThermoFisher Neon Electroporator to electroporate ES cells with PB vector and PB helper vector.
- 2 Add 10microliter buffer R to a sterile 1.5ml tube. Add 1.0mg of pAC150- TEX264-GFP, pAC150-TEX264deltaLIR-GFP, or pAC150- FAM134C-GFP vectors and 1.0mg of pCMV-hyPBBase hyperactive piggyBac vector. Pipet up and down to mix. Let it sit at RT for 10min. This is enough for 2 transfections (== one 6 well).
- 3 Individualize cells with Accutase. Neutralize Accutase with 5 times volume E8 with Rock inhibitor.
- 4 Count cells. You will need 2×10^5 for each transfection.
- 5 Spin down cells. Let it sit for a while so all the residue media can go down to the bottom of the tube. If the residue media is too much, take it out with a P200 pipet.
- 6 Resuspend cells to a concentration of 2×10^5 per 5 microliter (ie 4×10^7 per ml) using buffer R. You don't have to take all the residue media off but you will need to take into account the volume of residue media so you are not too much off.
- 7 Prepare a 6 well matrigel coated plate. Add 2mL of E8+ rock inhibitor (1:1000) to the wells you will use. Two transfections go into one well.
- 8 Wipe the Neon pipet station with EtOH and place it inside the hood.
- 9 Add 3ml of electrolytic buffer (buffer E) to the neon tube. Place the tube inside the station. You should feel a click before the tube is securely seated in the station.
- 10 Use program 13 from the optimization tab for electroporation parameter (Voltage: 1100. Pulse width: 20 Pulse number: 2). Program 9 should also work.
- 11 When everything is ready, mix 10-11microliters of resuspended cells with the plasmid containing R buffer. The final volume should be in the range of 21-22microliters.
- 12 Take up a neon tip, pipet 10microliters of the cell protein mix and electroporate with program 13. It is important to pipet slowly to avoid air bubble formation. It is also important to insert the pipet slowly into the station, especially during the end of the



insertion when you will feel a click. Help the pipet down slowly during the clicking so there is no sudden movement of the tip, which might create tiny air bubbles.

- 13 If you see air bubble in the tip, take it out, push everything out of the tip and re-pipet the mixture.
- 14 If you see sparking during the electroporation, your efficiency will reduce significantly.
- 15 Once electroporation is complete, push everything into one well of a 6 well plate. Do not pipet up and down with Neon tip.
- 16 Repeat the same procedure with the same tip and the left over cell mixture. Place the second electroporated mixture into the same well.
- 17 Disperse cells evenly in the well and place cells in a low O₂ incubator. 17. Put electroporated cells into low oxygen incubator for 2 days
- 18 Start selection of cells with 50mg/mL hygromycin B 4 days post-electroporation. Grow cells in selection medium for 7-10 days until there is no longer any cell death and every cell has integration. Cells can also be sorted for equal amount of GFP expression.

Differentiation of Stable Cell ES to induced neurons (iN) for imaging

- 19 Differentiation to induced neurons (iN) is done by following the protocol "Neural differentiation of AAVS1-TREG3-NGN2 pluripotent stem cells ([dx.doi.org/10.17504/protocols.io.br9em93e](https://doi.org/10.17504/protocols.io.br9em93e))
- 19.1 To help get spread out iNeurons for imaging, when cells get 90 percent confluent at any point in the day 5-7 range, cells are plated into onto 35 mm-glass bottom dishes (No. 1.5, 14 mm glass diameter, MatTek) at a low confluency (approximately 1×10^5 cells per 2mL dish).

Lentiviral iN transduction

- 20 LC3B tagged lentiviral pHAGE constructs (or if not using a stable piggyBac ER-phagy receptor line, co transduced with ER-phagy receptor lentiviral constructs) were packaged in HEK293T by cotransfection of pPAX2, pMD2 and the vector of interest in a 4:2:1 ratio using polyethelenimine (PEI). Normally a 6cm dish of 70-80% confluent HEK293T was transfected.
- 21 The next day the media on the transfected HEK293T cells is switched to ND2 no doxycycline media. If 6cm of HEK293T cells was transfected, 5mL of the ND2 no DOX



media is added.

- 22 The following day, virus containing supernatant (ND2 media with virus) was and filtered through a .22 micron syringe filter, aliquoted into epi tubes and frozen at -80C.
- 23 Target induced neurons are cultured to day10. After routine ND2 no doxycycline feeding cells were treated with varying amounts of the virus-containing media (for 20-40 percent iN confluent 2mL imaging dish, 300-500 microliters of ER-phagy receptor lentivirus media was added and 100-200 microliters of mCh-LC3B lentivirus media was added).

Analysis of autophagosomes and ER-phagy receptors in iN via confocal microscopy

- 24 Cells are imaged at various days in the differentiation live at 37C using a Yokogawa CSU-X1 spinning disk confocal on a Nikon Ti-E inverted microscope at the Nikon Imaging Center in Harvard Medical School. Nikon Perfect Focus System was used to maintain cell focus over time. The microscope is equipped with a Nikon Plan Apo 40x/1.30 N.A or 100x/1.40 N.A objective lens. 488nm (100mW) and 561nm (100mW) laser lines are controlled by AOTF. All images are collected with a Hamamatsu ORCA-ER cooled CCD camera (6.45 μm^2 photodiode) with Nikon Elements image acquisition software. Timelapses with various lengths and intervals were acquired.
- 25 Each channel time series are brightness and contrast adjusted equally and then converted to RGB for publication using FIJI software.

Protocol references

An H, Ordureau A, Paulo JA, Shoemaker CJ, Denic V, Harper JW. TEX264 Is an Endoplasmic Reticulum-Resident ATG8-Interacting Protein Critical for ER Remodeling during Nutrient Stress. *Mol Cell*. 2019 Jun 6;74(5):891-908.e10. doi: 10.1016/j.molcel.2019.03.034. Epub 2019 Apr 18. PMID: 31006537; PMCID: PMC6747008.