

Aug 11, 2023

ChAT Immunofluorescent Staining

DOI

dx.doi.org/10.17504/protocols.io.rm7vzxbz4gx1/v1

Haley Geertsma¹

¹University of Ottawa



Haley Geertsma

University of Ottawa

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzxbz4gx1/v1>

Protocol Citation: Haley Geertsma 2023. ChAT Immunofluorescent Staining. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.rm7vzxbz4gx1/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 11, 2023



Last Modified: August 11, 2023

Protocol Integer ID: 86350

Keywords: choline acetyltransferase, ChAT, day choline acetyltransferase, chat immunofluorescent, antibody, mouse brain, mice

Abstract

This protocol is designed for 3-day choline acetyltransferase (ChAT) staining using the EMD Millipore AB144P (RRID:AB_2079751) antibody. Tissue stained with this protocol include 40µm free-floating mouse brain and spinal cord sections, and 20µm pre-mounted retina sections. All tissue was from mice perfused with 10% buffered formalin.

Troubleshooting



ChAT Immunofluorescent Staining

2d 4h 35m

- 1 Wash tissue slices in 1X phosphate buffered saline (PBS) for 5 minutes. Repeat x2. 5m
- 2 Incubate in blocking buffer (10% horse serum, 0.5% gelatin, 0.1% triton X-100 in 1X PBS) for 2 hours at room temperature. 2h
- 3 Incubate in primary antibody (goat anti-ChAT (RRID:AB_2079751) at 1:200 dilution in blocking buffer) at 4°C for 48 hours. 2d
- 4 Wash tissue with 1X PBS for 5 minutes. Repeat x4. 5m
- 5 Incubate in secondary antibody for 2 hours at room temperature. 2h
- 6 Wash tissue with 1X PBS for 5 minutes. Repeat x4. 5m
- 7 Mount free-floating sections on SuperFrost+ slides (if staining free-floating tissue) and let dry at room temperature for 15 minutes. 15m
- 8 Coverslip with fluorescent mounting medium and #1.5 coverslips. Outline coverslip with clear nail polish and store at 4°C or -20°C depending on length of storage. 5m