



Jan 12, 2025

Chelex-based DNA extraction protocol for insect museum vouchers

DOI

<https://dx.doi.org/10.17504/protocols.io.bp2l6x54rlqe/v1>

Morgan Brown¹, Sara Ottati², Valeria Trivellone¹

¹University of Illinois at Urbana-Champaign, United States of America;

²Institute for Sustainable Plant Protection (IPSP), National Research Council, Italy



Valeria Trivellone

University of Illinois



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l6x54rlqe/v1>

Protocol Citation: Morgan Brown, Sara Ottati, Valeria Trivellone 2025. Chelex-based DNA extraction protocol for insect museum vouchers. [protocols.io https://dx.doi.org/10.17504/protocols.io.bp2l6x54rlqe/v1](https://dx.doi.org/10.17504/protocols.io.bp2l6x54rlqe/v1)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: January 31, 2024

Last Modified: January 12, 2025

Protocol Integer ID: 94500

Keywords: microorganisms, voucher, museum, acid nucleic extraction, DNA, next-generation sequencing, phytoplasma, bacteria, leafhoppers, bleaching, chelex, Column-based extraction, qPCR, insect, insect museum vouchers museum biorepository, insect museum, voucher specimen, taxonomic status, dna extraction, borne pathogens detection, taxonomic status of the associate, based dna extraction protocol, microbiome association, pathogens detection, museum specimen, tracking insect, insect, diverse community of microorganism, specimen, biodiversity, species, species on earth, phylogenetics, pathogen, high quality nucleic acid, dna, microorganism, described species, inferring evolutionary pathway, documenting association, nucleic acid, catalog record

Funders Acknowledgements:

U.S. National Science Foundation

Grant ID: DEB-2244871

Abstract

Museum biorepositories preserve precious material for scientific research in several fields such as biodiversity, phylogenetics, population genetics, toxicology, environmental monitoring and pathology.

About half of the 8.7 million described species on Earth are insects and harbor a diverse community of microorganisms representing different association types, such as parasitic, mutualistic, and commensal.

Voucher specimens may play an important role in tracking insect-microbiome associations.

Catalog records and voucher specimens are the primary source for: (1) tracking past associations and inferring evolutionary pathways; (2) documenting associations in no longer accessible collection sites; (3) revising the taxonomic status of the associates; (4) resampling the associations over time when new technologies become available.

Here we present a non-destructive, fast, inexpensive, non-toxic protocol for DNA extraction from insect museum specimens with the aim to obtain high quality nucleic acid and sufficient yield for vector-borne pathogens detection, quantification and characterization.

Image Attribution

Image source: Christopher H. Dietrich (University of Illinois at Urbana-Champaign)



Materials

⊗ UltraPure[®] SDS Solution, 10% **Thermo Fisher Catalog #24730020**

⊗ Chelex-100 resin **Bio-Rad Laboratories**

⊗ 95% EtOH

⊗ 1M Tris-HCl, pH 8.0 **Invitrogen - Thermo Fisher Catalog #15568025**

⊗ EDTA (0.5 M), pH 8.0 **Life Technologies Catalog #AM9260G**



Insect specimen selection and preparation

9m

- 1 Transfer the ethanol-preserved  insect specimen on a piece of paper towel to allow the ethanol to evaporate for a few minutes.

2m

Note

For dry insect voucher, re-hydrate the  insect specimen in a humidifying chamber

 Overnight

- 2 Place  insect specimen in a 1.5 mL centrifuge tube containing  500 μ L ethanol  90 % (v/v) for  00:02:00 .

2m

*

- 2.1 Discard the ethanol and add  500 μ L bleach  2.5 % (v/v) (NaOCl - The Clorox Co., Oakland, CA, USA) for  00:05:00 at  7.5 $^{\circ}$ C .

5m

Note

The insect body must be fully immersed in the bleach solution.

- 2.2 Discard bleach solution and rinse 3 X with  250 μ L of sterile water, discard the water each time.

10m

DNA Extraction

19h 19m

- 3 Place the  insect specimen into a fresh 1.5 mL centrifuge tube.
- 4 Pipet  50 μ L TES buffer and  2 μ L Proteinase K ( 10 mg/mL , Qiagen) into the tube.

5m

2m

**Note**

Prepare TES buffer (Tris-HCl [M] 20 millimolar (mM) pH 8.0 , EDTA

[M] 10 millimolar (mM) , SDS [M] 0.5 % (v/v))

For 10 mL TES add 0.2 mL EDTA [M] 0.5 Molarity (M) , 0.5 mL SDS

[M] 10 % (v/v) , 0.2 mL Tris-HCl [M] 1 Molarity (M) , 9.1 mL water (MilliQ)

- 5 Incubate 55 °C Overnight ;

18h



- 6 Centrifugate the sample for a few seconds to spin down the liquid. With a pipet, transfer 50 µL of the liquid sample into a fresh tube.

2m

Store the insect specimen as appropriate for further morphological studies.

- 7 Prepare [M] 25 Mass / % volume Chelex 100 suspension.

15m

Note

Add 2.5 g Chelex 100 resin into 10 mL water (MilliQ or NANOpure)

- 8 Place a flask containing a [M] 25 Mass / % volume Chelex 100 suspension on a magnetic agitator and mix.

2m

- 9 Pipet 40 µL of Chelex directly from the agitator into the tube with the sample.

2m

- 10 Incubate the tube at 55 °C for 00:30:00 . Do not shake the tube.

30m



- 11 Gently invert the tube 5 times.

2m



12 Allow the solution to cool down to  Room temperature .



15m

13 Centrifugate at  14000 rpm for  00:01:00 .



2m

14 Transfer the supernatant into a fresh 1.5 mL centrifuge tube.



2m

Safety information

The present protocol do not include toxic chemicals

Protocol references

1. Bleaching step

Arora, A. K., Pesko, K. N., Quintero-Hernández, V., Possani, L. D., Miller, T. A., & Durvasula, R. V. (2018). A paratransgenic strategy to block transmission of *Xylella fastidiosa* from the glassy-winged sharpshooter *Homalodisca vitripennis*. *BMC biotechnology*, 18, 1-10.

Greenstone, M. H., Weber, D. C., Coudron, T. A., Payton, M. E., & Hu, J. S. (2012). Removing external DNA contamination from arthropod predators destined for molecular gut-content analysis. *Molecular Ecology Resources*, 12(3), 464-469.

2. Chelex protocol modified from

Casquet, J., Thebaud, C., & Gillespie, R. G. (2012). Chelex without boiling, a rapid and easy technique to obtain stable amplifiable DNA from small amounts of ethanol-stored spiders. *Molecular ecology resources*, 12(1), 136-141.

Lienhard, A., & Schäffer, S. (2019). Extracting the invisible: obtaining high quality DNA is a challenging task in small arthropods. *PeerJ*, 7, e6753.